INSTRUCTIONS FOR USE FOR:

EXCLUDER[®]

ILIAC BRANCH ENDOPROSTHESIS

en

English

hu

Magyar

cs

Čeština

it

Italiano

da

Dansk

lt

Lietuvių

nl

Nederlands

no

Norsk

fi

Suomi

pl

Polski

fr

Français

pt

Português

de

Deutsch

es

Español

el

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
sv

Svenska

DIRECTORY TO LANGUAGE LOCATION

US ENGLISH (English).....	9	en
CZECH (Čeština).....	17	cs
DANISH (Dansk).....	25	da
DUTCH (Nederlands).....	33	nl
FINNISH (Suomi)	41	fi
FRENCH (Français)	49	fr
GERMAN (Deutsch).....	57	de
GREEK (Ελληνικά)	67	el
HUNGARIAN (Magyar)	77	hu
ITALIAN (Italiano)	85	it
LITHUANIAN (Lietuvių)	93	lt
NORWEGIAN (Norsk).....	101	no
POLISH (Polski)	109	pl
PORTUGUESE (Português)	119	pt
SPANISH (Español)	127	es
SWEDISH (Svenska)	135	sv

1

en - Figure 1: Iliac Branch Component (IBC); (Distal Iliac Diameters of 10, 12, or 14,5 mm)

Iliac Branch Component Radiopaque Markers:

- Two (2) short markers at proximal end
 - One (1) long and one (1) short marker at the endoprosthesis bifurcation level. The long mark denotes the internal iliac leg side location and orientation.
 - One (1) marker ring at the opening of the internal iliac leg hole.
 - One (1) short marker at the distal end of the external iliac leg.
1. Radiopaque Markers (2)
 2. Internal Iliac (Long) Radiopaque Marker
 3. Radiopaque Marker Ring
 4. Radiopaque Marker
 5. External Iliac (Short) Radiopaque Marker

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










TABLE OF CONTENTS

en

DESCRIPTION	10
INDICATIONS FOR USE	10
CONTRAINDICATIONS	10
WARNINGS AND PRECAUTIONS	10
PATIENT SELECTION AND TREATMENT	11
PATIENT COUNSELING INFORMATION	11
HOW SUPPLIED	11
CLINICAL USE INFORMATION	11
DIRECTIONS FOR USE	12
IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP	14
DEVICE-RELATED ADVERSE EVENT REPORTING	15
PATIENT TRACKING INFORMATION	15
DEFINITIONS	16

- **Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious surgical consequences or injury to the patient.**

INSTRUCTIONS FOR USE

GORE® EXCLUDER® ILIAC BRANCH ENDOPROSTHESIS

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious surgical consequences or injury to the patient.

DESCRIPTION

Iliac Branch Component and Internal Iliac Component

The GORE® EXCLUDER® Iliac Branch Endoprosthesis (IBE) provides endovascular treatment of common iliac artery aneurysms or aortoiliac aneurysms.

The GORE® EXCLUDER® Iliac Branch Endoprosthesis is comprised of two components, the Iliac Branch Component (IBC) (Figure 1) and the Internal Iliac Component (IIC) (Figure 2). The Iliac Branch Component (IBC) is a bifurcated iliac branch device with an external iliac leg and an internal iliac artery gate. An Internal Iliac Component (IIC) is used to extend into the internal iliac artery. The graft material is expanded polytetrafluoroethylene and fluorinated ethylene propylene (ePTFE and FEP), that is supported by nitinol (nickel titanium alloy) wire along its external surface. An ePTFE / FEP sleeve is used to constrain the endoprostheses on the leading end of their respective delivery catheters (Figures 3A and 3B).

The IBC is placed in the common iliac at a level wherein the internal iliac gate is at or proximal to the internal iliac artery. Deployment of the IBC initiates from the leading (aortic) end and proceeds toward the trailing end of the delivery catheter (Figure 3C). The IIC is introduced through the non-treatment femoral access side, (contralateral to the IBC delivery). The IIC delivery catheter (Figure 3D) is placed up-and-over the aortic bifurcation, through the internal iliac gate and into the internal iliac artery. Deployment of the IIC initiates from the trailing (common iliac) end and proceeds toward the leading (internal iliac) end of the delivery catheter. The ePTFE / FEP sleeve remains in situ between the endoprosthesis and the vessel wall.

The GORE® EXCLUDER® Iliac Branch Endoprosthesis is based on the design of the GORE® EXCLUDER® AAA Endoprosthesis; the graft materials are identical, maintaining the same luminal and abluminal ePTFE surfaces and materials specifications.

The GORE® EXCLUDER® Iliac Branch Endoprosthesis is to be used in conjunction with the GORE® EXCLUDER® AAA Endoprosthesis or GORE® EXCLUDER® Conformable Endoprosthesis and is not intended to be used on its own.

Trunk-Ipsilateral Leg Endoprosthesis and Contralateral Leg Endoprosthesis (Distal Iliac Diameters of 23 and 27 mm)

The GORE® EXCLUDER® or GORE® EXCLUDER® Conformable Trunk-Ipsilateral Leg and Contralateral Leg (distal diameters of 23 and 27 mm) Endoprostheses provide proximal seal and fixation for the GORE® EXCLUDER® Iliac Branch Endoprosthesis, and also provide endovascular treatment of infrarenal abdominal aortic aneurysms (AAAs). **For more information on use of these devices, see the GORE® EXCLUDER® AAA Endoprosthesis or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use.**

INDICATIONS FOR USE

Iliac Branch and Internal Iliac Components

The GORE® EXCLUDER® Iliac Branch Endoprosthesis (IBE) is intended to isolate the common iliac artery from systemic blood flow and preserve blood flow in the external iliac and internal iliac arteries in patients with a common iliac or aortoiliac aneurysm, who have appropriate anatomy, including:

1. Adequate iliac / femoral access
2. Minimum common iliac diameter of 17 mm at the proximal implantation zone of the IBE
3. External iliac artery treatment diameter range of 6.5 – 25 mm and seal zone length of at least 10 mm
4. Internal iliac artery treatment diameter range of 6.5 – 13.5 mm and seal zone length of at least 10 mm

Trunk-Ipsilateral Leg and Contralateral Leg Endoprosthesis Components

The Trunk-Ipsilateral Leg and Contralateral Leg Endoprostheses are intended to provide proximal seal and fixation to the GORE® EXCLUDER® Iliac Branch Endoprosthesis following deployment of the GORE® EXCLUDER® Iliac Branch Endoprosthesis. For more information on the Trunk-Ipsilateral Leg and Contralateral Leg Endoprosthesis Component indications for use and deployment, see the GORE® EXCLUDER® AAA Endoprosthesis and GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use.

Aortic Extender Endoprosthesis and Iliac Extender Endoprosthesis Components

The Aortic and Iliac Extender Endoprostheses can be used after deployment of the GORE® EXCLUDER® Iliac Branch and AAA Endoprostheses. These extensions are used when additional length and / or sealing for aneurysm exclusion is desired. For more information on Aortic Extender and Iliac Extender indications for use and deployment, see the GORE® EXCLUDER® AAA Endoprosthesis and GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use.

CONTRAINDICATIONS

The GORE® EXCLUDER® Iliac Branch Endoprosthesis is contraindicated in:

- Patients with known sensitivities or allergies to the device materials.
- Patients with a systemic infection who may be at increased risk of endovascular graft infection.

WARNINGS AND PRECAUTIONS

General

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious surgical consequences or injury to the patient.
- The GORE® Medical Device is designed for single use only; do not reuse device. Gore does not have data regarding reuse of this device. Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination. Reuse may result in infection, serious injury, or patient death.
- The long-term performance of stent-grafts has not been established. All patients should be advised this treatment modality requires long-term, regular follow-up to assess patients' health status and stent-graft performance. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should receive enhanced follow-up (see IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP).
- The GORE® EXCLUDER® Iliac Branch Endoprosthesis is not recommended in patients unable to undergo, or who will not be compliant with the necessary pre and post-operative imaging and follow-up described in IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP.
- Intervention or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms and / or persistent endoleak. An increase in aneurysm size and / or persistent endoleak may lead to aneurysm rupture.
- Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- The GORE® EXCLUDER® Iliac Branch Endoprosthesis should only be used by physicians experienced in vascular interventional techniques, and who have successfully completed the appropriate physician training program.

Patient Selection, Treatment, and Follow-Up

- The safety and effectiveness of the GORE® EXCLUDER® Iliac Branch Endoprosthesis have not been evaluated in the following patient populations:
 - traumatic aortic or iliac injury
 - leaking: pending rupture or ruptured aneurysms
 - mycotic aneurysms
 - pseudoaneurysms resulting from previous graft placement
 - revision of previously placed stent grafts
 - genetic connective tissue disease (e.g., Marfans or Ehlers-Danlos Syndromes)
 - concomitant thoracic aortic or thoracoabdominal aneurysms
 - inflammatory aneurysms
 - patients with active systemic infections
 - pregnant or nursing females
 - patients less than 21 years old
- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and / or tortuosity) should be compatible with vascular access techniques and the vascular introducer sheaths and accessories necessary to deliver the endoprostheses.
- Successful exclusion of the aneurysm(s) may be affected by significant thrombus and / or calcium at the distal iliac artery interfaces. Irregular calcium and / or plaque may compromise the fixation and sealing of the implantation sites.
- The GORE® EXCLUDER® Iliac Branch Endoprosthesis is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and post operative follow-up imaging.
- The GORE® EXCLUDER® Iliac Branch Endoprosthesis is not recommended in patients where weight and / or size compromises or prevents the necessary imaging requirements.
- The GORE® EXCLUDER® Iliac Branch Endoprosthesis is not recommended in patients with known sensitivities or allergies to the device materials.

Implant Procedure

- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not rotate any delivery catheters while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not rotate the Iliac Branch Component (IBC) delivery catheter beyond 360°. Delivery system damage and / or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

- **Do not rotate the Internal Iliac Component (IIC) delivery catheter during delivery, positioning or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.**
- Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheaths. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- **Do not attempt to reposition the endoprosthesis after complete deployment of the device. Vessel damage or device misplacement may result.**
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms.
- Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation and reintervention resulting in potential patient harms, see ADVERSE EVENTS.
- Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
- Incorrect deployment and / or movement of the endoprosthesis may require surgical intervention.
- Do not attempt to advance any GORE® EXCLUDER® Endoprostheses through smaller introducer sheaths than as recommended.

PATIENT SELECTION AND TREATMENT

(SEE WARNINGS AND PRECAUTIONS)

Individualization of Treatment

Gore recommends that the GORE® EXCLUDER® Iliac Branch Endoprosthesis sizing should follow the treatment range recommendations as described in **Tables 1 and 2**. The length to the internal iliac gate of the GORE® EXCLUDER® Iliac Branch Endoprosthesis, together with the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis and Contralateral Leg Endoprosthesis, should not exceed the distance from the most distal (lowest) major renal artery to the ostium of the internal iliac artery on the IBE treatment side. Lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters / lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

Additional considerations for patient selection include but are not limited to:

1. Patient's age and life expectancy
2. Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
3. Patient's suitability for open surgical repair
4. Patient's anatomical suitability for endovascular repair
5. The risk of aneurysm rupture compared to the risk of treatment with the GORE® EXCLUDER® Iliac Branch Endoprosthesis and GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis.
6. Ability to tolerate general, regional, or local anesthesia
7. Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and / or tortuosity) should be compatible with vascular access techniques and the vascular introducer sheaths and accessories necessary to deliver the endoprostheses.
8. Successful exclusion of the aneurysm(s) may be affected by significant thrombus and / or calcium at the distal iliac artery interfaces. Irregular calcium and / or plaque may compromise the fixation and sealing of the implantation sites.
9. For IBC, external iliac artery lengths of at least 30 mm of which at least 10 mm must be non-aneurysmal seal zone of 6.5 – 13.5 mm in diameter, or ≤ 25 mm in diameter if extending with Iliac Extender Endoprosthesis (see GORE® EXCLUDER® AAA Endoprosthesis Instructions For Use)..
10. For IIC, internal iliac artery length of at least 30 mm of which at least 10 mm must be non-aneurysmal seal zone of 6.5 – 13.5 mm in diameter.
11. Follow all required patient selection criteria set forth in the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instruction For Use.
12. Freedom from significant femoral / iliac artery occlusive disease that would impede inflow and/or outflow of stent-grafts.

The final treatment decision is at the discretion of the physician and patient.

PATIENT COUNSELING INFORMATION

The physician and patient should review the risks and benefits when discussing this endovascular device and procedure including:

1. Risks and differences between endovascular repair and open surgical repair
2. Potential advantages of traditional open surgical repair
3. Potential advantages of endovascular repair
4. The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair
5. Potential benefits of preservation of blood flow in the internal iliac artery.

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to post-operative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- The long-term safety and effectiveness of endovascular repair has not been established. Physicians should advise all patients that this treatment modality requires long-term, regular follow-up to assess patients' health status and stent-graft performance. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should receive additional follow-up. Patients should be counseled on the need for regular follow-up, even in the absence of obvious symptoms, e.g., pain, numbness, weakness (see IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP).
6. Regular follow-up including imaging of the device should be performed at least every 12 months for all patients and additional imaging for patients with known endoleaks or aneurysm enlargement for the duration of the implant (see IMAGING GUIDELINES AND POST- OPERATIVE FOLLOW UP).
 7. Physicians must advise all patients that it is important to seek prompt medical attention if he / she experiences signs or symptoms of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking, or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as pain, numbness, weakness in the legs; any back, chest, abdominal, or groin pain, dizziness, fainting, rapid heartbeat, or sudden weakness.
 8. Procedure related risks include cardiac, pulmonary, neurologic, bowel, and bleeding complications. Device related risks include occlusion, endoleak, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture and death (See POTENTIAL DEVICE OR PROCEDURE RELATED ADVERSE EVENTS).

HOW SUPPLIED

The GORE® EXCLUDER® Iliac Branch Endoprosthesis is preloaded on a delivery catheter and supplied sterile and non-pyrogenic.

Storage and Handling

1. Do not resterilize; for single use only.
2. Do not use if damaged or if sterile barrier has been compromised.
3. Do not use after the "use by" (expiration) date printed on the label.
4. Store in a cool, dry place.

CLINICAL USE INFORMATION

Physician Training Program

CAUTION: Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The GORE® EXCLUDER® Iliac Branch Endoprosthesis should only be used by physicians experienced in vascular interventional techniques, and who have successfully completed the appropriate physician training program.

The recommended skills / knowledge requirements for physicians using the GORE® EXCLUDER® Iliac Branch Endoprosthesis include certification of the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis, as well as those outlined below:

Patient selection:

1. Knowledge of the natural history of abdominal aortic (AAA), common iliac, and aortoiliac aneurysms, and co-morbidities associated with aneurysmal disease.
2. Knowledge of radiographic image interpretation, device selection and sizing

A multi-disciplinary team that has combined procedural experience with:

1. Femoral cutdown, arteriotomy, and repair
2. Percutaneous access and closure techniques
3. Non-selective and selective guidewire and catheter techniques
4. Fluoroscopic and angiographic image interpretation
5. Embolization
6. Angioplasty
7. Endovascular stent placement
8. Snare techniques
9. Appropriate use of radiographic contrast material
10. Techniques to minimize radiation exposure
11. Expertise in necessary patient follow-up modalities

Required Materials

1. Two 0.035" (0.89 mm) 'super stiff' guidewires, 145 cm or longer (180 cm is recommended)
2. 0.035" (0.89 mm) flexible guidewire, 160 cm or longer (260 cm is recommended)
3. Angiographic radiopaque marker catheter
4. Snare catheter
5. Contrast media
6. Syringe
7. Heparin and heparinized saline
8. 16 Fr introducer sheath (**Table 1**)
9. 12 Fr flexible reinforced introducer sheath (**Table 2**)
10. Large diameter, low pressure aortic balloon (monitor balloon volumes and pressures as recommended in balloon catheter Instructions For Use)
11. Percutaneous transluminal angioplasty (PTA) balloons (**Tables 1-4**)

12. Trunk-Ipsilateral Leg Endoprosthesis and Contralateral Leg Endoprosthesis and the required materials for implantation (see the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use)
13. GORE® EXCLUDER® Aortic Extender and Iliac Extender Endoprostheses as needed (see the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use)

Table 1: Iliac Branch Component Sizing Guide*

Part Number	Proximal IBC Diameter (mm)	Distal IBC Diameter ¹ (mm)	Overall Device Length (cm)	Length to Internal Iliac Gate (cm)	Intended External Iliac Vessel Diameter ¹ (mm)	Recommended Introducer Sheath ² (Fr)	Recommended Angioplasty Balloon Size (Distal) (mm x mm)
CEB231010	23	10	10	5.5	6.5 – 9	16	10 x 40
CEB231210	23	12	10	5.5	10 – 11	16	12 x 40
CEB231410	23	14.5	10	5.5	12 – 13.5	16	14 x 40

¹ Recommended endoprosthesis oversizing relative to the vessel diameter is approximately 7-35% in the external iliac vessel.

² GORE® Introducer Sheaths are recommended. The GORE® EXCLUDER® Iliac Branch Component is not compatible with introducer sheaths longer than 40 cm (total length including the hemostatic valve).

* Note: All dimensions are nominal.

Table 2: Internal Iliac Component Sizing Guide*

Part Number	IIC Distal Diameter ¹ (mm)	Overall Device Length ¹ (cm)	Intended Internal Iliac Vessel Diameter ² (mm)	Recommended Introducer Sheath ³ (Fr x cm)	Recommended Balloon Size for IBC-IIC Overlap (mm x mm)	Recommended Angioplasty Balloon Size (Distal) (mm x mm)
HGB161007	10	7	6.5 – 9	12 x 45	14 x 40	10 x 40
HGB161207	12	7	10 – 11	12 x 45	14 x 40	12 x 40
HGB161407	14.5	7	12 – 13.5	12 x 45	14 x 40	14 x 40

¹ 7 cm long Internal Iliac Component provides a maximum extension of 4.5 cm when placed into the Iliac Branch Component.

² Recommended endoprosthesis oversizing relative to the vessel diameter is approximately 7-35% in the internal iliac vessel.

³ Flexible Reinforced sheath.

* Note: All dimensions are nominal.

Table 3: Contralateral Leg Endoprosthesis* Sizing Guide for Bridging to the IBC

Common Iliac Artery Diameter at Proximal Landing Zone ¹ (mm)	Contralateral Leg Endoprosthesis Distal Diameter ² (mm)	Overall Device Lengths ³ (cm)	Recommended Angioplasty Balloon Size for IBC Overlap (mm x mm)
17 – 18	23	10, 12, 14	18 x 40
19 – 20	23	10, 12, 14	20 x 40
20 – 21.5	23	10, 12, 14	22 x 40
>21.5	27	10, 12, 14	24 x 40

¹ Treatment diameters reflect use of Contralateral Leg Endoprosthesis as bridging component to IBC only. For traditional use of Contralateral Leg Endoprosthesis to provide arterial apposition, see the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use.

² Recommended endoprosthesis oversizing relative to the IBC vessel diameter is 7-26%.

³ Labeled Contralateral Leg length includes 3 cm overlap in contralateral gate of Trunk-Ipsilateral Endoprosthesis, and 3 cm overlap in proximal end of Iliac Branch Endoprosthesis.

* Note: All dimensions are nominal. Please see GORE® EXCLUDER® AAA Endoprosthesis Instructions for Use.

DIRECTIONS FOR USE

Pre-Treatment Planning

- Determine accurate size of anatomy and proper size of Iliac Branch Component (Table 1), Internal Iliac Component (Table 2), Contralateral Leg Endoprostheses (Table 3), Trunk-Ipsilateral Leg Endoprosthesis and and, if necessary, Aortic and Iliac Extender Endoprostheses (see GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use).
- Use high resolution, non-contrast and contrast enhanced computerized tomography (CT / CTA) at ≤ 3 mm acquisition and reconstruction collimation.
- Use multiple-view Digital Subtraction Angiography (DSA) with a radiopaque marker catheter or spiral CT multi-planar reconstruction.
- For angiography, use correct imaging angulation (cranial-caudal, lateral-oblique) to accurately identify origin of branch vessel anatomy.
- Consider breath-hold technique to optimize Digital Subtraction Angiography (DSA) image quality.

Anatomical Requirements

- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and / or tortuosity) that is compatible with vascular access techniques and the labeled vascular introducer sheath size necessary to deliver the endoprostheses
- Successful exclusion of the aneurysm(s) may be affected by significant thrombus and / or calcium at the distal iliac artery interfaces. The US clinical studies for the GORE® EXCLUDER® AAA Endoprosthesis quantify significant thrombus as thrombus ≥ 2 mm in thickness and / or $\geq 25\%$ of the vessel circumference in the intended seal zone of the distal iliac artery. Irregular calcium and / or plaque may compromise the fixation and sealing of the implantation sites.
- For IBC, minimum common iliac diameter of ≥ 17 mm in diameter at the proximal IBC landing zone.
- For IBC, distal external iliac artery lengths of at least 30 mm of which at least 10 mm must be non-aneurysmal seal zone of 6.5 - 13.5 mm in diameter, or 6.5 – 25 mm in diameter if extending with Iliac Extender Endoprosthesis.
- For IIC, distal internal iliac artery length of at least 30 mm of which at least 10 mm must be non-aneurysmal seal zone of 6.5 - 13.5 mm in diameter.
- Freedom from significant femoral / iliac artery occlusive disease that would impede inflow or outflow of stent-grafts.
- Ability to tolerate general, regional, or local anesthesia.
- Patient's anatomical suitability for endovascular repair.

Arterial Access and Angiography

- Following standard practices, access the intended iliac treatment side via a percutaneous diagnostic sheath, and perform marker catheter Digital Subtraction Angiography (AP, oblique and lateral views as necessary) to confirm the correct device component sizing, and deployment locations. Consider breath-hold technique to optimize image quality. Leave marker catheter in place at the level of the renal arteries.
- Use an accurate radiopaque patient marking method to assure accurate device positioning and deployment locations.
- Following standard practices, perform percutaneous access and / or surgical exposure of the vessels selected to receive the Trunk-Ipsilateral Leg Endoprosthesis and IBE introducer sheaths.
- Following the manufacturer's Instructions For Use, advance through the IBE treatment side a 0.035" (0.89 mm) 'super stiff' guidewire or acceptable equivalent to the level of the renal arteries exchanging through an angiographic catheter.
- Following the manufacturer's Instructions For Use, prepare and advance the recommended 16 Fr introducer sheath over the guidewire, through the ilio-femoral anatomy and up to the level of the distal aortic neck according to standard practice.
- CAUTION:** Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Following the manufacturer's Instructions For Use, advance through the non-IBE treatment side a 0.035" (0.89 mm) 'super stiff' guidewire or acceptable equivalent to the level of the renal arteries exchanging through an angiographic catheter.
- Following the manufacturer's Instructions For Use, prepare and advance the recommended 12 Fr flexible reinforced introducer sheath over the guidewire, through the ilio-femoral anatomy and up to the level of the aortic bifurcation and remove the 0.035" (0.89 mm) 'super stiff' guidewire.
- Following the manufacturer's Instructions For Use, advance through the recommended 16 Fr introducer sheath a second wire, an 0.035" (0.89 mm) flexible guidewire of at least 180 cm in length (260 cm is recommended) to the level of the aortic bifurcation for snaring (**Figure 4**).
- Following the manufacturer's Instructions For Use, prepare and advance a snare catheter through, and snare and pull through the flexible guidewire, thus generating a flexible throughwire (**Figure 4**).
- Use standard heparinized saline to flush the sheath to prevent thrombus formation.

Device Preparation

- CAUTION:** Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Remove the appropriately sized device delivery catheters from their packaging and examine for possible damage.
- Remove protective packaging cover(s) from the leading end of the delivery catheters, and remove the stainless steel packaging mandrels from 1) the delivery catheter (for IBC and IIC) and 2) the removable guidewire tube (RGT) (IBC only). **CAUTION:** Do NOT remove the RGT at this point. The RGT provides pre-cannulation of the internal iliac gate, and removal of the RGT prior to cannulating with the throughwire (step 2 below) will remove this feature.
- Flush the flushing port on the trailing end of the delivery catheter with heparinized saline (**Figures 3A and 3B**).

Iliac Branched Component Positioning and Deployment

- Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
- Advance the Iliac Branch Component (IBC) device catheter over the 0.035" (0.89 mm) 'super stiff' guidewire, and cannulate the removable guidewire tube (RGT) with the 0.035" (0.89 mm) flexible throughwire (**Figure 5**). Remove the RGT leaving the throughwire in the constrained endoprosthesis (**Figure 6**).
- Advance the IBC delivery catheter over both guidewires, through the recommended 16 Fr introducer sheath into the common iliac artery to the approximate level of intended landing zone (**Figure 7**). **WARNING:** Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS. **WARNING:** Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS. **WARNING:** Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- While maintaining the delivery catheter in position, withdraw the introducer sheath (**Figure 7**) and visually verify that the leading end of the introducer sheath is not covering the Iliac Branch Component (IBC). **CAUTION:** Visually verify the throughwire is not wrapped around the IBC guidewire or delivery catheter. If wire wrap is observed, rotate the device and delivery catheter to resolve.
- Magnify and center the fluoroscope to the level of the common iliac artery. Reposition and rotate the Iliac Branch Component delivery catheter as necessary to properly position the device in the correct orientation in regard to the anatomy. The long radiopaque marker should be oriented toward the internal

iliac artery and the internal iliac leg radiopaque ring should be proximal to the origin of the internal iliac artery in order to obtain internal iliac access.

WARNING: Do not rotate the IBC delivery catheter beyond 360°. Delivery system damage and / or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

6. Stabilize the Iliac Branch Component (IBC) delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.
WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
7. Loosen the White Outer Deployment Knob (Figure 8A). Confirm final device position and orientation and deploy the IBC using a steady and continuous pull of the deployment knob to release the proximal endoprosthesis past the internal iliac leg hole of the IBC. Pull the deployment knob straight out and away from the catheter handle (Figure 8B). Deployment initiates from the leading (aortic) end toward the trailing (iliac) end. The external iliac leg will remain constrained on the delivery catheter (Figure 8C). If the Iliac Branch Component (IBC) is not properly positioned following removal of the White Outer Deployment Knob, continue to Section: Optional Repositioning of Iliac Branch Component Endoprosthesis.
WARNING: Do not attempt to remove a partially deployed Iliac Branch Component (IBC).
WARNING: Do not loosen and remove the Gray Inner Deployment Knob (Figure 3) until ready to proceed to External Iliac Leg Deployment of Iliac Branch Component (IBC) (Figure 12). Removal of the Gray Inner Deployment Knob will deploy the constrained external iliac leg of the Iliac Branch Component (IBC) and the device will no longer repositionable.

Optional Repositioning of Iliac Branched Component

- A. Magnify and center the fluoroscopic image on the IBC. Manually reposition by rotating the Iliac Branch Component (IBC) up to 90° in either direction as necessary to properly position the device to access the internal iliac artery (Figure 8C). The long radiopaque marker and internal iliac leg radiopaque ring should be oriented toward the internal iliac artery. The IBC can also be repositioned distally to bring the internal iliac gate in proximity to the internal iliac ostium.

CAUTION: Excessive repositioning of the IBC may result in tissue damage and/or thrombus dislodgment. Repositioning IBC in a proximal direction after partial deployment is not recommended. The internal iliac leg hole should be kept proximal to the internal iliac ostium.

Flexible Sheath Positioning and Internal Iliac Artery Cannulation through Internal Iliac Leg Hole

1. Use fluoroscopic visualization for all guidewire, sheath, and device catheter manipulations.
2. The internal iliac leg hole is pre-cannulated via the throughwire as per steps 1 through 3 above, Iliac Branch Component Positioning.
3. Following manufacturer's Instructions For Use, introduce a recommended 12Fr flexible reinforced sheath of 45 cm in length over the 0.035" (0.89 mm) throughwire via femoral access contralateral to the IBE treatment side.
4. Advance the flexible reinforced sheath along the throughwire up and over the aortic bifurcation, into and through the proximal opening of the Iliac Branch Component (IBC) to the level of the internal iliac leg hole (Figure 9).
5. Following manufacturer's Instructions For Use, introduce and advance an appropriate additional 0.035" (0.89 mm) guidewire (IIC guidewire) (and catheter if necessary) through the 12 Fr flexible reinforced sheath and internal iliac leg hole, and cannulate the internal iliac artery (Figure 10), and exchange for a 0.035" (0.89 mm) 'super stiff' guidewire if necessary.
NOTE: It is recommended to maintain the throughwire throughout IIC deployment to provide sheath stability.
CAUTION: Removal of the throughwire may result in loss of internal iliac leg hole cannulation.
6. Visually verify that the IIC guidewire is within the internal iliac artery by standard practice, such as angiography.
NOTE: It is recommended to maintain the throughwire throughout IIC deployment to provide sheath stability.
CAUTION: Removal of the throughwire may result in loss of internal iliac leg hole cannulation.

Internal Iliac Component Positioning and Deployment

1. Prepare the device delivery catheter as detailed in the "Device Preparation" section, and cannulate the Internal Iliac Component (IIC) device catheter lumen over the 0.035" (0.89 mm) IIC guidewire, and advance the delivery catheter through the 12 Fr flexible reinforced introducer sheath and into the internal iliac artery (Figure 11).
WARNING: Do not advance the device outside of the sheath while tracking it to the internal iliac artery. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not rotate the delivery device catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of the resistance. Vessel or catheter damage or premature

deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

2. Align the radiopaque marker of the anatomically proximal (catheter trailing) end of the IIC device with the long radiopaque marker on the partially deployed IBC. With the alignment of these markers, an approximate 2.5 cm overlap will be achieved.
3. While maintaining the delivery catheter in position, withdraw the introducer sheath and visually verify that the leading end of the introducer sheath is not covering the Internal Iliac Component (IIC).
WARNING: Do not rotate the Internal Iliac Component (IIC) delivery catheter during delivery, positioning, or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
4. Stabilize the Internal Iliac Component (IIC) delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.
WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
5. Loosen the deployment knob. Visually confirm final device position. Deploy the Internal Iliac Component (IIC) by using a steady, continuous pull of the deployment knob to release the endoprosthesis (Figure 12). Pull the deployment knob straight out and away from the catheter side-arm. Deployment initiates from the trailing (common iliac) end toward the leading (internal iliac) end.
WARNING: Do not attempt to reposition the endoprosthesis after deployment has been initiated. Vessel damage or device misplacement may result.
WARNING: Incorrect deployment or movement of the endoprosthesis may require surgical intervention.
6. Use fluoroscopic guidance during withdrawal of the delivery catheter to assure safe removal from the patient. If resistance is felt during removal of the delivery catheter through the introducer sheath, stop and assess cause of resistance. If necessary, withdraw delivery catheter and introducer sheath together.
WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation and reintervention resulting in potential patient harms, see ADVERSE EVENTS.
7. When withdrawing the device delivery system through the introducer sheath, verify all catheter components are intact.
WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
8. If extension of the Internal Iliac Component (IIC) is required for adequate seal, repeat steps 1 through 6 with an additional Internal Iliac Component (IIC). A minimum overlap of 3 cm is required between Iliac Extender Endoprostheses.
9. Following manufacturer's Instructions For Use, advance and inflate a 14 mm PTA balloon catheter to seat the proximal end of the Internal Iliac Component (IIC) within the internal iliac artery leg hole overlap region. Follow manufacturer's recommended method for preparation and use of iliac dilation balloons, carefully monitoring both volume and pressure to avoid complications.
10. With the 14 mm PTA balloon in position, carefully remove the throughwire.
11. Following manufacturer's Instructions For Use, advance and inflate the appropriate size PTA balloon (Table 2) to seat the distal end of the Internal Iliac Component (IIC) within the internal iliac artery. Follow manufacturer's recommended method for preparation and use of iliac dilation balloons, carefully monitoring both volume and pressure to avoid complications.

External Iliac Leg Deployment of Iliac Branched Component

1. Loosen the Gray Inner Deployment Knob by turning 90° counter-clockwise (Figure 13A). Deploy the external iliac leg of the IBC using a steady and continuous pull of the deployment knob straight out and away from the catheter handle to release the endoprosthesis (Figures 13B and 13C).
2. Use fluoroscopic guidance during withdrawal of the delivery catheter to assure safe removal from the patient. If resistance is felt during removal of the delivery catheter through the introducer sheath, stop and assess cause of resistance. If necessary, withdraw delivery catheter and introducer sheath together.
WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation and reintervention resulting in potential patient harms, see ADVERSE EVENTS.
3. When withdrawing the device delivery system through the introduce sheath, verify all catheter components are intact.

WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).

- Following manufacturer's Instructions For Use, advance and inflate the appropriate size PTA balloon (**Table 1**) to seat the distal end of the external iliac leg of the Iliac Branch Component within the external iliac artery. Following manufacturer's recommended method for size selection, preparation, and use of PTA balloons. Carefully inflate the balloon to avoid complications.

WARNING: Do not attempt to reposition the endoprosthesis after complete deployment of the device. Vessel damage or device misplacement may result.

WARNING: Incorrect deployment or movement of the endoprosthesis may require surgical intervention.

Trunk-Ipsilateral Leg Endoprosthesis Positioning and Deployment

NOTE: Please consult the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use (IFU) for information pertaining to the Trunk-Ipsilateral Leg Endoprosthesis.

- Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
- From the side contralateral to the GORE® EXCLUDER® Iliac Branch Endoprosthesis (IBE) treatment side, withdraw the IIC guidewire to the level of the native aortic bifurcation, then advance the IIC guidewire in the aorta past the renal arteries.
- Following the manufacturer's Instructions For Use, prepare and advance the appropriate introducer sheath over the guidewire, through the ilio-femoral anatomy and into the aortic neck according to standard practice.
- Position and deploy the Trunk-Ipsilateral Leg Endoprosthesis as described in GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use (**Figure 14**).
- Cannulate the contralateral leg hole as described in the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use (**Figure 14**).
- Visually verify the length from long radiopaque marker on the Trunk-Ipsilateral Leg Endoprosthesis to long radiopaque marker on the IBC before proceeding to Section: *Contralateral Leg Endoprosthesis as Bridge to IBC – Positioning and Deployment*.

Contralateral Leg Endoprosthesis as Bridge to IBC - Positioning and Deployment

NOTE: Please consult the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use (IFU) for information pertaining to the Contralateral Leg Endoprosthesis.

CAUTION: Only 23 mm or 27 mm distal diameter Contralateral Leg Endoprosthesis is compatible with the Iliac Branch Component (IBC). Use of a Contralateral Leg Endoprosthesis with a distal diameter <23 mm may result in potential adverse events.

WARNING: Incorrect deployment of the Endoprosthesis may require surgical intervention.

- Following manufacturer's Instructions For Use, prepare the device delivery catheter, and advance the delivery catheter through the appropriate introducer sheath to the level of the long radiopaque marker.

WARNING: Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

WARNING: Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- While maintaining the delivery catheter in position, withdraw the introducer sheath, and visually verify that the leading end of the introducer sheath is below the Contralateral Leg Endoprosthesis.
- Align the radiopaque marker of the proximal end of the Contralateral Leg Endoprosthesis with the long radiopaque marker on the Trunk-Ipsilateral Leg Endoprosthesis. With the alignment of these markers, an approximate 3 cm overlap will be achieved.
- Visually verify that the distal Contralateral Leg radiopaque marker is aligned with the proximal edge of the long and short radiopaque markers of the Iliac Branch Component. (IBC). Alignment of these markers will achieve an approximate 3 cm overlap.

WARNING: Deployment below these radiopaque markers may result in the Contralateral Leg Endoprosthesis deploying within the external iliac leg of the IBC, leading to insufficient blood flow into the internal iliac artery.
- Deploy the Contralateral Leg Endoprosthesis according to the GORE® EXCLUDER® AAA or GORE® EXCLUDER® Conformable Endoprosthesis Instructions For Use (IFU) (**Figure 15**).

- Following manufacturer's Instructions For Use, advance and inflate the appropriate size PTA balloon (**Table 3**) to seat the distal end of the Contralateral Leg Endoprosthesis within the proximal end of the Iliac Branch Component. Carefully inflate the balloon to avoid complications.

Completion of the Procedure

- Perform extended imaging angiography to confirm exclusion of the aneurysm(s). Consider breath-hold technique to optimize Digital Subtraction Angiography (DSA) image quality. Please consult the GORE® EXCLUDER® AAA Endoprosthesis Instructions For Use (IFU) for information pertaining to the use of Aortic and Iliac Extenders. These extensions can be used when additional length and / or sealing for aneurysmal exclusion is desired.
- Prior to removing wires and sheaths verify all catheter components have been withdrawn from the patient.

WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
- Close arterial access according to standard practice.
- Follow up patients as necessary to provide proper surveillance of the long-term performance of the endoprosthesis, procedure, and status of the aneurysm (See IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP).

IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

General

The long-term safety and effectiveness of endovascular repair has not been established. All patients should be advised this treatment modality requires long-term, regular follow-up to assess patients' health status and stent-graft performance. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should receive additional follow-up. Patients should be counseled on the need for regular follow-up, even in the absence of obvious symptoms (e.g., pain, numbness, weakness).

Regular and consistent follow-up is a critical part of ensuring continuing safety and efficacy of aortic endovascular repair. Physicians should tailor patient follow up to the needs and circumstances of each individual patient. At least one annual physician visit and the imaging schedule is recommended (**Table 4**)

Follow up modalities include CT / CTA, multi-view abdominal X-ray, MRI / MRA, and ultrasound. Data from these modalities is acquired and used to compare baseline and subsequent exams to review devices and morphological changes over time and their effects on exclusion of the aneurysm.

- CT / CTA imaging provides information on aneurysm size, vascular morphological changes, proximal device-Trunk fixation and migration, endoleak and patency / limb occlusion.
- Multi-view device X-ray imaging provides information on device wireform integrity (e.g., fracture, kinking) and relative component migration.
- MRI / MRA imaging provides information similar to CT / CTA and is often used as a surrogate for CT / CTA if patients are CT contrast intolerant.
- Ultrasound may be used to assess for endoleak and aneurysm size status but not for device integrity, specifically the wire form. Ultrasound is a less reliable and sensitive diagnostic method compared to CT.

Alternative imaging recommendations for patients with CT or angiography contrast intolerance issues include CO2 angiography, MRI-MRA with or without contrast, and ultrasound. These imaging and surveillance modalities may be less sensitive and difficult to compare with diagnostic findings from previous or subsequent follow-up exams.

Table 4. Recommended Schedule for Patient Imaging Follow Up

Visit	Schedule for Patient Imaging Follow-Up		
	Angiogram	Abdominal X-Ray ¹	CT Pre-Contrast and Contrast
Pre-Treatment	X ²		X ²
Treatment (Pre and Post Deployment)	X		
Discharge		X	
1 Month			X
3 Month			X ³
6 Month		X	X
12 Month (Annually Thereafter)		X	X

¹ Recommended if wire fracture is suspected

² Imaging should be performed ≤ six months prior to procedure

³ Recommended if endoleak and/or aneurysm enlargement reported at one month

Angiographic Imaging

Angiographic images are recommended pre-treatment to evaluate the length and tortuosity of abdominal aorta, iliac and common femoral arteries.

- Images should include an angiographic marker catheter with incremental one centimeter radiopaque marks over a 10 to 20 cm length.
- The following views are recommended for optimal evaluation and case planning:
 - Abdominal aorta; Supine-AP, Lateral
 - Pelvis (to include bilateral common femorals); AP, both obliques

Angiographic images are recommended during the treatment procedure both pre and post-deployment to evaluate device placement and orientation. Selective angiography during subsequent follow-up exams may provide useful device position and device integrity information.

CT / CTA Images

- Film sets should include all sequential images at lowest possible slice thickness (≤ 2 mm). Do NOT perform large slice thickness (> 3 mm) and / or omission of CT images / film sets (non-consecutive) as it prevents precise anatomical and device comparisons over time.
- If an endoleak is suspected or there is aneurysm enlargement, it is recommended that pre-contrast and contrast runs be performed.
- Pre-contrast and contrast run slice thickness and interval must match.
- DO NOT change patient orientation or re-landmark patient between non contrast and contrast runs.

Non-contrast and contrast enhanced baseline and follow-up exams are important for optimal patient surveillance. The following CT / CTA imaging guidelines are recommended (Table 5).

Table 5. CT / CTA Imaging Guidelines

	Pre-Contrast	CT / CTA
IV Contrast	No	Yes
Injection Volume (ml)	NA	150
Injection Rate (cc / sec)	NA	≥ 2.5
Delay	NA	Smart-Prep*, CARE or equivalent
Start Position	Diaphragm	1 cm above Celiac Axis
End Position	Proximal Femur	Femoral Bifurcation
Scan FOV	Large	Large
DFOV	32 cm	32 cm
Scan Type	Helical	Helical
Rotation Speed	0.8	0.8
Slice Thickness (mm)	≤ 2.0 mm	≤ 2.0 mm
Scan Mode	HS	HS
Table Speed (mm / rot)	15	15
Interval (mm)	2.0	2.0
KV / mA	120 / 300	120 / 300
Reconstruction / Algorithm	≤ 3.0 mm Soft	≤ 3.0 mm Soft
* Smart Prep	ROI Loc: 1 cm Sup. to Celiac Axis Scan Phase: 3 Sec MA: 40	Monitor Delay: 6 Sec Monitor ISD: 3 Sec Enhance Thres: 100 HU

Abdominal X-Ray Film Series (plain film)

If there is any concern regarding device integrity (e.g., kinking, stent-wire breaks, relative component migration), a chest X-ray film series may be acquired and evaluated by the attending physician. The following chest X-ray views are recommended for optimal visualization of the endoprosthesis. Magnified views (2-4x) may aid in evaluation of device integrity.

- Supine – frontal (AP)
- Lateral
- 45° LPO
- 45° RPO

Ensure the entire device is captured on each single image format.

MRI Information MR Conditional

Non-clinical testing has demonstrated that the GORE® EXCLUDER® Iliac Branch Endoprosthesis is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 1.5 or 3.0 Tesla
- Highest spatial gradient magnetic field of less than or equal to 3000 Gauss / cm
- Whole-body-averaged specific absorption rate (SAR) of 2.0 W / kg for 15 minutes of scanning (i.e., per pulse sequence)
- Normal Operating Mode for MR system

MR Related Heating:

In non-clinical testing, the IBE produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numeris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and during MRI performed for 15-min of scanning 3-Tesla/128-MHz (Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	1.5 Tesla	3 Tesla
MR system reported, whole body averaged SAR	2.9 W/kg	2.9 W/kg
Calorimetry measured values, whole body averaged SAR	2.1 W/kg	2.7 W/kg
Highest recorded temperature change	1.9°C	2.3°C
Calculated Temperature Change Normal Operating Mode SAR = 2.0 W/kg	1.3°C	1.6°C
Calculated Temperature Change First Level Control Mode SAR = 4.0 W/kg	2.6°C	3.2°C

Image Artifact:

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the IBE. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 5 mm relative to the size and shape of this implant.

Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement ≥ 5 mm of maximum diameter (regardless of endoleak status)
- Branch vessel stenosis or occlusion

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy, and the patient's personal choices. Patients should be counseled as to the possibility of subsequent reinterventions including catheter-based and open surgical conversion.

Potential Device or Procedure Related Adverse Events

Adverse events that may occur and / or require intervention or additional procedure time include, but are not limited to:

- allergic reaction and / or anaphylactoid response to x-ray contrast dye, anti-platelet therapy, device materials
- amputation
- anesthetic complications
- aneurysm enlargement
- aneurysm rupture and death
- arterial or venous thrombosis and / or pseudoaneurysm
- arteriovenous fistula
- bleeding, hematoma, or coagulopathy
- bowel (e.g., ileus, transient ischemia, infarction, necrosis)
- cardiac (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension or hypertension)
- claudication (e.g., buttock, lower limb)
- death
- dissection, perforation, or rupture of the aortic vessel and surrounding vasculature
- edema
- embolization (micro and macro) with transient or permanent ischemia
- endoleak
- endoprosthesis or delivery system: improper component placement; incomplete component deployment; unintentional/premature component deployment; leading end catheter component retention; component migration; separation of graft material from stent; occlusion; infection; stent fracture; graft material failure, dilatation, erosion, puncture, perigraft flow
- fever and localized inflammation
- genitourinary (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- hepatic failure
- impotence
- infection (e.g., aneurysm, device or access sites)
- lymph fistula / complications
- multi-system organ failure
- neurologic damage, local or systemic (e.g., stroke, paraplegia, paraparesis)
- occlusion of device or native vessel
- post-implant syndrome
- pulmonary complications (e.g., pneumonia, respiratory failure)
- radiation injury, late malignancy
- renal (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- surgical cut down, bypass or conversion
- tissue necrosis
- wound (e.g., infection, dehiscence)
- vascular spasm or vascular trauma (e.g., ilio-femoral vessel dissection, bleeding, rupture, death)

DEVICE-RELATED ADVERSE EVENT REPORTING

Any adverse event involving the GORE® EXCLUDER® Iliac Branch Endoprosthesis should be reported to W. L. Gore & Associates immediately. To report an event in the US, call 800.437.8181 (as required by US Federal Regulation). Outside the US, contact your local Gore Representative.

PATIENT TRACKING INFORMATION

In addition to these Instructions For Use, the GORE® EXCLUDER® Iliac Branch Endoprosthesis is packaged with a Device Tracking Form which US hospital staff are required to complete and forward to W. L. Gore & Associates for the purposes of tracking all patients who receive a GORE® EXCLUDER® Iliac Branch Endoprosthesis product (as required by US Federal Regulation).

DEFINITIONS

 Authorised Representative in the European Community

 Batch Code

 Caution

 **Only** CAUTION: USA Federal Law restricts the sale, distribution, or use of this device to, by, or on the order of a physician.

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© 2013-2015, 2017, 2019, 2020 W. L. Gore & Associates, Inc.
GORE®, EXCLUDER®, and designs are trademarks of W. L. Gore & Associates.
All other trademarks are the property of their respective owners.
W. L. GORE & ASSOCIATES B.V.