INSTRUCTIONS FOR USE FOR:

EXCLUDER[®]

A A A E N D O P R O S T H E S I S

us

US English

Low Permeability Design

en

English

Low Permeability Design

cs

Čeština

Model s nízkou propustností

da

Dansk

Design med lav permeabilitet

nl

Nederlands

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fi

Suomi

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fr

Français

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de

Deutsch

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

1A • 1B

us en - Figure 1A. Trunk-Ipsilateral Leg Endoprosthesis;
(Aortic Diameters of 23, 26 or 28.5 mm)

us en - Figure 1B. Trunk-Ipsilateral Leg Endoprosthesis;
(Aortic Diameter of 31 mm or 35 mm)

Trunk-Ipsilateral Leg Endoprosthesis Radiopaque Markers:

- Three (3) short markers at the aortic end.
- One (1) long and one (1) short marker at the endoprosthesis bifurcation level. The long marker denotes the contralateral leg side location and orientation.
- One (1) marker ring at the opening of the contralateral leg hole.
- One (1) short marker at the iliac end of the ipsilateral leg.

1. Radiopaque Markers (3)
2. Anchors
3. Ipsilateral (Short) Radiopaque Marker
4. Radiopaque Marker Ring
5. Radiopaque Marker
6. Sealing Cuff
7. Contralateral (Long) Radiopaque Marker

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
2B

us en - Figure 2B. Contralateral Leg Endoprosthesis;
(Distal Iliac Diameters of 16, 18, 20, 23, and 27 mm)

Contralateral Leg Endoprosthesis Radiopaque Markers:

- One (1) marker at each end
- One (1) marker located 3 cm below the proximal end
- One (1) marker located 4 cm above the distal end (16, 18, and 20 mm Contralateral Leg devices only)
- One (1) marker located 5 cm above the distal end (23 mm Contralateral Leg device only)
- One (1) marker located 6 cm above the distal end (27 mm Contralateral Leg device only) *

1. Leading End
2. Trailing End
3. Radiopaque Marker

* Note: The 27 mm x 10 cm does not have a radiopaque marker 6 cm above the distal end.

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
us en - Figure 3E. Constrained GORE® EXCLUDER® AAA Endoprosthesis (Contralateral Distal Iliac Diameters of 16, 18, 20, 23, and 27 mm) on Delivery Catheter with Radiopaque Markers

1. Leading End
2. Radiopaque Marker
3. Radiopaque Marker on Delivery Catheter

* Note: the 27 mm x 10 cm Contralateral Leg Endoprosthesis has three (3) radiopaque markers; One (1) at each end and one (1) located 3 cm below the proximal end.

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
7A • 7B • 7C

us en - Figure 7. Contralateral Leg Endoprosthesis as an Iliac Extender (Distal Iliac Diameters of 16, 18, 20, 23, and 27 mm)

1. A minimum overlap of 3 cm with the Ipsilateral Leg is required.
2. If the Contralateral Leg and Iliac Extender diameters are identical, the taper zone can be deployed inside the previously deployed Contralateral Leg.
3. A minimum overlap of 3 cm with the Contralateral Leg is required. This overlap should be achieved prior to the beginning of the distal taper zone of the 18, 20, 23, and 27 mm Contralateral Leg.

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TABLE OF CONTENTS

us

DESCRIPTION.....	22
INDICATIONS FOR USE	22
CONTRAINDICATIONS.....	22
WARNINGS AND PRECAUTIONS.....	22
CLINICAL STUDY OUTCOMES* (98-03 and 99-04)	23
CLINICAL STUDY OUTCOMES (04-04).....	27
CLINICAL STUDY OUTCOMES (03-02).....	28
PATIENT SELECTION AND TREATMENT.....	28
PATIENT COUNSELING INFORMATION	29
HOW SUPPLIED.....	29
CLINICAL USE INFORMATION	29
DIRECTIONS FOR USE	30
IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP.....	33
DEVICE-RELATED ADVERSE EVENT REPORTING.....	34
PATIENT TRACKING INFORMATION	34
DEFINITIONS	34

- **Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious surgical consequences or injury to the patient.**

INSTRUCTIONS FOR USE

GORE® EXCLUDER® AAA ENDOPROSTHESIS

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious surgical consequences or injury to the patient.

DESCRIPTION

Trunk-Ipsilateral Leg Endoprosthesis and Contralateral Leg Endoprosthesis

The GORE® EXCLUDER® AAA Endoprosthesis provides endovascular treatment of infrarenal abdominal aortic aneurysms (AAAs).

The GORE® EXCLUDER® AAA Endoprosthesis is comprised of two components, the Trunk-Ipsilateral Leg Endoprosthesis (Trunk) (Figures 1A and 1B) and the Contralateral Leg Endoprosthesis (Figures 2A and 2B). The graft material is expanded polytetrafluoroethylene and fluorinated ethylene propylene (ePTFE and FEP), that is supported by nitinol (nickel titanium alloy) wire along its external surface. Nitinol anchors and an ePTFE / FEP sealing cuff are located at the aortic end of the trunk (Figures 1A and 1B). An ePTFE / FEP sleeve is used to constrain the endoprosthesis on the leading end of the delivery catheter (Figures 3A, 3B, 3C, 3D and 3E).

Deployment of both endoprosthesis components initiates from the leading (aortic) end and proceeds toward the trailing (iliac) end of the delivery catheter (Figures 3A, 3B, 3C, 3D and 3E). The ePTFE / FEP sleeve remains in situ between the endoprosthesis and the vessel wall.

The low permeability GORE® EXCLUDER® Endoprosthesis is the only available design, and differs from the original GORE® EXCLUDER® Endoprosthesis by the addition of film layers that decrease the overall permeability of the graft. The luminal and abluminal ePTFE surface materials, microstructure, and characteristics of the low permeability GORE® EXCLUDER® Endoprosthesis meet the same specifications as the original GORE® EXCLUDER® Endoprosthesis.

Aortic Extender Endoprosthesis and Iliac Extender Endoprosthesis Components**Aortic Extender Endoprosthesis**

The Aortic Extender Endoprosthesis (Aortic Extender) provides an extension of approximately 1.6 cm or 2.2 cm of the leading (proximal) end of the Trunk-Ipsilateral Leg Endoprosthesis (Trunk). The Aortic Extender can be placed at variable extension lengths from 0 cm to 1.6 or 2.2 cm of the leading (proximal) end of the Trunk-Ipsilateral Leg Endoprosthesis (Trunk), allowing customization of extender treatment length based on patient anatomy and physician preference. This extension requires a minimum of approximately 1.6 cm or 2.2 cm overlap with the Trunk. The graft material is expanded polytetrafluoroethylene and fluorinated ethylene propylene (ePTFE and FEP), and is supported by nitinol (nickel titanium alloy) wire along its external surface. An ePTFE / FEP sealing cuff is located near the proximal end of the endoprosthesis (Figure 4). An ePTFE / FEP sleeve is used to constrain the endoprosthesis on the leading end of the delivery catheter (Figures 5A and 5B). Deployment of the Aortic Extender initiates from the trailing (trunk) end and proceeds toward the leading (aortic) end of the endoprosthesis and delivery catheter. Following deployment, the ePTFE / FEP sleeve remains *in situ* between the endoprosthesis and the vessel wall.

Iliac Extender Endoprosthesis**(Distal Iliac Diameters of 10, 12, and 14.5 mm)**

The Iliac Extender Endoprosthesis (Iliac Extender) provides an extension of up to 4 cm of either the ipsilateral or contralateral limb. The extender component can be placed at variable extension lengths from 0 cm to 4 cm allowing customization of extender treatment length based on patient anatomy and physician preference. The graft material is ePTFE / FEP, and is supported by nitinol wire along its external surface. A radiopaque marker is located 3 cm from the proximal or top end (Figure 6A). This marker denotes the recommended minimum overlap with the ipsilateral or contralateral limb of the GORE® EXCLUDER® AAA Endoprosthesis. An ePTFE / FEP sleeve is used to constrain the endoprosthesis on the leading end of the delivery catheter (Figure 6B). Deployment of the Iliac Extender initiates from the leading (aortic) end and proceeds toward the trailing (iliac) end of the delivery catheter. Following deployment, the ePTFE / FEP sleeve remains *in situ* between the endoprosthesis and the vessel wall.

Contralateral Leg Endoprosthesis Used as Iliac Extender Endoprosthesis (Distal Iliac Diameters of 16, 18, 20, 23, and 27 mm)

Only Contralateral Leg Endoprostheses described in Figure 2B may also be utilized as Iliac Extenders as depicted in Figure 7. A minimum overlap of 3 cm between the original Ipsilateral Leg (Figure 7A) or Contralateral Leg (Figure 7B and 7C) is required. A radiopaque marker located 3 cm below the proximal end is provided to confirm this overlap. This overlap should be achieved prior to the beginning of the distal taper zone of the 18, 20, 23, and 27 mm Contralateral Leg. Further, the distal end including the taper zone should not be deployed inside the previously deployed Ipsilateral Leg or Contralateral Leg of the GORE® EXCLUDER® AAA Endoprosthesis if the distal end diameters are 12 or 14.5 mm.

The 16, 18, and 20 mm Contralateral Legs have a radiopaque marker located 4 cm above the distal end which defines the recommended minimal extension required of the Contralateral Leg as an Iliac Extender. The 23 and 27 mm Contralateral Legs have a radiopaque marker located 5 and 6 cm, respectively,

above the distal end which defines the recommended minimal extension required of the Contralateral Leg as an Iliac Extender. * However, when the Contralateral Leg and Iliac Extender diameters are identical, the taper zone can be deployed inside the previously deployed Contralateral Leg (Figure 7B)

- * **Note: The 27 mm x 10 cm Contralateral Leg Endoprosthesis has only three (3) radiopaque markers; One (1) located at each end and one (1) located 3 cm below the proximal end.**

INDICATIONS FOR USE

Trunk-Ipsilateral Leg Endoprosthesis and Contralateral Leg Endoprosthesis Components

The GORE® EXCLUDER® AAA Endoprosthesis is intended to exclude the aneurysm from the blood circulation in patients diagnosed with infrarenal abdominal aortic aneurysm (AAA) disease and who have appropriate anatomy as described below:

- Adequate iliac / femoral access
- Infrarenal aortic neck treatment diameter range of 19 – 32 mm and a minimum aortic neck length of 15 mm
- Proximal aortic neck angulation $\leq 60^\circ$
- Iliac artery treatment diameter range of 8 – 25 mm and iliac distal vessel seal zone length of at least 10 mm

Aortic Extender Endoprosthesis and Iliac Extender Endoprosthesis Components

The Aortic and Iliac Extender Endoprostheses are intended to be used after deployment of the GORE® EXCLUDER® AAA Endoprosthesis. These extensions are intended to be used when additional length and / or sealing for aneurysmal exclusion is desired.

CONTRAINDICATIONS

The GORE® EXCLUDER® AAA Endoprosthesis is contraindicated in:

- Patients with known sensitivities or allergies to the device materials
- Patients with a systemic infection who may be at increased risk of endovascular graft infection.

WARNINGS AND PRECAUTIONS

General

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious surgical consequences or injury to the patient.
- The GORE® Medical Device is designed for single use only; do not reuse device. Gore does not have data regarding reuse of this device. Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination. Reuse may result in infection, serious injury, or patient death.
- The long-term performance of stent-grafts has not been established.** All patients should be advised this treatment modality requires long-term, regular follow-up to assess patients' health status and stent-graft performance. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should receive enhanced follow-up (See IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP).
- The GORE® EXCLUDER® AAA Endoprosthesis should only be used by physicians experienced in vascular interventional techniques, and who have successfully completed the appropriate physician training program.
- The GORE® EXCLUDER® AAA Endoprosthesis is not recommended in patients unable to undergo, or who will not be compliant with the necessary pre and post-operative imaging and follow-up described in IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP.
- Intervention or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms and / or persistent endoleak. An increase in aneurysm size and / or persistent endoleak may lead to aneurysm rupture.
- Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

Patient Selection, Treatment, and Follow-Up

- The safety and effectiveness of the GORE® EXCLUDER® AAA Endoprosthesis have not been evaluated in the following patient populations:
 - traumatic aortic injury
 - leaking: pending rupture or ruptured aneurysms
 - mycotic aneurysms
 - pseudoaneurysms resulting from previous graft placement
 - revision of previously placed stent grafts
 - genetic connective tissue disease (e.g., Marfans or Ehlers-Danlos Syndromes)
 - concomitant thoracic aortic or thoracoabdominal aneurysms
 - inflammatory aneurysms
 - patients with active systemic infections
 - pregnant or nursing females
 - morbidly obese patients
 - patients less than 21 years old
 - patients with less than 15 mm in length or $> 60^\circ$ angulation of the proximal aortic neck

- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and / or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 12 Fr - 18 Fr vascular introducer sheath.
- Key anatomic elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation, short proximal aortic neck and significant thrombus and / or calcium at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. The US clinical studies quantify significant thrombus as thrombus ≥ 2 mm in thickness and / or $\geq 25\%$ of the vessel circumference in the intended seal zone of the aortic neck. Irregular calcium and / or plaque may compromise the fixation and sealing of the implantation sites.
- The GORE® EXCLUDER® AAA Endoprosthesis is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and post-operative follow-up imaging.
- The GORE® EXCLUDER® AAA Endoprosthesis is not recommended in patients exceeding weight and / or size limits which compromise or prevent the necessary imaging requirements.
- The GORE® EXCLUDER® AAA Endoprosthesis is not recommended in patients with known sensitivities or allergies to ePTFE, FEP, nickel, or titanium.

Implant Procedure

- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not attempt to advance the Aortic Extender through the 12 Fr introducer sheath. The Aortic Extender is designed for a 16 Fr, 17 Fr or 18 Fr sheath.
- Do not rotate the Trunk, Contralateral Leg, Iliac Extender, or Aortic Extender delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not rotate the Trunk delivery catheter beyond 360°. Delivery system damage and / or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not rotate the Contralateral Leg, Iliac Extender, or Aortic Extender delivery catheter during delivery, positioning or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- **Do not attempt to reposition the endoprosthesis after complete deployment of the device. Vessel damage or device misplacement may result.**
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- **Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patients harms, see ADVERSE EVENTS.**
- Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
- Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.
- While using 16, 18, 20, 23, or 27 mm Contralateral Legs as an Iliac Extender, ensure that the distal end including the taper zone will not be deployed inside the previously deployed Ipsilateral Leg or Contralateral Leg of the GORE® EXCLUDER® AAA Endoprosthesis. However, when the Contralateral Leg and Iliac Extender diameters are identical, the taper zone can be deployed inside the previously deployed Contralateral Leg (Figure 7B).
- While using 16, 18, 20, 23, or 27 mm Contralateral Legs as an Iliac Extender, the 3 cm mandatory overlap must be achieved prior to the beginning of the distal taper zone of the 18, 20, 23, and 27 mm Contralateral Leg. Inadequate sealing may lead to endoleak.

MRI Information MR Conditional

Non-clinical testing has demonstrated that the GORE® EXCLUDER® AAA Endoprosthesis is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Highest spatial gradient magnetic field of 3000-Gauss/cm or less
- Whole body-averaged specific absorption rate (SAR) of 2.0-W/kg in the Normal Operating Mode (i.e., the mode of operation of the MR equipment in which none of the outputs have a value that cause physiological stress to patients) for 15 minutes of scanning (i.e., per pulse sequence).

MR Related Heating:

In non-clinical testing, the GORE® EXCLUDER® AAA Endoprosthesis produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numeris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and during MRI performed for 15-min of scanning 3-Tesla/128-MHz (Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

Table 1. MRI Information

	1.5 Tesla	3 Tesla
MR system reported, whole body averaged SAR	2.9 W/kg	2.9 W/kg
Calorimetry measured values, whole body averaged SAR	2.1 W/kg	2.7 W/kg
Highest recorded temperature change	2.2°C	2.7°C
Calculated Temperature Change Normal Operating Mode SAR = 2.0 W/kg	1.5°C	1.9°C
Calculated Temperature Change First Level Control Mode SAR = 4.0 W/kg	3.0°C	3.7°C

Image Artifact:

For each vascular device and assembly, the artifacts that appeared on the MR images were shown as localized signal voids (i.e., signal loss) that were minor in size relative to the size and shape of these implants. The gradient echo pulse sequence produced larger artifacts than the T1 – weighted, spin echo pulse sequence for the GORE® EXCLUDER® AAA Endoprosthesis. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10 mm relative to the size and shape of this implant. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the GORE® EXCLUDER® AAA Endoprosthesis. Therefore, it may be necessary to optimize the MR imaging parameters to compensate for the presence of this implant.

CLINICAL STUDY OUTCOMES* (98-03 AND 99-04)

Two US clinical studies were conducted to evaluate the safety and efficacy of the GORE® EXCLUDER® AAA Endoprosthesis. The first, referred to as 98-03, evaluated the original design of the GORE® EXCLUDER® AAA Endoprosthesis. This study enrolled 235 test subjects and 99 control subjects; 49 test subjects were enrolled in Continued Access. The second US clinical study, referred to as 99-04, evaluated a modified version of the GORE® EXCLUDER® AAA Endoprosthesis. These modifications were minor design changes that did not affect the risk benefit profile or clinical performance of the device. This study enrolled 193 test subjects; 88 test subjects were enrolled in Continued Access. Collectively, 565 test subjects were treated with the GORE® EXCLUDER® AAA Endoprosthesis in these studies and are described herein as the "Combined IDE Cohort". Safety outcomes of the Combined IDE Cohort are compared to the 99 control subjects treated with open surgery in the 98-03 study and are described herein as "98-03 Controls."

* Please note that the Contralateral Leg components with distal iliac diameters of 16, 18, 20, 23, and 27 mm were not evaluated under any clinical studies.

Baseline Characteristics

Overall, baseline characteristics were comparable in the Combined IDE Cohort vs. 98-03 Controls (Tables 2 and 3). The Combined IDE Cohort was an average of three years older and had a higher proportion of males vs. the 98-03 Controls. Smoking history, coronary artery disease, and chronic obstructive pulmonary disease were the most common comorbidities in each treatment group. Subjects with aneurysm diameter between 50-59 mm were most commonly treated, regardless of treatment group.

Table 2. Subject Demographics and Medical History

	98-03 Control	Combined IDE Cohort
Number of Subjects Enrolled	99	565
Age at Treatment (Years)		
n	99	565
Mean (Std Dev)	70.6 (8.5)	73.9 (7.9)
Median	71.8	74.3
Range	(51.8, 87.6)	(48.9, 93.0)
Gender		
Male	73 (73.7%)	456 (80.7%)
Female	26 (26.3%)	109 (19.3%)
Ethnic Background		
Caucasian	99 (100.0%)	553 (97.9%)
Black	0 (0.0%)	5 (0.9%)
Hispanic	0 (0.0%)	5 (0.9%)
Asian	0 (0.0%)	2 (0.4%)
NYHA Functional Class		
1	64 (64.6%)	289 (51.2%)
2	24 (24.2%)	201 (35.6%)
3	11 (11.1%)	48 (8.5%)
Missing	0 (0.0%)	27 (4.8%)
ASA Functional Class		
1	2 (2.0%)	8 (1.4%)
2	21 (21.2%)	111 (19.6%)
3	60 (60.6%)	373 (66.0%)
4	15 (15.2%)	73 (12.9%)
Missing	1 (1.0%)	0 (0.0%)
SVS Summary Risk Score		
n	99	565
Mean (Std Dev)	4.6 (2.7)	4.9 (2.7)
Median	4.0	5.0
Range	(0.0, 12.0)	(0.0, 13.3)
Aneurysm Symptomatic	15 (15.2%)	38 (6.7%)
Arrhythmia	21 (21.2%)	135 (23.9%)
Bleeding Disorder	1 (1.0%)	16 (2.8%)
Cancer	19 (19.2%)	150 (26.5%)
Congestive Heart Failure	8 (8.1%)	58 (10.3%)
COPD	25 (25.3%)	167 (29.6%)
Coronary Artery Disease	53 (53.5%)	344 (60.9%)
Erectile Dysfunction (Males Only)	10 (13.7%)	78 (17.1%)
Family History of AAA	9 (9.1%)	42 (7.4%)
Hepatic Dysfunction	1 (1.0%)	9 (1.6%)
Inflammatory AAA	1 (1.0%)	7 (1.2%)
Long-Term Use of Steroids	1 (1.0%)	23 (4.1%)
Other Concomitant Aneurysms	13 (13.1%)	47 (8.3%)
Peripheral Arterial Occlusive Disease	14 (14.1%)	108 (19.1%)
Paraplegia	0 (0.0%)	2 (0.4%)
Prior Vascular Intervention	10 (10.1%)	63 (11.2%)
Renal Dialysis	0 (0.0%)	2 (0.4%)
Smoking History	85 (85.9%)	484 (85.7%)
Stroke	10 (10.1%)	59 (10.4%)
Thromboembolic Event	4 (4.0%)	34 (6.0%)
Valvular Heart Disease	7 (7.1%)	45 (8.0%)

Table 3. Distribution of Pre-Treatment Aneurysm Diameter

	98-03 Control	Combined IDE Cohort
Subjects Enrolled	99	565
Diameter Range		
<30 mm	0 (0.0%)	0 (0.0%)
30-39 mm	0 (0.0%)	0 (0.0%)
40-49 mm	15 (15.2%)	122 (21.6%)
50-59 mm	46 (46.5%)	298 (52.7%)
60-69 mm	21 (21.2%)	99 (17.5%)
70-79 mm	10 (10.1%)	31 (5.5%)
80-89 mm	5 (5.1%)	10 (1.8%)
≥90 mm	1 (1.0%)	5 (0.9%)
Missing ¹	1 (1.0%)	0 (0.0%)

1 Missing indicates no pre-treatment diameter measurement was available.

Subject Compliance and Disposition

Subjects in each group were followed for five years post-treatment. Of 565 subjects in the Combined IDE Cohort, 255 returned for the five-year follow-up visit, 149 died, and 177 were discontinued from the study (Table 4). The number of "Subjects with Visit in Window" (compliance) is an exclusive category from the "Death," "Discontinued," and "Not Due for Next Follow-Up" (disposition) categories. A subject may have had a visit within a follow-up window and died or discontinued from the study later in the same window. Therefore, adding these entries will not necessarily equal the number of subjects "Eligible for Follow-Up."

Follow-up visit and imaging compliance remained high throughout the five-year follow-up period. Of 99 subjects in the 98-03 Control group, 14 died, and 46 were discontinued from the five-year follow-up visit. Abdominal x-rays were not required for the 98-03 Controls.

Table 4. Subject Compliance and Disposition by Study Interval: Combined IDE Cohort

		Follow-Up Compliance				Events Prior to Next Interval		
Study Period	Eligible for Follow-Up	Subjects with Visit in Window ¹	With CT ¹	With XRay ¹	With ABI ¹	Death ¹	Discontinued ¹	Not Due for Next Follow-Up ¹
Procedure	565	—	—	—	—	2 (0.4%)	1 (0.2%)	0 (0.0%)
Post-Procedure	562	—	—	—	—	1 (0.2%)	1 (0.2%)	0 (0.0%)
1 Month	560	545 (97.3%)	522 (93.2%)	103 (18.4%)	461 (82.3%)	6 (1.1%)	1 (0.2%)	0 (0.0%)
3 Months	553	94 (17.0%)	85 (15.4%)	11 (2.0%)	37 (6.7%)	6 (1.1%)	0 (0.0%)	0 (0.0%)
6 Months	547	507 (92.7%)	486 (88.8%)	411 (75.1%)	422 (77.1%)	7 (1.3%)	7 (1.3%)	0 (0.0%)
12 Months	533	501 (94.0%)	484 (90.8%)	412 (77.3%)	414 (77.7%)	37 (6.9%)	13 (2.4%)	0 (0.0%)
24 Months	483	441 (91.3%)	412 (85.3%)	360 (74.5%)	350 (72.5%)	25 (5.2%)	31 (6.4%)	0 (0.0%)
36 Months	427	340 (79.6%)	319 (74.7%)	274 (64.2%)	265 (62.6%)	28 (6.6%)	38 (8.9%)	0 (0.0%)
48 Months	361	306 (84.8%)	285 (78.3%)	240 (66.5%)	219 (60.7%)	26 (7.2%)	44 (12.2%)	0 (0.0%)
60 Months	291	255 (87.6%)	236 (81.1%)	185 (63.6%)	178 (61.2%)	11 (3.8%)	41 (14.1%)	—

1 Denominators are based on the number of subjects eligible for follow-up
Study period definitions: Procedure (0-0 days), Post-Procedure (1-14 days), 1 Month (15-60 days), 3 Months (61-120 days), 6 Months (121-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days).

Table 5. Subject Compliance and Disposition by Study Interval: 98-03 Controls

		Follow-Up Compliance				Events Prior to Next Interval		
Study Period	Eligible for Follow-Up	Subjects with Visit in Window ¹	With CT ¹	With XRay ¹	With ABI ¹	Death ¹	Discontinued ¹	Not Due for Next Follow-Up ¹
Procedure	99	—	—	—	—	0 (0.0%)	0 (0.0%)	0 (0.0%)
Post-Procedure	99	—	—	—	—	0 (0.0%)	1 (1.0%)	0 (0.0%)
1 Month	98	89 (90.8%)	0 (0.0%)	0 (0.0%)	67 (68.4%)	1 (1.0%)	1 (1.0%)	0 (0.0%)
3 Months	96	11 (11.5%)	0 (0.0%)	0 (0.0%)	5 (5.2%)	1 (1.0%)	1 (1.0%)	0 (0.0%)
6 Months	94	81 (86.2%)	4 (4.3%)	0 (0.0%)	66 (70.2%)	2 (2.1%)	5 (5.3%)	0 (0.0%)
12 Months	87	83 (95.4%)	74 (85.1%)	4 (4.6%)	67 (77.0%)	1 (1.1%)	5 (5.7%)	0 (0.0%)
24 Months	81	70 (86.4%)	67 (82.7%)	6 (7.4%)	56 (69.1%)	2 (2.5%)	9 (11.1%)	0 (0.0%)
36 Months	70	52 (74.3%)	46 (65.7%)	4 (5.7%)	46 (65.7%)	3 (4.3%)	9 (12.9%)	0 (0.0%)
48 Months	58	47 (81.0%)	40 (69.0%)	5 (8.6%)	37 (63.8%)	3 (5.2%)	9 (15.5%)	0 (0.0%)
60 Months	46	40 (87.0%)	34 (73.9%)	7 (15.2%)	33 (71.7%)	1 (2.2%)	6 (13.0%)	—

1 Denominators are based on the number of subjects eligible for follow-up
Study period definitions: Procedure (0-0 days), Post-Procedure (1-14 days), 1 Month (15-60 days), 3 Months (61-120 days), 6 Months (121-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days).

Procedure and Recovery

Endovascular access was successful in 99.6% of cases in the Combined IDE Cohort (Table 6). The Combined IDE Cohort experienced less procedural blood loss, required fewer transfusions, had a shorter treatment time, and had shorter convalescence compared to 98-03 Controls (Table 7).

Table 6. Device Use at Initial Procedure

	Combined IDE Cohort
Number of Subjects Enrolled	565
Number of Subjects with Devices Implanted	563
Subjects with Trunks Implanted	563 (100.0%)
Subjects with Contralateral Legs Implanted	556 (98.8%)
Subjects with Aortic Extender(s) Implanted	66 (11.7%)
Subjects with Iliac Extender(s) Implanted	152 (27.0%)

Table 7. Procedure and Recovery

	98-03 Control	Combined IDE Cohort
Number of Subjects Enrolled	99	565
Blood Loss During Procedure (mL)		
n	98	557
Mean (Std Dev)	1589 (1226)	324 (331)
Median	1175	200
Range	(100, 7000)	(0, 3000)
Procedure Transfusion	88 (88.9%)	77 (13.6%)
Procedure Time (Minutes)		
n	95	563
Mean (Std Dev)	196 (78)	143 (57)
Median	175	132
Range	(67, 420)	(45, 509)
ICU Stay	86 (86.9%)	114 (20.2%)
Length of Hospital Stay (Days)		
n	99	565
Mean (Std Dev)	10 (14)	2 (1)
Median	7	2
Range	(3, 114)	(0, 15)
Time to First Oral Intake (Days)		
n	98	561
Mean (Std Dev)	2.7 (1.7)	0.4 (0.6)
Median	3.0	0.0
Range	(0.0, 10.0)	(0.0, 7.0)
Time to Ambulation (Days)		
n	98	559
Mean (Std Dev)	2.6 (2.3)	1.0 (0.7)
Median	2.0	1.0
Range	(0.0, 18.0)	(0.0, 8.0)
Time to Return to Normal Daily Activities (Days)		
n	90	530
Mean (Std Dev)	101.1 (137.4)	40.7 (60.9)
Median	41.0	26.0
Range	(1.0, 828.0)	(0.0, 527.0)

Safety Outcomes

Adverse events were classified as major or minor¹. An adverse event classified as major required therapy and short hospitalization (24-48 hours), major therapy, unplanned increase in level of care, prolonged hospitalization (> 48 hours), or resulted in permanent adverse sequelae, or death. Freedom from a major adverse event was statistically greater ($p < 0.001$) in the Combined IDE Cohort (39%) vs. 98-03 Controls (15%) through five years post-treatment (Figure 10, Table 8). Cardiac and pulmonary complications were the most-commonly reported major adverse events in the Combined IDE Cohort (Table 9); blood loss (> 1 L) and cardiac events were the most commonly reported major adverse events in 98-03 Controls (Table 10). Cumulative major adverse events per subject were statistically lower in the Combined IDE Cohort (1.6 / subject) vs. 98-03 Controls (2.8 / subject) through five years post-treatment (Figure 11, Table 11).

¹ Sacks D, Marinelli DL, Martin LG, et al. Reporting standards for clinical evaluation of new peripheral arterial revascularization devices. *Journal of Vascular and Interventional Radiology* 1997; 8: 137-149.

Figure 10. Freedom From Major Adverse Event Through Five Years Post-Treatment


Table 8. Freedom From Major Adverse Event Through Five Years Post-Treatment

	Months After Procedure							
	Day 0	1 month	6 months	1 year	2 years	3 years	4 years	5 years
98-03 Control								
Subjects at Risk	99	32	30	26	20	16	12	7
% Free From Major Adverse Event	47%	31%	30%	26%	24%	23%	21%	15%
Combined IDE								
Subjects at Risk	565	477	428	371	302	243	190	110
% Free From Major Adverse Event	94%	84%	76%	67%	57%	50%	44%	39%

Table 9. Major Adverse Events Through Five Years Post-Treatment by Visit Interval: Combined IDE Cohort

	1 Month	6 Months	12 Months	24 Months	36 Months	48 Months	60 Months	Overall
Subjects Available at Beginning of Interval ¹	565	553	533	483	427	361	291	565
Subjects with Any Major Safety Adverse Event ²	105 (18.6%)	72 (13.0%)	91 (17.1%)	81 (16.8%)	64 (15.0%)	53 (14.7%)	22 (7.6%)	326 (57.7%)
Bleeding Complications	27 (4.8%)	2 (0.4%)	1 (0.2%)	3 (0.6%)	1 (0.2%)	5 (1.4%)	0 (0.0%)	38 (6.7%)
Pulmonary Complications	13 (2.3%)	13 (2.4%)	17 (3.2%)	10 (2.1%)	24 (5.6%)	10 (2.8%)	7 (2.4%)	81 (14.3%)
Cardiac Complications	22 (3.9%)	21 (3.8%)	24 (4.5%)	23 (4.8%)	23 (5.4%)	21 (5.8%)	10 (3.4%)	124 (21.9%)
Renal Function Complications	9 (1.6%)	6 (1.1%)	5 (0.9%)	6 (1.2%)	6 (1.4%)	1 (0.3%)	3 (1.0%)	30 (5.3%)
Wound Complications	22 (3.9%)	9 (1.6%)	3 (0.6%)	4 (0.8%)	1 (0.2%)	1 (0.3%)	0 (0.0%)	37 (6.5%)
Bowel Complications	12 (2.1%)	5 (0.9%)	10 (1.9%)	10 (2.1%)	4 (0.9%)	2 (0.6%)	2 (0.7%)	43 (7.6%)
Vascular Complications	16 (2.8%)	10 (1.8%)	9 (1.7%)	8 (1.7%)	5 (1.2%)	2 (0.6%)	2 (0.7%)	44 (7.8%)
Neurologic Complications	7 (1.2%)	5 (0.9%)	12 (2.3%)	9 (1.9%)	8 (1.9%)	3 (0.8%)	0 (0.0%)	38 (6.7%)
Other Complications	8 (1.4%)	12 (2.2%)	26 (4.9%)	19 (3.9%)	8 (1.9%)	11 (3.0%)	2 (0.7%)	77 (13.6%)
Genitourinary	6 (1.1%)	5 (0.9%)	0 (0.0%)	1 (0.2%)	0 (0.0%)	3 (0.8%)	1 (0.3%)	16 (2.8%)
Sepsis	1 (0.2%)	2 (0.4%)	0 (0.0%)	2 (0.4%)	0 (0.0%)	2 (0.6%)	0 (0.0%)	7 (1.2%)
Neoplasm	3 (0.5%)	8 (1.4%)	16 (3.0%)	10 (2.1%)	11 (2.6%)	14 (3.9%)	1 (0.3%)	61 (10.8%)
Death of Unknown Cause	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	1 (0.3%)	0 (0.0%)	1 (0.2%)
Other	8 (1.4%)	12 (2.2%)	26 (4.9%)	19 (3.9%)	8 (1.9%)	11 (3.0%)	2 (0.7%)	77 (13.6%)

¹ Number of subjects with study day of last contact $\geq$ lower limit of specified time window. This is the denominator for all cells.
² Entries represent number of subjects with the event (percentage). Time frames for each interval are as follows: 1 Month (0-60 days), 6 Months (61-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days), Overall (0-2006 days).

Table 10. Major Adverse Events Through Five Years Post-Treatment by Visit Interval: 98-03 Controls

	1 Month	6 Months	12 Months	24 Months	36 Months	48 Months	60 Months	Overall
Subjects Available at Beginning of Interval ¹	99	96	87	81	70	58	46	99
Subjects with Any Major Safety Adverse Event ²	68 (68.7%)	15 (15.6%)	10 (11.5%)	11 (13.6%)	12 (17.1%)	13 (22.4%)	3 (6.5%)	80 (80.8%)
Bleeding Complications	45 (45.5%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	1 (1.7%)	0 (0.0%)	45 (45.5%)
Pulmonary Complications	12 (12.1%)	3 (3.1%)	1 (1.1%)	2 (2.5%)	1 (1.4%)	4 (6.9%)	0 (0.0%)	19 (19.2%)
Cardiac Complications	16 (16.2%)	7 (7.3%)	7 (8.0%)	6 (7.4%)	6 (8.6%)	6 (10.3%)	1 (2.2%)	36 (36.4%)
Renal Function Complications	3 (3.0%)	0 (0.0%)	1 (1.1%)	0 (0.0%)	1 (1.4%)	3 (5.2%)	0 (0.0%)	7 (7.1%)
Wound Complications	4 (4.0%)	3 (3.1%)	0 (0.0%)	1 (1.2%)	2 (2.9%)	1 (1.7%)	0 (0.0%)	10 (10.1%)
Bowel Complications	17 (17.2%)	2 (2.1%)	1 (1.1%)	1 (1.2%)	2 (2.9%)	1 (1.7%)	1 (2.2%)	22 (22.2%)
Vascular Complications	7 (7.1%)	1 (1.0%)	3 (3.4%)	1 (1.2%)	0 (0.0%)	0 (0.0%)	1 (2.2%)	12 (12.1%)
Neurologic Complications	3 (3.0%)	3 (3.1%)	2 (2.3%)	0 (0.0%)	0 (0.0%)	2 (3.4%)	0 (0.0%)	9 (9.1%)
Other Complications	1 (1.0%)	2 (2.1%)	2 (2.3%)	2 (2.5%)	2 (2.9%)	1 (1.7%)	0 (0.0%)	9 (9.1%)
Genitourinary	1 (1.0%)	0 (0.0%)	2 (2.3%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	3 (3.0%)
Sepsis	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Neoplasm	0 (0.0%)	1 (1.0%)	2 (2.3%)	0 (0.0%)	3 (4.3%)	0 (0.0%)	0 (0.0%)	6 (6.1%)
Death of Unknown Cause	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Other	1 (1.0%)	2 (2.1%)	2 (2.3%)	2 (2.5%)	2 (2.9%)	1 (1.7%)	0 (0.0%)	9 (9.1%)

¹ Number of subjects with study day of last contact $\geq$ lower limit of specified time window. This is the denominator for all cells.
² Entries represent number of subjects with the event (percentage). Time frames for each interval are as follows: 1 Month (0-60 days), 6 Months (61-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days), Overall (0-2006 days).

Figure 11. Cumulative Major Adverse Events Through Five Years Post-Treatment


Table 11. Cumulative Major Adverse Events Through Five Years Post-Treatment

	Months After Procedure							
	Day 0	1 month	6 months	1 year	2 years	3 years	4 years	5 years
98-03 Control								
Subjects at Risk	99	97	91	86	79	64	52	35
Cumulative Major AE	0.60	1.17	1.46	1.64	1.88	2.06	2.55	2.83
Combined IDE								
Subjects at Risk	565	556	545	512	457	402	332	196
Cumulative Major AE	0.08	0.24	0.43	0.63	0.91	1.18	1.42	1.64

Mortality

Freedom from aneurysm-related death was similar in the Combined IDE Cohort (98%) vs. 98-03 Controls (98%) through five years post-treatment (Figure 11, Table 12). Survival was lower ($p < 0.05$) in the Combined IDE Cohort (70%) vs. 98-03 Controls (81%) through five years post-treatment (Figure 12, Table 13). However, the Combined IDE Cohort was an average of three years older vs. the 98-03 Controls (73.9 vs. 70.6 years), and after controlling for age, no difference between groups was observed in survival through five years post-treatment.

Figure 12. Freedom From Aneurysm-Related Death Through Five Years Post-Treatment


Table 12. Freedom From Aneurysm-Related Death Through Five Years Post-Treatment

	Months After Procedure							
	Day 0	1 month	6 months	1 year	2 years	3 years	4 years	5 years
98-03 Control								
Subjects at Risk	99	97	91	86	79	64	52	35
% Free of Aneurysm-Related Death	100%	100%	98%	98%	98%	98%	98%	98%
Combined IDE								
Subjects at Risk	565	556	545	512	457	401	332	196
% Free of Aneurysm-Related Death	100%	99%	99%	99%	99%	98%	98%	98%

Figure 13. Survival Through Five Years Post-Treatment


Table 13. Survival Through Five Years Post-Treatment

	Months After Procedure							
	Day 0	1 month	6 months	1 year	2 years	3 years	4 years	5 years
98-03 Control								
Subjects at Risk	99	97	91	86	79	64	52	35
% Survival	100%	100%	97%	95%	94%	92%	85%	81%
Combined IDE								
Subjects at Risk	565	556	545	512	457	401	332	196
% Survival	100%	99%	97%	93%	87%	82%	76%	70%

Device Efficacy Outcomes

A major device event was observed in 14.3% of the Combined IDE Cohort through five years post-treatment (Table 14); endoleak and aneurysm diameter increase were the most commonly reported major device events. Aneurysm rupture was noted in one subject during the five-year follow-up period. The rupture was non-hemorrhagic and the subject was successfully converted to open surgical repair during an elective procedure. Endoleak incidence ranged from 14 – 27% at each follow-up period through five years post-treatment (Table 15). New onset of endoleak became less frequent over time. Type II endoleak was most common. Incidence of Type I or III endoleak ranged from 0 – 2.6% at all follow-up visits through five years post-treatment. The incidence of aneurysm diameter increase was 34% at the five-year follow-up visit in the Combined IDE Cohort (Table 16). However, most (78%, 62 / 79) of these subjects were free from a detectable endoleak. Reintervention for endoleak or aneurysm diameter increase was performed in 1.9 – 5.3% of subjects at each follow-up visit (Table 17). Common treatments for aneurysm enlargement were explant ($n = 12$), endovascular stent graft ($n = 10$), and embolization ($n = 8$); endoleak treatments were embolization ($n = 61$), stent graft ($n = 10$), and explant ($n = 3$). Treatments for these events are not additive as a single treatment may be used to treat aneurysm enlargement and endoleak.

Table 14. Major Device Events Through Five Years Post-Treatment by Visit Interval: Combined IDE Cohort

	1 Month	6 Months	12 Months	24 Months	36 Months	48 Months	60 Months	Overall
Subjects Available at Beginning of Interval ¹	565	553	533	483	427	361	291	565
Subjects with Any Major Efficacy Adverse Event ²	32 (5.7%)	16 (2.9%)	12 (2.3%)	18 (3.7%)	8 (1.9%)	9 (2.5%)	6 (2.1%)	81 (14.3%)
Branch Vessel Occlusion	4 (0.7%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	4 (0.7%)
Endoleak	24 (4.2%)	15 (2.7%)	7 (1.3%)	10 (2.1%)	3 (0.7%)	4 (1.1%)	1 (0.3%)	55 (9.7%)
Access Failure	2 (0.4%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	2 (0.4%)
Deployment Failure	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Lumen Obstruction	1 (0.2%)	1 (0.2%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	2 (0.4%)
Prosthesis Material Failure	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Aneurysm Rupture	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	1 (0.3%)	1 (0.2%)
Other Device-Related Event	1 (0.2%)	0 (0.0%)	0 (0.0%)	2 (0.4%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	3 (0.5%)
Extrusion / Erosion	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Prosthesis Migration	1 (0.2%)	0 (0.0%)	2 (0.4%)	0 (0.0%)	1 (0.2%)	0 (0.0%)	0 (0.0%)	4 (0.7%)
Prosthesis Realignment	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Aneurysm Increase with Reintervention	2 (0.4%)	1 (0.2%)	3 (0.6%)	8 (1.7%)	4 (0.9%)	7 (1.9%)	5 (1.7%)	29 (5.1%)

¹ Number of subjects with study day of last contact $\geq$ lower limit of specified time window. This is the denominator for all cells.
² Entries represent number of subjects with the event (percentage). Time frames for each interval are as follows: 1 Month (0-60 days), 6 Months (61-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days), Overall (0-2006 days).

Table 15. Endoleaks Through Five Years Post-Treatment by Visit Interval: Combined IDE Cohort

	Treatment	1 Month	6 Months	12 Months	24 Months	36 Months	48 Months	60 Months
Subjects Available at Beginning of Interval	565	562	553	533	483	427	361	291
Subjects with Endoleak Evaluation or Ongoing Endoleak	565	538	509	489	425	335	297	244
Subjects With One or More Endoleak Adverse Events Ongoing in Window	79 (14.0%)	147 (27.3%)	130 (25.5%)	110 (22.5%)	93 (21.9%)	74 (22.1%)	61 (20.5%)	40 (16.4%)
New	79 (14.0%)	83 (15.4%)	46 (9.0%)	26 (5.3%)	33 (7.8%)	11 (3.3%)	13 (4.4%)	6 (2.5%)
Ongoing	—	77 (14.3%)	90 (17.7%)	91 (18.6%)	72 (16.9%)	65 (19.4%)	51 (17.2%)	35 (14.3%)
Type I	11 (1.9%)	14 (2.6%)	6 (1.2%)	2 (0.4%)	4 (0.9%)	4 (1.2%)	3 (1.0%)	2 (0.8%)
New	11 (1.9%)	5 (0.9%)	1 (0.2%)	2 (0.4%)	2 (0.5%)	1 (0.3%)	0 (0.0%)	1 (0.4%)
Ongoing	—	9 (1.7%)	5 (1.0%)	1 (0.2%)	2 (0.5%)	3 (0.9%)	3 (1.0%)	1 (0.4%)
Type II	66 (11.7%)	123 (22.9%)	106 (20.8%)	99 (20.2%)	85 (20.0%)	65 (19.4%)	50 (16.8%)	30 (12.3%)
New	66 (11.7%)	61 (11.3%)	37 (7.3%)	24 (4.9%)	25 (5.9%)	8 (2.4%)	7 (2.4%)	4 (1.6%)
Ongoing	—	66 (12.3%)	72 (14.1%)	79 (16.2%)	66 (15.5%)	58 (17.3%)	44 (14.8%)	27 (11.1%)
Type III	1 (0.2%)	3 (0.6%)	6 (1.2%)	4 (0.8%)	2 (0.5%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
New	1 (0.2%)	2 (0.4%)	3 (0.6%)	2 (0.4%)	1 (0.2%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Ongoing	—	1 (0.2%)	3 (0.6%)	2 (0.4%)	1 (0.2%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Other	3 (0.5%)	8 (1.5%)	6 (1.2%)	4 (0.8%)	0 (0.0%)	1 (0.3%)	3 (1.0%)	2 (0.8%)
New	3 (0.5%)	6 (1.1%)	1 (0.2%)	1 (0.2%)	0 (0.0%)	1 (0.3%)	2 (0.7%)	0 (0.0%)
Ongoing	—	2 (0.4%)	5 (1.0%)	3 (0.6%)	0 (0.0%)	0 (0.0%)	1 (0.3%)	2 (0.8%)
Indeterminate	3 (0.5%)	14 (2.6%)	15 (2.9%)	10 (2.0%)	13 (3.1%)	10 (3.0%)	11 (3.7%)	7 (2.9%)
New	3 (0.5%)	11 (2.0%)	5 (1.0%)	2 (0.4%)	7 (1.6%)	1 (0.3%)	4 (1.3%)	1 (0.4%)
Ongoing	—	3 (0.6%)	11 (2.2%)	8 (1.6%)	7 (1.6%)	9 (2.7%)	8 (2.7%)	6 (2.5%)
Subjects With No Endoleak Adverse Events Ongoing in Window	486 (86.0%)	391 (72.7%)	379 (74.5%)	379 (77.5%)	332 (78.1%)	261 (77.9%)	236 (79.5%)	204 (83.6%)

Note: Denominators are the number of subjects with endoleak evaluation or ongoing endoleak at each interval.
Time frames for each interval are as follows: Treatment (0-0 days), 1 Month (1-60 days), 6 Months (61-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days).

Table 16. Aneurysm Diameter Change Through Five Years Post-Treatment by Visit Interval: Combined IDE Cohort

	6 Months	12 Months	24 Months	36 Months	48 Months	60 Months
Subjects Available at Beginning of Interval ¹	553	533	483	427	361	291
Subjects with Aneurysm Evaluation ²	476 (86.1%)	458 (85.9%)	397 (82.2%)	311 (72.8%)	279 (77.3%)	231 (79.4%)
$\geq$ 5 mm Aneurysm Size Decrease ⁴	66 (13.9%)	100 (21.8%)	100 (25.2%)	79 (25.4%)	65 (23.3%)	53 (22.9%)
No Change in Aneurysm Diameter	401 (84.2%)	335 (73.1%)	251 (63.2%)	174 (55.9%)	134 (48.0%)	99 (42.9%)
$\geq$ 5 mm Aneurysm Size Increase ^{3,4}	9 (1.9%)	23 (5.0%)	46 (11.6%)	58 (18.6%)	80 (28.7%)	79 (34.2%)
Subjects with Aneurysm Increase and an Endoleak	4	13	20	16	22	17
Type I	0	1	1	0	0	1
Type II	3	12	19	15	18	13
Type III	0	0	1	0	0	0
Indeterminate	1	0	2	1	4	3

¹ The number eligible for evaluation is based on the number of subjects that enter the beginning of the follow-up interval alive and in the study.
² The denominator used to calculate percentages at each interval is based on the number of subjects eligible for evaluation.
³ The denominator used to calculate percentages at each interval is based on the number of subjects with aneurysm evaluation at each timepoint.
⁴ Increase is based on a change of greater than or equal to 5 mm from baseline.
Time frames for each interval are as follows: 6 Months (61-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days).

Table 17. Reinterventions Due to Aneurysm Diameter Increase or Endoleak: Combined IDE Cohort

	1 Month	6 Months	12 Months	24 Months	36 Months	48 Months	60 Months
Subjects Available at Beginning of Interval	565	553	533	483	427	361	291
Subjects with Endoleak or Aneurysm Increase AEs Requiring an Intervention	30 (5.3%)	19 (3.4%)	17 (3.2%)	19 (3.9%)	8 (1.9%)	9 (2.5%)	6 (2.1%)
Aneurysm Increase Interventions							
Embolization	1	0	2	3	1	1	0
Endovascular Stent Graft	0	0	0	2	3	3	2
Explant	2	0	1	3	1	3	2
Other	0	1	0	0	1	0	1
Endoleak Interventions							
Embolization	19	17	12	8	2	2	1
Endovascular Stent Graft	5	0	1	3	1	0	0
Explant	2	0	0	0	0	1	0
Other	5	3	1	1	0	0	0

Note: Time frames for each interval are as follows: 1 Month (0-60 days), 6 Months (61-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days).

Imaging Core Laboratory Outcomes

An independent imaging core laboratory assessed follow-up CTs and X-rays for device events through five years post-treatment in the Combined IDE Cohort (Table 18). Endoleak incidence ranged from 6.6 – 23.5% at each follow-up visit through five years post-treatment; most endoleaks were Type II. Aneurysm enlargement was observed in 36.1% of subjects at the five-year follow-up visit. Most subjects with aneurysm enlargement had no detectable endoleak.

Table 18. Imaging Core Laboratory Results Through Five Years Post-Treatment by Visit Interval: Combined IDE Cohort

	1 Month	6 Months	12 Months	24 Months	36 Months	48 Months	60 Months
Subjects Available at Beginning of Interval	565	553	533	483	427	361	291
Subjects with CT Imaging (sent to Core Lab)	443	409	408	339	249	211	170
Subjects with X-Ray Imaging (sent to Core Lab)	423	327	325	295	213	186	139
Trunk Migration	—	0 / 370 (0.0%)	0 / 382 (0.0%)	0 / 320 (0.0%)	0 / 226 (0.0%)	0 / 197 (0.0%)	0 / 155 (0.0%)
Component Migration	—	1 / 254 (0.4%)	2 / 301 (0.7%)	0 / 281 (0.0%)	0 / 206 (0.0%)	0 / 184 (0.0%)	0 / 135 (0.0%)
Fracture	1 / 352 (0.3%)	1 / 304 (0.3%)	1 / 301 (0.3%)	1 / 278 (0.4%)	0 / 209 (0.0%)	0 / 184 (0.0%)	0 / 133 (0.0%)
Narrowing of Flow Channel	8 / 427 (1.9%)	3 / 391 (0.8%)	4 / 389 (1.0%)	4 / 326 (1.2%)	3 / 238 (1.3%)	2 / 197 (1.0%)	1 / 155 (0.6%)
Endoleak	68 / 351 (19.4%)	76 / 324 (23.5%)	60 / 310 (19.4%)	41 / 271 (15.1%)	26 / 198 (13.1%)	14 / 179 (7.8%)	9 / 137 (6.6%)
Type I	7	9	4	3	0	0	0
Type II	41	45	40	31	18	10	7
Type III	1	2	5	0	1	1	0
Indeterminate	20	29	16	8	9	3	1
Missing	0	1	0	1	0	0	1
Aneurysm Size Change From Baseline							
$\geq$ 5 mm Aneurysm Size Decrease	—	44 / 372 (11.8%)	58 / 374 (15.5%)	65 / 314 (20.7%)	45 / 226 (19.9%)	47 / 193 (24.4%)	35 / 158 (22.2%)
No Change in Aneurysm Diameter	—	315 / 372 (84.7%)	292 / 374 (78.1%)	205 / 314 (65.3%)	119 / 226 (52.7%)	87 / 193 (45.1%)	66 / 158 (41.8%)
$\geq$ 5 mm Aneurysm Size Increase	—	13 / 372 (3.5%)	24 / 374 (6.4%)	44 / 314 (14.0%)	62 / 226 (27.4%)	59 / 193 (30.6%)	57 / 158 (36.1%)
Subjects with Aneurysm Increase and an Endoleak	—	4	8	12	15	10	5
Type I	—	1	0	1	0	0	0
Type II	—	2	5	7	10	7	4
Type III	—	0	0	0	1	1	0
Indeterminate	—	2	3	4	5	2	1
Missing	—	0	0	1	0	0	0

Note: Denominators for each variable are based on the number of evaluable images at each follow-up visit.
Time frames for each interval are as follows: 1 Month (0-60 days), 6 Months (61-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days).
Note: Migration and Aneurysm Size Change use one month as baseline, and therefore are not assessed at the one month time point.

CLINICAL STUDY OUTCOMES (04-04)

In response to the aneurysm enlargement rates associated with the original GORE® EXCLUDER® Device (shown in Table 16), Gore enhanced the design of the device by making changes to the graft material that decreased the overall permeability of the device. The low permeability GORE® EXCLUDER® AAA Endoprosthesis was approved by the FDA and released in June 2004. Gore has conducted a post-approval clinical study (04-04) of 139 subjects to evaluate the clinical performance of the low permeability GORE® EXCLUDER® AAA Endoprosthesis. These subjects are described herein as the "04-04 Test Group." Controls in this study are 120 randomly selected subjects treated with the GORE® EXCLUDER® AAA Endoprosthesis from the 98-03 clinical study.

Subjects in the 04-04 Test Group have completed a two-year follow-up (Table 19). Seven subjects died during the study. Causes of death included myocardial infarction (2), lung cancer, esophageal cancer, pneumonia (2), and sepsis / bowel perforation.

Aneurysm enlargement assessed by diameter and volume was less frequent in the 04-04 Test Group vs. Controls (Tables 20 and 21). Furthermore, aneurysm shrinkage was more frequent in the 04-04 Test Group vs. Controls. Outcomes of the 04-04 study suggest that the design enhancement of the low permeability GORE® EXCLUDER® AAA Endoprosthesis mitigate the risk of aneurysm enlargement without endoleak.

Table 19. Subject Compliance and Disposition by Study Interval

Study Period	Eligible for Follow-Up	Follow-Up Compliance				Events Prior to Next Interval		
		Subjects with Visit in Window ¹	Site Reported CT Scan Read ¹	Site Reported X-Ray Read ²	Corelab CT Scan Read ¹	Death ¹	Discontinued ¹	Not Due for Next Follow-Up ¹
Procedure	139	139 (100.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Post-Procedure	139	139 (100.0%)	0 (0.0%)	1 (0.7%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
1 Month	139	134 (96.4%)	130 (93.5%)	60 (43.2%)	130 (93.5%)	0 (0.0%)	4 (2.9%)	0 (0.0%)
6 Months	135	115 (82.2%)	111 (82.2%)	35 (25.9%)	111 (82.2%)	3 (2.2%)	3 (2.2%)	0 (0.0%)
12 Months	129	123 (95.3%)	121 (93.8%)	30 (23.3%)	121 (93.8%)	3 (2.3%)	8 (6.2%)	0 (0.0%)
24 Months	118	115 (97.5%)	110 (93.2%)	27 (22.9%)	110 (93.2%)	0 (0.0%)	2 (1.7%)	3 (2.5%)

¹ Denominators are based on number of subjects eligible for follow-up.

² X-rays were not required as part of the 04-04 Protocol.

Time frames for each interval are as follows: Procedure: (days 0-0), Post-Procedure: (days 1-14), 1 Month: (days 15-60), 6 Months: (days 61-244), 12 Months: (days 245-548) and 24 Months: (days 549-912).

Table 20. Aneurysm Diameter Change Through Two Years

Post-Treatment: Imaging Core Laboratory: 04-04 Test Subjects

	6 Months		12 Months		24 Months	
	Test	Control	Test	Control	Test	Control
Number of Subjects Enrolled	139	120	139	120	139	120
Subjects with Available Data ¹	103	88	114	94	104	100
≥5 mm Increase	2 (1.9%)	2 (2.3%)	0 (0.0%)	3 (3.2%)	3 (2.9%)	19 (19.0%)
No Change	76 (73.8%)	75 (85.2%)	69 (60.5%)	76 (80.9%)	50 (48.1%)	60 (60.0%)
≥5 mm Decrease	25 (24.3%)	11 (12.5%)	45 (39.5%)	15 (16.0%)	23 (46.0%)	21 (21.0%)

¹ Subjects are considered to be available for evaluation if they had both a baseline (one month) and post-baseline measurement.

Table 21. Aneurysm Volume Change Through Two Years

Post-Treatment: Imaging Core Laboratory: 04-04 Test Subjects

	6 Months		12 Months		24 Months	
	Test	Control	Test	Control	Test	Control
Number of Subjects Enrolled	139	120	139	120	139	120
Overall Aneurysm Change ¹						
Subjects with Available Data ²	102	86	112	92	104	101
≥5% Increase	6 (5.9%)	20 (23.3%)	9 (8.0%)	31 (33.7%)	17 (16.3%)	47 (46.5%)
No Change	29 (28.4%)	32 (37.2%)	26 (23.2%)	26 (28.3%)	15 (14.4%)	23 (22.8%)
≥5% Decrease	67 (65.7%)	34 (39.5%)	77 (68.8%)	35 (38.0%)	72 (69.2%)	31 (30.7%)

¹ Overall volume increase is computed as follows: if either distal renal to aortic bifurcation or right hypogastric to aortic bifurcation volume measurement increases ≥ 5% then the subject is considered an increaser. If either measure decreases ≥ 5% then the subject is considered a decreaser. Other non-missing combinations are considered No Change.

² Subjects are considered to be available for evaluation if they had both a baseline (one month) and post-baseline measurement.

CLINICAL STUDY OUTCOMES (03-02)

The GORE® EXCLUDER® Device product line has been expanded to include a 31 mm Trunk-Ipsilateral Leg and 32 mm Aortic Extender to treat patients with larger proximal aortic neck diameters. The 03-02 study was initiated to evaluate the 31 mm GORE® EXCLUDER® Device compared to open surgical repair and to compare device performance to the Original GORE® EXCLUDER® Device in the treatment of AAAs with large proximal neck diameters. Enrollment of Test patients into this study was completed in January 2008 with a total of 35 patients enrolled. The larger size GORE® EXCLUDER® Devices were approved by the FDA in March 2009.

Potential Device or Procedure Related Adverse Events

Adverse events that may occur and / or require intervention or additional intraoperative procedure time include, but are not limited to:

- allergic reaction and / or anaphylactoid response to x-ray contrast dye, anti-platelet therapy, device materials
- amputation
- anesthetic complications
- aneurysm enlargement
- aneurysm rupture and death
- arterial or venous thrombosis and / or pseudoaneurysm
- arteriovenous fistula
- bleeding, hematoma, or coagulopathy
- bowel (e.g., ileus, transient ischemia, infarction, necrosis)
- cardiac (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension or hypertension)
- claudication (e.g., buttock, lower limb)
- death
- dissection, perforation, or ruptures of the aortic vessel and surrounding vasculature
- edema
- embolization (micro and macro) with transient or permanent ischemia
- endoleak
- endoprosthesis or delivery system: improper component placement; incomplete component deployment; unintentional/premature component deployment; leading end catheter component retention; component migration; separation of graft material from stent; occlusion; infection; stent fracture; graft material failure, dilatation, erosion, puncture, perigraft flow
- fever and localized inflammation
- genitourinary (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- hepatic failure
- impotence
- infection (e.g., aneurysm, device or access sites)
- lymph fistula / complications
- multi-system organ failure
- neurologic damage, local or systemic (e.g., stroke, paraplegia, paraparesis)
- occlusion of device or native vessel
- post-implant syndrome
- pulmonary complications (e.g., pneumonia, respiratory failure)
- radiation injury, late malignancy
- renal (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- surgical cut down, bypass or conversion
- tissue necrosis
- wound (e.g., infection, dehiscence)
- vascular spasm or vascular trauma (e.g., ilio-femoral vessel dissection, bleeding, rupture, death)

Device Related Adverse Event Reporting

Any adverse event involving the GORE® EXCLUDER® AAA Endoprosthesis should be reported to W. L. Gore & Associates immediately. To report an event in the US, call 800.437.8181. Outside the US, contact your local technical representative.

PATIENT SELECTION AND TREATMENT

(SEE WARNINGS AND PRECAUTIONS)

Individualization of Treatment

Gore recommends that the GORE® EXCLUDER® AAA Endoprosthesis diameter be at least 2 mm larger than the aortic inner diameter (10 – 21% oversizing) and 1 mm larger than the iliac inner diameter (7 – 25% oversizing) as described in Tables 22 and 23. The length of the GORE® EXCLUDER® AAA Endoprosthesis should be sufficient to reach from just inferior to the most distal (lowest) major renal artery to non-aneurysmal tissue in the common or external iliac arteries. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters / lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

Additional considerations for patient selection include but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- Patient's anatomical suitability for endovascular repair
- The risk of aneurysm rupture compared to the risk of treatment with the GORE® EXCLUDER® AAA Endoprosthesis
- Ability to tolerate general, regional, or local anesthesia
- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and / or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 12 Fr - 18 Fr vascular introducer sheath
- An infrarenal non-aneurysmal aortic neck length of at least 15 mm and a diameter of no greater than 32 mm
- Proximal neck angulation ≤ 60° with minimal thrombus or calcification

- Distal segment iliac vessel lengths of at least 30 mm of which at least 10 mm must be less than or equal to 25 mm in diameter
- Freedom from significant femoral / iliac artery occlusive disease that would impede inflow or outflow of stent-grafts

The final treatment decision is at the discretion of the physician and patient.

PATIENT COUNSELING INFORMATION

The physician and patient should review the risks and benefits when discussing this endovascular device and procedure including:

- Risks and differences between endovascular repair and open surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to post-operative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term safety and effectiveness of endovascular repair has not been established.** Physicians should advise all patients that this treatment modality requires long-term, regular follow-up to assess patients' health status and stent-graft performance. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should receive enhanced follow-up. Patients should be counseled on the need for regular follow-up, even in the absence of obvious symptoms, e.g., pain, numbness, weakness (see IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP).
- Regular follow-up including imaging of the device should be performed at least every 12 months for all patients and at least every 6 to 12 months for patients with known endoleaks or aneurysm enlargement for the duration of the implant (see IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP).
- Physicians must advise all patients that it is important to seek prompt medical attention if he / she experiences signs of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking, or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as pain, numbness, weakness in the legs; any back, chest, abdominal, or groin pain, dizziness, fainting, rapid heartbeat, or sudden weakness.

US physicians are encouraged to refer the patient to the Patient Brochure regarding risks occurring during or after implantation of the device. Procedure related risks include cardiac, pulmonary, neurologic, bowel, and bleeding complications. Device related risks include occlusion, endoleak, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture and death (See POTENTIAL DEVICE OR PROCEDURE RELATED ADVERSE EVENTS).

US physicians are encouraged to complete the Patient Wallet Card and give it to the patient so that he / she can carry it with them at all times. The patient should refer to the wallet card any time they visit additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

HOW SUPPLIED

The GORE® EXCLUDER® AAA Endoprosthesis is preloaded on a delivery catheter and supplied sterile and non-pyrogenic.

Storage and Handling

- Do not resterilize; for single use only.
- Do not use if damaged or if sterile barrier has been compromised.
- Do not use after the "use by" (expiration) date printed on the label.
- Store in a cool, dry place.

CLINICAL USE INFORMATION

Physician Training Program

CAUTION: Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The GORE® EXCLUDER® AAA Endoprosthesis should only be used by physicians experienced in vascular interventional techniques, and who have successfully completed the appropriate physician training program.

The recommended skill / knowledge requirements for physicians using the GORE® EXCLUDER® AAA Endoprosthesis are outlined below:

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair
- Knowledge of radiographic image interpretation, device selection and sizing

A multi-disciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy, and repair
- Percutaneous access and closure techniques
- Non-selective and selective guidewire and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement

- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

Recommended Materials

- 0.035" (0.89 mm) 'super stiff' guidewire, 145 cm or longer
- Angiographic radiopaque marker catheter
- Contrast media
- Snare Catheter
- Syringe
- Heparin and heparinized saline
- **Trunk-Ipsilateral Leg Endoprosthesis and Contralateral Leg Endoprosthesis:**
 - 18 Fr and 12 Fr, 14 Fr or 15 Fr introducer sheaths (Tables 22 and 23)
 - Large diameter, low pressure aortic balloon (monitor balloon volumes and pressures as recommended in balloon catheter Instructions for Use)
 - Percutaneous transluminal angioplasty (PTA) balloons (Table 23)
- **Aortic Extender Endoprosthesis:**
 - 16 Fr, 17 Fr or 18 Fr introducer sheath
 - Large diameter, low pressure aortic balloon (monitor balloon volumes and pressures as recommended in balloon catheter Instructions for Use)
- **Iliac Extender Endoprosthesis:**
 - 12 Fr, 14 Fr and 15 Fr introducer sheaths
 - PTA balloon catheters, 10 mm x 40 mm, 12 mm x 40 mm, 14 mm x 40 mm, 16 mm x 40 mm, 18 mm x 40 mm, 20 mm x 40 mm, 24 mm x 40 mm, and 28 mm x 40 mm (Table 25)

Table 22. Trunk-Ipsilateral Leg Endoprosthesis Sizing Guide*

Intended Aortic Vessel Diameter (mm)	Aortic Endoprosthesis Diameter ¹ (mm)	Intended Iliac Vessel Diameter (mm)	Iliac Endoprosthesis Diameter ² (mm)	Overall Device Lengths (cm)	Recommended Introducer Sheath ³ (Fr)
19 - 21	23	10 - 11	12	12, 14, 16, 18	18
		12 - 13.5	14.5		
22 - 23	26	10 - 11	12	12, 14, 16, 18	18
		12 - 13.5	14.5		
24 - 26	28.5	10 - 11	12	12, 14, 16, 18	18
		12 - 13.5	14.5		
27 - 29	31	12 - 13.5	14.5	13, 15, 17	18
30 - 32	35	12 - 13.5	14.5	14, 16, 18	18

1 Recommended endoprosthesis oversizing relative to the aortic vessel is approximately 10-21%, and for the iliac vessel approximately 7-25%.

2 Recommended angioplasty balloon size is 12 mm and 14 mm respectively.

3 GORE® Introducer Sheaths are recommended. The GORE® EXCLUDER® AAA Endoprosthesis is not compatible with introducer sheaths longer than 40 cm (total length including the hemostatic valve).

* Note: All dimensions are nominal.

Table 23. Contralateral Leg Endoprosthesis Sizing Guide*

Intended Iliac Vessel Diameter (mm)	Iliac Endoprosthesis Diameter ¹ (mm)	Overall Device Lengths ^{2,3} (cm)	Recommended Contralateral Introducer Sheath ⁴ (Fr)	Recommended Angioplasty Balloon Size (mm x mm)
10 - 11	12	10, 12, 14	12	12 x 40
12 - 13.5	14.5	10, 12, 14	12	14 x 40
13.5 - 14.5	16	9.5, 11.5, 13.5	12	16 x 40
14.5 - 16.5	18	9.5, 11.5, 13.5	12	18 x 40
16.5 - 18.5	20	9.5, 11.5, 13.5	12	20 x 40
18.5 - 21.5	23	10, 12, 14	14	24 x 40
21.5 - 25	27	10, 12, 14	15	28 x 40

1 Recommended endoprosthesis oversizing relative to the vessel is approximately 7-25%.

2 Total treatable lengths include 4 cm of the trunk region of the Trunk-Ipsilateral Leg Endoprosthesis.

3 Labeled Contralateral Leg length includes 3 cm overlap.

4 GORE® Introducer Sheaths are recommended. The GORE® EXCLUDER® AAA Endoprosthesis is not compatible with introducersheaths longer than 40 cm (total length including the hemostatic valve).

* Note: All dimensions are nominal.

Table 24. Aortic Extender Endoprosthesis Sizing Guide*

Intended Aortic Vessel Diameter (mm)	Aortic Extender Diameter ¹ (mm)	Endoprosthesis Length (cm)	Recommended Introducer Sheath ² (Fr)
19 - 21	23	3.3	16
22 - 23	26	3.3	16
24 - 26	28.5	3.3	16
27 - 29	32	4.5	17
30 - 32	36	4.5	18

1 Recommended endoprosthesis oversizing relative to the vessel diameter is approximately 10-21%.

2 GORE® Introducer Sheaths are recommended. The GORE® EXCLUDER® AAA Endoprosthesis is not compatible with introducer sheaths longer than 40 cm (total length including the hemostatic valve).

* Note: All dimensions are nominal.

Table 25. Iliac Extender Endoprosthesis Sizing Guide*

Intended Iliac Vessel Diameter (mm)	Distal Iliac Extender Diameter ¹ (mm)	Endoprosthesis Length ² (cm)	Recommended Introducer Sheath ⁴ (Fr)	Recommended Balloon Size (Proximal) (mm)	Recommended Angioplasty Balloon Size (Distal) (mm x mm)
8 - 9	10	7	12	14	10 x 40
10 - 11	12	7	12	14	12 x 40
12 - 13.5	14.5	7	12	14	14 x 40
13.5 - 14.5	16 ^{3,5}	9.5, 11.5, 13.5	12	14	16 x 40
14.5 - 16.5	18 ^{3,5}	9.5, 11.5, 13.5	12	14	18 x 40
16.5 - 18.5	20 ^{3,5}	9.5, 11.5, 13.5	12	14	20 x 40
18.5 - 21.5	23 ^{3,5}	10, 12, 14	14	14	24 x 40
21.5 - 25	27 ^{3,5}	10, 12, 14	15	14	28 x 40

1 Recommended endoprosthesis oversizing relative to the vessel diameter is approximately 7-25%.

2 7 cm long Iliac Extender Endoprosthesis provides a maximum extension of 4 cm when placed in the Trunk-Ipsilateral or Contralateral Leg Endoprosthesis; labeled length includes 3 cm overlap.

3 16, 18, 20, 23, and 27 mm Contralateral Legs may be used as Iliac Extenders.

4 GORE® Introducer Sheaths are recommended. The GORE® EXCLUDER® AAA Endoprosthesis is not compatible with introducer sheaths longer than 40 cm (total length including the hemostatic valve).

5 If extension of a 16, 18, 20, 23, or 27 mm Contralateral Leg is necessary, a 16 mm angioplasty balloon size is recommended for the proximal end of the extension.

* Note: All dimensions are nominal.

DIRECTIONS FOR USE

Pre-Treatment Planning

- Determine accurate size of anatomy and proper size of Trunk-Ipsilateral and Contralateral Endoprostheses (Tables 22 and 23) and Aortic and Iliac Extender Endoprostheses (Tables 24 and 25).
- Use high resolution, non-contrast and contrast enhanced computerized tomography (CT / CTA) at ≤ 3 mm acquisition and reconstruction collimation.
- Use multiple view, digital subtraction angiography with a radiopaque marker catheter or spiral CT multi-planar reconstruction.
- For angiography, use correct imaging angulation (cranial-caudal, lateral-oblique) to accurately identify origin of branch vessel anatomy.
- Consider breath-hold technique to optimize digital subtraction angiography image quality.

Anatomical Requirements

- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and / or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 12 Fr - 18 Fr vascular introducer sheath.
- An infrarenal, non-aneurysmal aortic neck length of at least 15 mm and an infrarenal aortic neck treatment diameter range of 19-32 mm.
- For Trunk-Ipsilateral Leg Endoprosthesis and Aortic Extender Endoprosthesis: Proximal aortic neck angulation $\leq 60^\circ$ with minimal thrombus and / or calcification.
- Key anatomic elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation, short proximal aortic neck and significant thrombus and / or calcium at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. The US clinical studies quantify significant thrombus as thrombus ≥ 2 mm in thickness and / or $\geq 25\%$ of the vessel circumference in the intended seal zone of the aortic neck. Irregular calcium and / or plaque may compromise the fixation and sealing of the implantation sites.
- Distal segment iliac vessel lengths of at least 30 mm of which at least 10 mm must be less than or equal to 25 mm in diameter for Iliac Extender Endoprosthesis: Non-aneurysmal iliac artery length ≥ 10 mm of appropriate diameter.
- Freedom from significant femoral / iliac artery occlusive disease that would impede inflow or outflow of stent-grafts.
- Ability to tolerate general, regional, or local anesthesia.
- Patient's anatomical suitability for endovascular repair.

Arterial Access and Angiography

- Following standard practices, access the intended contralateral side via a percutaneous diagnostic sheath, and perform marker catheter digital subtraction angiography (AP, oblique and lateral views as necessary) to confirm the correct device component sizing, and deployment locations. Consider breath-hold technique to optimize image quality. Leave marker catheter in place at the level of the renal arteries.
- Following standard practices, perform percutaneous access and / or surgical exposure of the vessels selected to receive the Trunk-Ipsilateral and Contralateral side introducer sheaths.
- Following the manufacturer's instructions for use, advance a 0.035" (0.89 mm) 'super stiff' guidewire, or acceptable equivalent to the level of the renal arteries.
- Following the manufacturer's instructions for use, prepare and advance the recommended introducer sheath (Tables 22 - 25) over the guidewire, through the ilio-femoral anatomy, aortic aneurysm and up to the level of the proximal aortic neck according to standard practice.
- CAUTION:** Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Use standard heparinized saline, pressure flush system technique to prevent thrombus formation in the introducer sheaths.
- Use an accurate radiopaque patient marking method to assure accurate device positioning and deployment locations.

Catheter Preparation

- Use new, sterile gloves when preparing device.

CAUTION: Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.

- Remove the appropriately sized device delivery catheters from their packaging and examine for possible damage.
- Remove protective packaging mandrel and packaging sheath(s) from the leading end of the delivery catheters (Figures 3A - 3F).
- Flush with heparinized saline through the flushing port on the trailing end of the delivery catheter (Figures 3A - 3F).
- Follow the manufacturer's recommended method for size selection, preparation and use of aortic and iliac dilation balloons. Carefully inflate the balloon to avoid complications.

Trunk-Ipsilateral Leg Endoprosthesis Positioning

- Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
- Advance the Trunk delivery catheter over a 0.035" (0.89 mm) 'super stiff' guidewire, through the 18 Fr introducer sheath into the aorta to the approximate level of intended positioning.
WARNING: Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- While maintaining the delivery catheter in position, withdraw the introducer sheath (Table 22) and visually verify that the leading end of the introducer sheath is below the distal Trunk-Ipsilateral Leg radiopaque marker.
- Magnify and center the fluoroscopic image on the proximal trunk. Reposition and rotate the Trunk-Ipsilateral delivery catheter as necessary to properly position the proximal device marker as well as orient the long contralateral, and short ipsilateral radiopaque markers and device position on the appropriate side of the anatomy. Maximize the separation between these two markers to achieve maximum lateral positioning of the iliac legs of the device. The long marker should be oriented toward the contralateral side (Figures 1A and 1B).
WARNING: Do not rotate the Trunk delivery catheter beyond 360°. Delivery system damage and / or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- It is recommended to view and confirm the distal position of the iliac end of the device relative to the internal iliac artery to ensure accurate and desired deployment position of the distal aspect of the device.
- If clinically acceptable, lower the patient's mean arterial pressure to 60 - 70 mm Hg during Trunk deployment and aortic balloon inflation to decrease blood flow and reduce the risk of endoprosthesis movement.
- Maintain a contralateral access side sheath, catheter or guidewire in position across the distal, native bifurcation to protect and ensure that contralateral access is maintained into the aneurysm sac and contralateral leg hole of the device during Trunk-Ipsilateral component deployment.
- Re-center and magnify the image on the proximal Trunk of the device to assure final desired position of proximal device relative to anatomy. Stabilize the Trunk delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.
WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- When using the Original Delivery System (Figure 3A), please refer to the deployment instructions entitled, ORIGINAL DELIVERY SYSTEM - TRUNK IPSILATERAL ENDOPROSTHESIS DEPLOYMENT.
- When using the GORE® C3® Delivery System (Figure 3B), please refer to the deployment instructions entitled, GORE® C3® DELIVERY SYSTEM - TRUNK-IPSILATERAL ENDOPROSTHESIS DEPLOYMENT.

ORIGINAL DELIVERY SYSTEM - TRUNK IPSILATERAL ENDOPROSTHESIS DEPLOYMENT

- Loosen the deployment knob. Confirm final device position and orientation and deploy the Trunk using a steady and continuous pull of the deployment knob to release the endoprosthesis. Pull the deployment knob straight out from the catheter handle. Deployment initiates from the leading end toward the trailing end.
WARNING: Do not attempt to reposition the endoprosthesis after deployment has been initiated. Vessel damage or device misplacement may result.
WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.

CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.

2. Use fluoroscopic guidance during the withdrawal of the delivery catheter to assure safe removal from, and to avoid catching on, the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop and withdraw delivery catheter and introducer sheath together.


WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.

3. When withdrawing the device delivery system through the introducer sheath, verify all catheter components are intact.

WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred, see ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).

4. Position the aortic balloon inside the proximal region of the trunk. Avoid balloon contact with the flow splitter which is aligned with the long and short radiopaque markers. Inflate and deflate the balloon quickly with dilute contrast solution to seat the aortic end of the endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of aortic and iliac dilation balloons, carefully monitoring both volume and pressure to avoid complications.
5. Use fluoroscopic guidance to ensure the balloon is completely deflated and is safely removed from the endoprosthesis.
6. Advance and inflate the appropriate size PTA balloon catheter to seat the iliac end of the endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of PTA balloons. Carefully inflate the balloon to avoid complication.

Gore® C3® DELIVERY SYSTEM - TRUNK IPSI LATERAL ENDOPROSTHESIS DEPLOYMENT

1. Loosen the White Outer Deployment Knob by turning it 90° counter-clockwise (Figure 8A). Confirm final device position and orientation and deploy the Trunk using a steady and continuous pull of the deployment knob to release the proximal endoprosthesis past the contralateral leg hole of the Trunk. Pull the deployment knob straight out from the catheter handle. Deployment initiates from the leading end toward the trailing end. The ipsilateral leg will remain constrained on the delivery catheter. If the Trunk is not properly positioned following removal of the white deployment knob, continue with Optional Repositioning of Trunk-Ipsilateral Leg Endoprosthesis, described below.

WARNING: Do not attempt to reposition the endoprosthesis after complete deployment of the device. Vessel damage or device misplacement may result.

WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.

WARNING: Do not attempt to remove a partially deployed Trunk component.

CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.

BACKUP DEPLOYMENT MECHANISM: In the event that the removal of the White Outer Deployment Knob does not initiate the deployment of the proximal end of the Trunk device, and the problem is suspected to have occurred in the deployment handle, the deployment handle has a design feature that is intended to protect the patient against this problem becoming a hazard. Remove the Deployment Line Access Hatch Lid (Figure 9) and then pull the first deployment line (labeled "1") with an appropriate dull or blunt-edged tool. Fully remove the deployment line in a steady, continuous motion. This action will cause the proximal Trunk to deploy past the contralateral leg hole.

WARNING: Do not cut the first deployment line (labeled "1"). Patient injury, catheter damage, and / or device damage may occur.

Optional Repositioning of Trunk-Ipsilateral Leg Endoprosthesis

- A. To constrain the proximal end of the Trunk, turn the Gray Constraining Dial clockwise until it stops (Figure 8B). The Black Nut within the Transparent Knob housing will move toward the trailing end of the catheter when the Trunk end of the device is constrained. Verify the distance between the three proximal radiopaque markers is reduced (Figures 1A and 1B).
- B. Magnify and center the fluoroscopic image on the proximal Trunk. Manually reposition and / or rotate the Trunk-Ipsilateral component delivery catheter as necessary to properly position the Trunk device. Maximize the separation between the long and short radiopaque markers to achieve maximum lateral iliac leg positioning of the device. The long marker should be oriented toward the contralateral side (Figures 1A and 1B).
- C. Repeat steps 4 through 8 listed above to properly position the endoprosthesis while the Trunk end of the endoprosthesis is constrained.
- D. When the Trunk is properly positioned, turn the Gray Constraining Dial counter-clockwise until it stops to reopen the proximal end of the Trunk (Figure 8C). The Black Nut within the Transparent Knob housing will move toward the leading end of the catheter when the proximal end of the Trunk is reopened. Verify that the distance between the three proximal radiopaque markers is maximized (Figures 1A and 1B).

WARNING: The constraining mechanism cannot be removed unless the proximal end of the Trunk device is fully opened.

WARNING: Do not constrain / reopen the Trunk device more than two times during a procedure. Device and / or catheter damage may occur.

WARNING: In the event that the complete counter-clockwise rotation of the Gray Constraining Dial does not initiate the reopening of the proximal end of the Trunk device, and the problem is suspected to have occurred in the deployment handle, the deployment handle has a design feature that is intended to protect the patient against this problem becoming a hazard. Verify the final desired Trunk device position relative to anatomy. Stabilize the Trunk delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site. Disengage the constraining mechanism by sliding the Red Safety Lock back while rotating the entire Transparent Knob 90° counter-clockwise (Figure 8D). Remove the constraining mechanism by pulling the Transparent Knob straight out from the catheter handle (Figure 8E). These actions will cause the constraining mechanism to be removed and the proximal end of the Trunk device to reopen.

Contralateral Leg Hole Guidewire Cannulation

1. Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
2. Following manufacturer's instructions for use, advance a 0.035" (0.89 mm) 'super stiff' guidewire into the contralateral leg hole of the Trunk according to standard practice.
3. Verify that the guidewire is within the contralateral leg hole of the Trunk by standard practice used to verify guidewire location. If multiple cannulation attempts fail, the Trunk-Ipsilateral Leg Endoprosthesis may be repositioned for better contralateral gate access by following steps A through D above.
WARNING: Do not constrain / reopen the Trunk device more than two times during a procedure. Device and / or catheter damage may occur.
4. Following manufacturer's instructions for use, introduce the recommended introducer sheath (Table 23). Advance the sheath over the guidewire and through the contralateral leg hole of the Trunk.

Ipsilateral Leg Deployment

1. Verify the proximal end of the Trunk is fully open by turning the Gray Constraining Dial counter-clockwise until it stops (Figure 8C). The Black Nut within the Transparent Knob housing will move toward the leading end of the catheter when the proximal end of the Trunk is opened. Verify that the distance between the three proximal radiopaque markers is maximized (Figures 1A and 1B).
2. Disengage the constraining mechanism by sliding the Red Safety Lock back while rotating the entire Transparent Knob 90° counter-clockwise (Figure 8D). Remove the constraining mechanism by pulling the Transparent Knob straight out from the catheter handle (Figure 8E).

WARNING: The constraining mechanism cannot be removed unless the proximal end of the Trunk device is fully opened. Verify the Gray Constraining Dial is rotated counter-clockwise until it stops, prior to removing the constraining mechanism. Attempting to remove the constraining mechanism while the Trunk is constrained may result in misplacement of the device, patient injury, and / or device damage.

BACKUP DEPLOYMENT MECHANISM: In the event, that the removal of the Transparent Knob does not initiate the removal of the constraining mechanism and the problem is suspected to have occurred in the deployment handle, the deployment handle has a design feature that is intended to protect the patient against this problem becoming a hazard. Remove the Deployment Line Access Hatch Lid (Figure 9) and then cut the silver Lock Pin line (labeled "2") with an appropriate tool. Grasp, pull, and fully remove the silver Lock Pin line (labeled "2") in a steady, continuous motion. Cut, grasp, pull, and fully remove the gray Constraining Loop line (labeled "3") in a steady, continuous motion. These actions will cause the constraining mechanism to be removed.

3. Loosen the gray deployment knob (Figure 8F). Deploy the Ipsilateral Leg of the Trunk device using a steady and continuous pull of the gray deployment knob to release the endoprosthesis. Pull the deployment line straight out from the catheter handle.

BACKUP DEPLOYMENT MECHANISM: In the event that the removal of the gray deployment knob does not initiate the deployment of the ipsilateral leg and the problem is suspected to have occurred in the deployment handle, the deployment handle has a design feature that is intended to protect the patient against this problem becoming a hazard. Remove the Deployment Line Access Hatch Lid (Figure 9) and then pull the second deployment line (labeled "4") with an appropriate dull or blunt-edged tool. Fully remove the deployment line in a steady, continuous motion. This action will cause the ipsilateral leg to deploy.

WARNING: Do not cut the second deployment line (labeled "4"). Patient injury, catheter damage, and / or device damage may occur.

4. Use fluoroscopic guidance during the withdrawal of the delivery catheter to assure safe removal from, and to avoid catching on, the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop, visually assess the cause of the resistance, and withdraw delivery catheter and introducer sheath together if appropriate.
WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.
5. When withdrawing the device delivery system through the introduce sheath, verify all catheter components are intact.

WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred, see ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).

- Position the aortic balloon inside the proximal region of the trunk. Avoid balloon contact with the flow splitter which is aligned with the long and short radiopaque markers. Inflate and deflate the balloon quickly with dilute contrast solution to seat the aortic end of the endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of aortic and iliac dilation balloons, carefully monitoring both volume and pressure to avoid complications.
- Use fluoroscopic guidance to ensure the balloon is completely deflated and is safely removed from the endoprosthesis.
- Advance and inflate the appropriate size PTA balloon catheter to seat the iliac end of the endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of PTA balloons. Carefully inflate the balloon to avoid complication.

Contralateral Leg Endoprosthesis Positioning and Deployment

- Advance the prepped Contralateral Leg Endoprosthesis delivery catheter to the level of the long radiopaque marker (Figure 1).
WARNING: Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment may occur have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Align the radiopaque marker at the proximal end of the Contralateral Leg Endoprosthesis with the long contralateral radiopaque marker on the Trunk-Ipsilateral Leg Endoprosthesis. With the alignment of these markers, an approximate 3 cm overlap will be achieved.
- While maintaining the delivery catheter in position, withdraw the introducer sheath (Table 23) and visually verify that the leading end of the introducer sheath is below the distal Contralateral Leg radiopaque marker.
WARNING: Do not rotate the Contralateral Leg delivery catheter during delivery, positioning or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Stabilize the Contralateral Leg Endoprosthesis delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.
WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Loosen the deployment knob. Confirm final device position. Deploy the Contralateral Leg Endoprosthesis by using a steady, continuous pull of the deployment knob to release the endoprosthesis. Pull the deployment knob straight out from the catheter side-arm. Deployment initiates from the leading (aortic) end toward the trailing (iliac) end.
WARNING: Do not attempt to reposition the endoprosthesis after deployment has been initiated. Vessel damage or device misplacement may result.
WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.
- Use fluoroscopic guidance during withdrawal of the delivery catheter to assure safe removal from the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop and withdraw delivery catheter and introducer sheath together.
WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.
- When withdrawing the device delivery system through the introduce sheath, verify all catheter components are intact.
WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred, see ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
- Following manufacturer's instructions for use, advance and inflate a 14 mm PTA balloon catheter to seat the proximal end of the Contralateral Leg Endoprosthesis within the contralateral leg hole overlap region. Follow the manufacturer's

recommended method for size selection, preparation and use of aortic and iliac dilation balloons, carefully monitoring both volume and pressure to avoid complications.

- Following manufacturer's instructions for use, advance and inflate the appropriate size PTA balloon to seat the iliac end of the Contralateral Leg Endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of PTA balloons. Carefully inflate the balloon to avoid complications.

Aortic Extender Endoprosthesis Positioning and Deployment

- Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
- Advance the Aortic Extender Endoprosthesis delivery catheter over a 0.035" (0.89 mm) 'super stiff' guidewire, through the 16 Fr, 17 Fr or 18 Fr introducer sheath into the aorta, just proximal to the level of intended device positioning.
WARNING: Do not attempt to advance the Aortic Extender through the 12 Fr introducer sheath. The Aortic Extender is designed for a 16 Fr, 17 Fr, or 18 Fr sheath.
WARNING: Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment may occur have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- While maintaining the delivery catheter in position, withdraw the introducer sheath (Table 24) and visually verify that the leading end of the introducer sheath is below the distal Aortic Extender radiopaque marker.
WARNING: Do not rotate the Aortic Extender delivery catheter during delivery, positioning or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Magnify and center the fluoroscopic image on the proximal Aortic Extender endoprosthesis. Reposition the Aortic Endoprosthesis delivery catheter as necessary to position the proximal and distal radiopaque markers in appropriate position. The maximum recommended extension with each Aortic Extender component is approximately one-half of the Extender length inside (1.6 cm or 2.2 cm) and one-half, outside (1.6 cm or 2.2 cm), or proximal to the Trunk or Aortic Extender host component. The proximal three (3) and distal one (1) markers are visible relative to host device and anatomy pre and post-deployment (Figures 4 and 5B).
- If clinically acceptable, lower the patient's mean arterial pressure to 60 – 70 mm Hg during Aortic Extender deployment and aortic balloon inflation to decrease blood flow and reduce the risk of endoprosthesis movement.
- Stabilize the Extender delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.
WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Loosen the deployment knob. Using fluoroscopy, confirm final device position and deploy the Aortic Extender using a steady and continuous pull of the deployment knob to release the endoprosthesis. Pull the deployment knob straight out from the catheter side-arm. Deployment initiates from the trailing end of the device toward the leading end of the device.
WARNING: Do not attempt to reposition the endoprosthesis after deployment has been initiated. Vessel damage or device misplacement may result.
WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.
- Use fluoroscopic guidance during the withdrawal of the delivery catheter to assure safe removal from, and to avoid catching on, the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop and withdraw delivery catheter and introducer sheath together.
WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.
- When withdrawing the device delivery system through the introduce sheath, verify all catheter components are intact.
WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred, see ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
- Advance the aortic dilation balloon until it is centered relative to the endoprosthesis. Inflate and deflate the balloon quickly with dilute contrast

solution to seat the Aortic Extender Endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of aortic dilation balloons. Carefully inflate the balloon to avoid complications.

- Use fluoroscopic guidance to ensure the balloon is completely deflated and to assure safe removal from, and to avoid catching on, the endoprosthesis.

Iliac Extender Endoprosthesis Positioning and Deployment

- Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
- Advance the Iliac Extender Endoprosthesis delivery catheter into the distal end of the host device, via the recommended introducer sheath (Table 25).
WARNING: Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- For maximum extension, align the radiopaque marker at the iliac (distal) end of the host device with the marker located 3 cm below the proximal end of the Extender component (Figures 2B, 3D, 6A and 6B).
WARNING: While using 16, 18, 20, 23, or 27 mm Contralateral Legs as an Iliac Extender, ensure that the distal end including the taper zone will not be deployed inside the previously deployed Ipsilateral Leg or Contralateral Leg of the GORE® EXCLUDER® AAA Endoprosthesis. However, when the Contralateral Leg and Iliac Extender diameters are identical, the taper zone can be deployed inside the previously deployed Contralateral Leg (Figure 7B).
WARNING: While using 16, 18, 20, 23, or 27 mm Contralateral Legs as an Iliac Extender, the 3 cm mandatory overlap must be achieved prior to the beginning of the distal taper zone of the 18, 20, 23, and 27 mm Contralateral Leg. Inadequate sealing may lead to endoleak.
- While maintaining the delivery catheter in position, withdraw the introducer sheath (Table 25) and visually verify that the leading end of the introducer sheath is below the distal Iliac Extender radiopaque marker.
WARNING: Do not rotate the Iliac Extender delivery catheter during delivery, positioning or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Stabilize the Iliac Extender delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.
WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Loosen the deployment knob. Confirm final device position. Using fluoroscopy, deploy the Iliac Extender Endoprosthesis by using a steady, continuous pull of the deployment knob to release the endoprosthesis. Pull the deployment knob straight out from the catheter side-arm. The device deploys from the leading (proximal) end toward the trailing (distal) end.
WARNING: Do not attempt to reposition the endoprosthesis after deployment has been initiated. Vessel damage or device misplacement may result.
WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.
- Use fluoroscopic guidance during the withdrawal of the delivery catheter to assure safe removal from, and avoid catching on, the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop and withdraw delivery catheter and introducer sheath together.
WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.
- When withdrawing the device delivery system through the introduce sheath, verify all catheter components are intact.
WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
- Advance and inflate an appropriate size PTA balloon catheter to seat the proximal overlap end and the distal end of the Iliac Extender Endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of PTA balloons. Carefully inflate the balloon to avoid complications.
- Use fluoroscopic guidance to ensure the balloon is completely deflated and to assure safe removal from, and to avoid catching on, the endoprosthesis.

Completion of the Procedure

- Perform extended imaging angiography to confirm exclusion of the aneurysm. Consider breath-hold technique to optimize digital subtraction angiography image quality. Consider use of GORE® EXCLUDER® AAA Endoprosthesis Extender components as necessary. For Aortic Extenders, a minimum overlap of 1.6 cm or 2.2 cm is required, offering a maximum of 1.6 cm or 2.2 cm of extension; for Iliac Extenders, a minimum overlap of 3 cm is required.
- Prior to removing wires and sheaths verify all catheter components have been withdrawn from the patient.
WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
- Close arterial access according to standard practice.
- Prior to removing patients as necessary to provide proper surveillance of the long-term performance of the endoprosthesis, procedure and status of the aneurysm. Minimally, annual CTs, multiple view X-rays, and ultrasound may be used for such surveillance.

IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

General

The long-term safety and effectiveness of endovascular repair has not been established. All patients should be advised this treatment modality requires long-term, regular follow-up to assess patients' health status and stent-graft performance. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should receive enhanced follow-up. Patients should be counseled on the need for regular follow-up, even in the absence of obvious symptoms (e.g., pain, numbness, weakness).

Regular and consistent follow-up is a critical part of ensuring continuing safety and efficacy of aortic endovascular repair. Physicians should tailor patient follow-up to the needs and circumstances of each individual patient. In the US clinical studies, at least one annual physician visit and the imaging schedule (Table 26) were employed.

Follow-up modalities include CT / CTA, multi-view abdominal X-ray, MRI / MRA, and ultrasound. Data from these modalities is acquired and used to compare baseline and subsequent exams to review devices and morphological changes over time and their effects on exclusion of the aneurysm.

- CT / CTA imaging provides information on aneurysm size, vascular morphological changes, proximal device-trunk fixation and migration, endoleak and patency / limb occlusion.
- Multi-view device X-ray imaging provides information on device wireform integrity (e.g., fracture, kinking) and relative component migration.
- MRI / MRA imaging provides information similar to CT / CTA and is often used as a surrogate for CT / CTA if patients are CT contrast intolerant.
- Ultrasound may be used to assess for endoleak and aneurysm size status but not for device integrity, specifically the wire form. Ultrasound is a less reliable and sensitive diagnostic method compared to CT.

Alternative imaging recommendations for patients with CT or angiography contrast intolerance issues include CO₂ angiography, MRI-MRA with or without contrast, and ultrasound. These imaging and surveillance modalities may be less sensitive and difficult to compare with diagnostic findings from previous or subsequent follow-up exams.

Table 26. Recommended Schedule for Patient Imaging Follow-Up

Schedule for Patient Imaging Follow-Up			
Visit	Angiogram	Abdominal X-ray	CT Pre-Contrast and Contrast
Pre-Treatment	X ¹		X ¹
Treatment (Pre and Post Deployment)	X		
Discharge		X	
1 Month			X
3 Month			X ²
6 Month		X	X
12 Month (Annually Thereafter)		X	X

¹ Imaging should be performed ≤ six months prior to the procedure

² Recommended if endoleak reported at one month

Angiographic Imaging

Angiographic images are recommended pre-treatment to evaluate the length and tortuosity of abdominal aorta, iliac and common femoral arteries.

- Images should include an angiographic marker catheter with incremental one centimeter markers over a 10 – 20 cm length.
- The following views are recommended for optimal evaluation and case planning:
 - Abdominal aorta; Supine-AP, Lateral
 - Pelvis (to include bilateral common femorals); AP, both Obliques

Angiographic images are recommended during the treatment procedure both pre and post-deployment to evaluate device placement and orientation. Selective angiography during subsequent follow-up exams may provide useful device position and device integrity information.

CT / CTA Images

- Film sets should include all sequential images at lowest possible slice thickness (≤ 3 mm). Do NOT perform large slice thickness (> 3 mm) and / or omission of CT images / film sets (non-consecutive) as it prevents precise anatomical and device comparisons over time.
- All images should include a scale for each film / image. Images should be arranged no smaller than 20:1 images on 14" x 17" sheets if film is used.
- If an endoleak is suspected or there is aneurysm enlargement, it is recommended that pre-contrast and contrast runs be performed.**
- Pre-contrast and contrast run slice thickness and interval must match.
- DO NOT change patient orientation or re-landmark patient between non-contrast and contrast runs.

Non-contrast and contrast enhanced baseline and follow-up exams are important for optimal patient surveillance. For the best results, use the following CT / CTA imaging guidelines listed in Table 27.

Table 27. CT / CTA Imaging Guidelines

	Pre-Contrast	CT / CTA
IV Contrast	No	Yes
Injection Volume (ml)	NA	150
Injection Rate (cc / sec)	NA	≥ 2.5
Delay	NA	Smart-Prep [®] , CARE or equivalent
Start Position	Diaphragm	1 cm above Celiac Axis
End Position	Proximal Femur	Femoral Bifurcation
Scan FOV	Large	Large
DFOV	32 cm	32 cm
Scan Type	Helical	Helical
Rotation Speed	0.8	0.8
Slice Thickness (mm)	≤ 3.0 mm	≤ 3.0 mm
Scan Mode	HS	HS
Table Speed (mm / rot)	15	15
Interval (mm)	2.0	2.0
KV / mA	120 / 300	120 / 300
Reconstruction / Algorithm	≤ 3.0 mm Soft	≤ 3.0 mm Soft
* Smart Prep ROI Loc: 1 cm Sup. to Celiac Axis Scan Phase: 3 Sec MA: 40 Monitor Delay: 6 Sec Monitor ISO: 3 Sec Enhance Thres: 100 HU		

Abdominal X-ray Film Series (plain film)

The following abdominal X-ray views are recommended for optimal visualization of the engraft:

- Supine – Frontal (AP)
- Lateral
- 30 degree LPO
- 30 degree RPO

Ensure entire device is captured on each single image format lengthwise.

If there is any concern about the device integrity (e.g., kinking, stent-wire breaks, relative component migration), it is recommended to use magnified views. The attending physician should evaluate films for device integrity (entire device length including components) using 2 – 4x magnification.

MRI Information MR Conditional

Non-clinical testing has demonstrated that the GORE® EXCLUDER® AAA Endoprosthesis is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Highest spatial gradient magnetic field of 3000-Gauss/cm or less
- Whole body-averaged specific absorption rate (SAR) of 2.0-W/kg in the Normal Operating Mode (i.e., the mode of operation of the MR equipment in which none of the outputs have a value that cause physiological stress to patients) for 15 minutes of scanning (i.e., per pulse sequence).

MR Related Heating:

In non-clinical testing, the GORE® EXCLUDER® AAA Endoprosthesis produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numeris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and during MRI performed for 15-min of scanning 3-Tesla/128-MHz (Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

Please see TABLE 1. MRI Information

Image Artifact:

For each vascular device and assembly, the artifacts that appeared on the MR images were shown as localized signal voids (i.e., signal loss) that were minor in size relative to the size and shape of these implants. The gradient echo pulse sequence produced larger artifacts than the T1 – weighted, spin echo pulse sequence for the GORE® EXCLUDER® AAA Endoprosthesis. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10 mm relative to the size and shape of this implant. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the GORE® EXCLUDER® AAA Endoprosthesis. Therefore, it may be necessary to optimize the MR imaging parameters to compensate for the presence of this implant.

Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement ≥ 5 mm of maximum diameter (regardless of endoleak status)

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy, and the patient's personal choices. Patients should be counseled as to the possibility of subsequent reinterventions including catheter-based and open surgical conversion.

DEVICE-RELATED ADVERSE EVENT REPORTING

Any adverse event involving the GORE® EXCLUDER® AAA Endoprosthesis should be reported to W. L. Gore & Associates immediately. To report an event in the US, call 800.437.8181. Outside the US, contact your local technical representative.

PATIENT TRACKING INFORMATION

In addition to these Instructions for Use, the GORE® EXCLUDER® AAA Endoprosthesis is packaged with a Device Tracking Form which US hospital staff are required to complete and forward to Gore for the purposes of tracking all patients who receive a GORE® EXCLUDER® AAA Endoprosthesis product (as required by US Federal Regulation).

DEFINITIONS

 Use By

 Caution

 Consult Instructions for Use

 Do Not Resterilize

 Do Not Reuse

 Catalogue Number

 Batch Code

 Serial Number

 Authorised Representative in the European Community

 MR Conditional

 **CAUTION:** USA Federal Law restricts the sale, distribution, or use of this device to, by, or on the order of a physician.

 Sterile

 Sterilized using Ethylene Oxide

 Do Not Use if Package is Damaged

 Keep Dry

 Store in a Cool Place

 Catheter Working Length

 Delivery Profile

 Guidewire Compatibility

 Manufacturer

 Date of Manufacture

TABLE OF CONTENTS

DESCRIPTION.....	36
INDICATIONS FOR USE	36
CONTRAINDICATIONS.....	36
WARNINGS AND PRECAUTIONS.....	36
CLINICAL STUDY OUTCOMES* (98-03 and 99-04)	37
CLINICAL STUDY OUTCOMES (03-02).....	38
PATIENT SELECTION AND TREATMENT.....	39
PATIENT COUNSELING INFORMATION	39
HOW SUPPLIED.....	39
CLINICAL USE INFORMATION	39
DIRECTIONS FOR USE.....	40
IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP.....	44
DEVICE-RELATED ADVERSE EVENT REPORTING.....	45
PATIENT TRACKING INFORMATION	45
DEFINITIONS	45

- **Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious surgical consequences or injury to the patient.**

INSTRUCTIONS FOR USE

GORE® EXCLUDER® AAA ENDOPROSTHESIS

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious surgical consequences or injury to the patient.

DESCRIPTION

Trunk-Ipsilateral Leg Endoprosthesis and Contralateral Leg Endoprosthesis

The GORE® EXCLUDER® AAA Endoprosthesis provides endovascular treatment of infrarenal abdominal aortic aneurysms (AAAs).

The GORE® EXCLUDER® AAA Endoprosthesis is comprised of two components, the Trunk-Ipsilateral Leg Endoprosthesis (Trunk) (Figures 1A and 1B) and the Contralateral Leg Endoprosthesis (Figures 2A and 2B). The graft material is expanded polytetrafluoroethylene and fluorinated ethylene propylene (ePTFE and FEP), that is supported by nitinol (nickel titanium alloy) wire along its external surface. Nitinol anchors and an ePTFE / FEP sealing cuff are located at the aortic end of the trunk (Figures 1A and 1B). An ePTFE / FEP sleeve is used to constrain the endoprostheses on the leading end of the delivery catheter (Figures 3A, 3B, 3C, 3D and 3E).

Deployment of both endoprosthesis components initiates from the leading (aortic) end and proceeds toward the trailing (iliac) end of the delivery catheter (Figures 3A, 3B, 3C, 3D and 3E). The ePTFE / FEP sleeve remains in situ between the endoprosthesis and the vessel wall.

The low permeability GORE® EXCLUDER® Endoprosthesis is the only available design, and differs from the original GORE® EXCLUDER® Endoprosthesis by the addition of film layers that decrease the overall permeability of the graft. The luminal and abluminal ePTFE surface materials, microstructure, and characteristics of the low permeability GORE® EXCLUDER® Endoprosthesis meet the same specifications as the original GORE® EXCLUDER® Endoprosthesis.

Aortic Extender Endoprosthesis and Iliac Extender Endoprostheses Components

Aortic Extender Endoprosthesis

The Aortic Extender Endoprosthesis (Aortic Extender) provides an extension of approximately 1.6 cm or 2.2 cm of the leading (proximal) end of the Trunk-Ipsilateral Leg Endoprosthesis (Trunk). The Aortic Extender can be placed at variable extension lengths from 0 cm to 1.6 or 2.2 cm of the leading (proximal) end of the Trunk-Ipsilateral Leg Endoprosthesis (Trunk), allowing customization of extender treatment length based on patient anatomy and physician preference. This extension requires a minimum of approximately 1.6 cm or 2.2 cm overlap with the Trunk. The graft material is expanded polytetrafluoroethylene and fluorinated ethylene propylene (ePTFE and FEP), and is supported by nitinol (nickel titanium alloy) wire along its external surface. An ePTFE / FEP sealing cuff is located near the proximal end of the endoprosthesis (Figure 4). An ePTFE / FEP sleeve is used to constrain the endoprosthesis on the leading end of the delivery catheter (Figures 5A and 5B). Deployment of the Aortic Extender initiates from the trailing (trunk) end and proceeds toward the leading (aortic) end of the endoprosthesis and delivery catheter. Following deployment, the ePTFE / FEP sleeve remains *in situ* between the endoprosthesis and the vessel wall.

Iliac Extender Endoprosthesis

(Distal Iliac Diameters of 10, 12, and 14.5 mm)

The Iliac Extender Endoprosthesis (Iliac Extender) provides an extension of up to 4 cm of either the ipsilateral or contralateral limb. The extender component can be placed at variable extension lengths from 0 cm to 4 cm allowing customization of extender treatment length based on patient anatomy and physician preference. The graft material is ePTFE / FEP, and is supported by nitinol wire along its external surface. A radiopaque marker is located 3 cm from the proximal or top end (Figure 6A). This marker denotes the recommended minimum overlap with the ipsilateral or contralateral limb of the GORE® EXCLUDER® AAA Endoprosthesis. An ePTFE / FEP sleeve is used to constrain the endoprosthesis on the leading end of the delivery catheter (Figure 6B). Deployment of the Iliac Extender initiates from the leading (aortic) end and proceeds toward the trailing (iliac) end of the delivery catheter. Following deployment, the ePTFE / FEP sleeve remains *in situ* between the endoprosthesis and the vessel wall.

Contralateral Leg Endoprosthesis Used as Iliac Extender Endoprosthesis (Distal Iliac Diameters of 16, 18, 20, 23, and 27 mm)

Only Contralateral Leg Endoprostheses described in Figure 2B may also be utilized as Iliac Extenders as depicted in Figure 7. A minimum overlap of 3 cm between the original Ipsilateral Leg (Figure 7A) or Contralateral Leg (Figures 7B and 7C) is required. A radiopaque marker located 3 cm below the proximal end is provided to confirm this overlap. This overlap should be achieved prior to the beginning of the distal taper zone of the 18, 20, 23, and 27 mm Contralateral Leg. Further, the distal end including the taper zone should not be deployed inside the previously deployed Ipsilateral Leg or Contralateral Leg of the GORE® EXCLUDER® AAA Endoprosthesis if the distal end diameters are 12 or 14.5 mm.

The 16, 18, and 20 mm Contralateral Legs have a radiopaque marker located 4 cm above the distal end which defines the recommended minimal extension required of the Contralateral Leg as an Iliac Extender. The 23 and 27 mm Contralateral Legs have a radiopaque marker located 5 and 6 cm, respectively,

above the distal end which defines the recommended minimal extension required of the Contralateral Leg as an Iliac Extender. * However, when the Contralateral Leg and Iliac Extender diameters are identical, the taper zone can be deployed inside the previously deployed Contralateral Leg (Figure 7B)

- * **Note: The 27 mm x 10 cm Contralateral Leg Endoprosthesis has only three (3) radiopaque markers; One (1) located at each end and one (1) located 3 cm below the proximal end.**

INDICATIONS FOR USE

Trunk-Ipsilateral Leg Endoprosthesis and Contralateral Leg Endoprosthesis Components

The GORE® EXCLUDER® AAA Endoprosthesis is intended to exclude the aneurysm from the blood circulation in patients diagnosed with infrarenal abdominal aortic aneurysm (AAA) disease and who have appropriate anatomy as described below:

- Adequate iliac / femoral access
- Infrarenal aortic neck treatment diameter range of 19 – 32 mm and a minimum aortic neck length of 15 mm
- Proximal aortic neck angulation $\leq 60^\circ$
- Iliac artery treatment diameter range of 8 – 25 mm and iliac distal vessel seal zone length of at least 10 mm

Aortic Extender Endoprosthesis and Iliac Extender Endoprostheses Components

The Aortic and Iliac Extender Endoprostheses are intended to be used after deployment of the GORE® EXCLUDER® AAA Endoprosthesis. These extensions are intended to be used when additional length and / or sealing for aneurysmal exclusion is desired.

CONTRAINDICATIONS

The GORE® EXCLUDER® AAA Endoprosthesis is contraindicated in:

- Patients with known sensitivities or allergies to the device materials
- Patients with a systemic infection who may be at increased risk of endovascular graft infection.

WARNINGS AND PRECAUTIONS

General

- Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious surgical consequences or injury to the patient.
- The GORE® Medical Device is designed for single use only; do not reuse device. Gore does not have data regarding reuse of this device. Reuse may cause device failure or procedural complications including device damage, compromised device biocompatibility, and device contamination. Reuse may result in infection, serious injury, or patient death.
- **The long-term performance of stent-grafts has not been established.** All patients should be advised this treatment modality requires long-term, regular follow-up to assess patients' health status and stent-graft performance. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should receive enhanced follow-up (See IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP).
- The GORE® EXCLUDER® AAA Endoprosthesis should only be used by physicians experienced in vascular interventional techniques, and who have successfully completed the appropriate physician training program.
- The GORE® EXCLUDER® AAA Endoprosthesis is not recommended in patients unable to undergo, or who will not be compliant with the necessary pre and post-operative imaging and follow-up described in IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP.
- Intervention or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms and / or persistent endoleak. An increase in aneurysm size and / or persistent endoleak may lead to aneurysm rupture.
- Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

Patient Selection, Treatment, and Follow-Up

- The safety and effectiveness of the GORE® EXCLUDER® AAA Endoprosthesis have not been evaluated in the following patient populations:
 - traumatic aortic injury
 - leaking: pending rupture or ruptured aneurysms
 - mycotic aneurysms
 - pseudoaneurysms resulting from previous graft placement
 - revision of previously placed stent grafts
 - genetic connective tissue disease (e.g., Marfans or Ehlers-Danlos Syndromes)
 - concomitant thoracic aortic or thoracoabdominal aneurysms
 - inflammatory aneurysms
 - patients with active systemic infections
 - pregnant or nursing females
 - morbidly obese patients
 - patients less than 21 years old
 - patients with less than 15 mm in length or $> 60^\circ$ angulation of the proximal aortic neck

- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and / or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 12 Fr – 18 Fr vascular introducer sheath.
- Key anatomic elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation, short proximal aortic neck and significant thrombus and / or calcium at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. The US clinical studies quantify significant thrombus as thrombus ≥ 2 mm in thickness and / or $\geq 25\%$ of the vessel circumference in the intended seal zone of the aortic neck. Irregular calcium and / or plaque may compromise the fixation and sealing of the implantation sites.
- The GORE® EXCLUDER® AAA Endoprosthesis is not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and post-operative follow-up imaging.
- The GORE® EXCLUDER® AAA Endoprosthesis is not recommended in patients exceeding weight and / or size limits which compromise or prevent the necessary imaging requirements.
- The GORE® EXCLUDER® AAA Endoprosthesis is not recommended in patients with known sensitivities or allergies to ePTFE, FEP, nickel, or titanium.

Implant Procedure

- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not attempt to advance the Aortic Extender through the 12 Fr introducer sheath. The Aortic Extender is designed for a 16 Fr, 17 Fr or 18 Fr sheath.
- Do not rotate the Trunk, Contralateral Leg, Iliac Extender or Aortic Extender delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not rotate the Trunk delivery catheter beyond 360°. Delivery system damage and / or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not rotate the Contralateral Leg, Iliac Extender or Aortic Extender delivery catheter during delivery, positioning or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- **Do not attempt to reposition the endoprosthesis after complete deployment of the device. Vessel damage or device misplacement may result.**
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
- **Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.**
- Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
- Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.
- While using 16, 18, 20, 23, or 27 mm Contralateral Legs as an Iliac Extender, ensure that the distal end including the taper zone will not be deployed inside the previously deployed Ipsilateral Leg or Contralateral Leg of the GORE® EXCLUDER® AAA Endoprosthesis. However, when the Contralateral Leg and Iliac Extender diameters are identical, the taper zone can be deployed inside the previously deployed Contralateral Leg (Figure 7B).
- While using 16, 18, 20, 23, or 27 mm Contralateral Legs as an Iliac Extender, the 3 cm mandatory overlap must be achieved prior to the beginning of the distal taper zone of the 18, 20, 23, and 27 mm Contralateral Leg. Inadequate sealing may lead to endoleak.

MRI Information MR Conditional

Non-clinical testing has demonstrated that the GORE® EXCLUDER® AAA Endoprosthesis is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Highest spatial gradient magnetic field of 3000-Gauss/cm or less
- Whole body-averaged specific absorption rate (SAR) of 2.0-W/kg in the Normal Operating Mode (i.e., the mode of operation of the MR equipment in which none of the outputs have a value that cause physiological stress to patients) for 15 minutes of scanning (i.e., per pulse sequence).

MR Related Heating:

In non-clinical testing, the GORE® EXCLUDER® AAA Endoprosthesis produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numeris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and during MRI performed for 15-min of scanning 3-Tesla/128-MHz (Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

Table 1. MRI Information

	1.5 Tesla	3 Tesla
MR system reported, whole body averaged SAR	2.9 W/kg	2.9 W/kg
Calorimetry measured values, whole body averaged SAR	2.1 W/kg	2.7 W/kg
Highest recorded temperature change	2.2°C	2.7°C
Calculated Temperature Change Normal Operating Mode SAR = 2.0 W/kg	1.5°C	1.9°C
Calculated Temperature Change First Level Control Mode SAR = 4.0 W/kg	3.0°C	3.7°C

Image Artifact:

For each vascular device and assembly, the artifacts that appeared on the MR images were shown as localized signal voids (i.e., signal loss) that were minor in size relative to the size and shape of these implants. The gradient echo pulse sequence produced larger artifacts than the T1 – weighted, spin echo pulse sequence for the GORE® EXCLUDER® AAA Endoprosthesis. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10 mm relative to the size and shape of this implant. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the GORE® EXCLUDER® AAA Endoprosthesis. Therefore, it may be necessary to optimize the MR imaging parameters to compensate for the presence of this implant.

CLINICAL STUDY OUTCOMES* (98-03 AND 99-04)

Two US clinical studies were conducted to evaluate the safety and efficacy of the GORE® EXCLUDER® AAA Endoprosthesis. The first, referred to as 98-03, evaluated the original design of the GORE® EXCLUDER® AAA Endoprosthesis. This study enrolled 235 test subjects and 99 control subjects; 49 test subjects were enrolled in Continued Access. The second US clinical study, referred to as 99-04, evaluated a modified version of the GORE® EXCLUDER® AAA Endoprosthesis. These modifications were minor design changes that did not affect the risk benefit profile or clinical performance of the device. This study enrolled 193 test subjects; 88 test subjects were enrolled in Continued Access. Collectively, 565 test subjects were treated with the GORE® EXCLUDER® AAA Endoprosthesis in these studies and are described herein as the "Combined IDE Cohort". Safety outcomes of the Combined IDE Cohort are compared to the 99 control subjects treated with open surgery in the 98-03 study and are described herein as "98-03 Controls."

* Please note that the Contralateral Leg components with distal iliac diameters of 16, 18, 20, 23, and 27 mm were not evaluated under any clinical studies.

Table 2. Major Adverse Events Through Five Years Post-Treatment by Visit Interval: Combined IDE Cohort

	1 Month	6 Months	12 Months	24 Months	36 Months	48 Months	60 Months	Overall
Subjects Available at Beginning of Interval ¹	565	553	533	483	427	361	291	565
Subjects with Any Major Safety Adverse Event ²	105 (18.6%)	72 (13.0%)	91 (17.1%)	81 (16.8%)	64 (15.0%)	53 (14.7%)	22 (7.6%)	326 (57.7%)
Bleeding Complications	27 (4.8%)	2 (0.4%)	1 (0.2%)	3 (0.6%)	1 (0.2%)	5 (1.4%)	0 (0.0%)	38 (6.7%)
Pulmonary Complications	13 (2.3%)	13 (2.4%)	17 (3.2%)	10 (2.1%)	24 (5.6%)	10 (2.8%)	7 (2.4%)	81 (14.3%)
Cardiac Complications	22 (3.9%)	21 (3.8%)	24 (4.5%)	23 (4.8%)	23 (5.4%)	21 (5.8%)	10 (3.4%)	124 (21.9%)
Renal Function Complications	9 (1.6%)	6 (1.1%)	5 (0.9%)	6 (1.2%)	6 (1.4%)	1 (0.3%)	3 (1.0%)	30 (5.3%)
Wound Complications	22 (3.9%)	9 (1.6%)	3 (0.6%)	4 (0.8%)	1 (0.2%)	1 (0.3%)	0 (0.0%)	37 (6.5%)
Bowel Complications	12 (2.1%)	5 (0.9%)	10 (1.9%)	10 (2.1%)	4 (0.9%)	2 (0.6%)	2 (0.7%)	43 (7.6%)
Vascular Complications	16 (2.8%)	10 (1.8%)	9 (1.7%)	8 (1.7%)	5 (1.2%)	2 (0.6%)	2 (0.7%)	44 (7.8%)
Neurologic Complications	7 (1.2%)	5 (0.9%)	12 (2.3%)	9 (1.9%)	8 (1.9%)	3 (0.8%)	0 (0.0%)	38 (6.7%)
Other Complications	8 (1.4%)	12 (2.2%)	26 (4.9%)	19 (3.9%)	8 (1.9%)	11 (3.0%)	2 (0.7%)	77 (13.6%)
Genitourinary	6 (1.1%)	5 (0.9%)	0 (0.0%)	1 (0.2%)	0 (0.0%)	3 (0.8%)	1 (0.3%)	16 (2.8%)
Sepsis	1 (0.2%)	2 (0.4%)	0 (0.0%)	2 (0.4%)	0 (0.0%)	2 (0.6%)	0 (0.0%)	7 (1.2%)
Neoplasm	3 (0.5%)	8 (1.4%)	16 (3.0%)	10 (2.1%)	11 (2.6%)	14 (3.9%)	1 (0.3%)	61 (10.8%)
Death of Unknown Cause	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	1 (0.3%)	0 (0.0%)	1 (0.2%)
Other	8 (1.4%)	12 (2.2%)	26 (4.9%)	19 (3.9%)	8 (1.9%)	11 (3.0%)	2 (0.7%)	77 (13.6%)

1 Number of subjects with study day of last contact ≥ lower limit of specified time window. This is the denominator for all cells.
2 Entries represent number of subjects with the event (percentage). Time frames for each interval are as follows: 1 Month (0-60 days), 6 Months (61-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days), Overall (0-2006 days).

Table 3. Major Adverse Events Through Five Years Post-Treatment by Visit Interval: 98-03 Controls

	1 Month	6 Months	12 Months	24 Months	36 Months	48 Months	60 Months	Overall
Subjects Available at Beginning of Interval ¹	99	96	87	81	70	58	46	99
Subjects with Any Major Safety Adverse Event ²	68 (68.7%)	15 (15.6%)	10 (11.5%)	11 (13.6%)	12 (17.1%)	13 (22.4%)	3 (6.5%)	80 (80.8%)
Bleeding Complications	45 (45.5%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	1 (1.7%)	0 (0.0%)	45 (45.5%)
Pulmonary Complications	12 (12.1%)	3 (3.1%)	1 (1.1%)	2 (2.5%)	1 (1.4%)	4 (6.9%)	0 (0.0%)	19 (19.2%)
Cardiac Complications	16 (16.2%)	7 (7.3%)	7 (8.0%)	6 (7.4%)	6 (8.6%)	6 (10.3%)	1 (2.2%)	36 (36.4%)
Renal Function Complications	3 (3.0%)	0 (0.0%)	1 (1.1%)	0 (0.0%)	1 (1.4%)	3 (5.2%)	0 (0.0%)	7 (7.1%)
Wound Complications	4 (4.0%)	3 (3.1%)	0 (0.0%)	1 (1.2%)	2 (2.9%)	1 (1.7%)	0 (0.0%)	10 (10.1%)
Bowel Complications	17 (17.2%)	2 (2.1%)	1 (1.1%)	1 (1.2%)	2 (2.9%)	1 (1.7%)	1 (2.2%)	22 (22.2%)
Vascular Complications	7 (7.1%)	1 (1.0%)	3 (3.4%)	1 (1.2%)	0 (0.0%)	0 (0.0%)	1 (2.2%)	12 (12.1%)
Neurologic Complications	3 (3.0%)	3 (3.1%)	2 (2.3%)	0 (0.0%)	0 (0.0%)	2 (3.4%)	0 (0.0%)	9 (9.1%)
Other Complications	1 (1.0%)	2 (2.1%)	2 (2.3%)	2 (2.5%)	2 (2.9%)	1 (1.7%)	0 (0.0%)	9 (9.1%)
Genitourinary	1 (1.0%)	0 (0.0%)	2 (2.3%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	3 (3.0%)
Sepsis	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Neoplasm	0 (0.0%)	1 (1.0%)	2 (2.3%)	0 (0.0%)	3 (4.3%)	0 (0.0%)	0 (0.0%)	6 (6.1%)
Death of Unknown Cause	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Other	1 (1.0%)	2 (2.1%)	2 (2.3%)	2 (2.5%)	2 (2.9%)	1 (1.7%)	0 (0.0%)	9 (9.1%)

1 Number of subjects with study day of last contact ≥ lower limit of specified time window. This is the denominator for all cells.
2 Entries represent number of subjects with the event (percentage). Time frames for each interval are as follows: 1 Month (0-60 days), 6 Months (61-242 days), 12 Months (243-546 days), 24 Months (547-911 days), 36 Months (912-1275 days), 48 Months (1276-1640 days), 60 Months (1641-2006 days), Overall (0-2006 days).

CLINICAL STUDY OUTCOMES (04-04)

In response to the aneurysm enlargement rates associated with the original GORE® EXCLUDER® Device, Gore enhanced the design of the device by making changes to the graft material that decreased the overall permeability of the device. The low permeability GORE® EXCLUDER® AAA Endoprosthesis was approved by the FDA and released in June 2004. Gore has conducted a post-approval clinical study (04-04) of 139 subjects to evaluate the clinical performance of the low permeability GORE® EXCLUDER® AAA Endoprosthesis. These subjects are described herein as the "04-04 Test Group." Controls in this study are 120 randomly selected subjects treated with the GORE® EXCLUDER® AAA Endoprosthesis from the 98-03 clinical study.

Subjects in the 04-04 Test Group have completed a two-year follow-up (Table 4). Seven subjects died during the study. Causes of death included myocardial infarction (2), lung cancer, esophageal cancer, pneumonia (2), and sepsis / bowel perforation.

Aneurysm enlargement assessed by diameter and volume was less frequent in the 04-04 Test Group vs. Controls (Tables 5 and 6). Furthermore, aneurysm shrinkage was more frequent in the 04-04 Test Group vs. Controls. Outcomes of the 04-04 study suggest that the design enhancement of the low permeability

GORE® EXCLUDER® AAA Endoprosthesis mitigate the risk of aneurysm enlargement without endoleak.

Table 4. Subject Compliance and Disposition by Study Interval

		Follow-Up Compliance				Events Prior to Next Interval		
Study Period	Eligible for Follow-Up	Subjects with Visit in Window ¹	Site Reported CT Scan Read ¹	Site Reported X-Ray Read ²	Corelab CT Scan Read ¹	Death ¹	Discontinued ¹	Not Due for Next Follow-Up ¹
Procedure	139	139 (100.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Post-Procedure	139	139 (100.0%)	0 (0.0%)	1 (0.7%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
1 Month	139	134 (96.4%)	130 (93.5%)	60 (43.2%)	130 (93.5%)	0 (0.0%)	4 (2.9%)	0 (0.0%)
6 Months	135	115 (82.2%)	111 (82.2%)	35 (25.9%)	111 (82.2%)	3 (2.2%)	3 (2.2%)	0 (0.0%)
12 Months	129	123 (95.3%)	121 (93.8%)	30 (23.3%)	121 (93.8%)	3 (2.3%)	8 (6.2%)	0 (0.0%)
24 Months	118	115 (97.5%)	110 (93.2%)	27 (22.9%)	110 (93.2%)	0 (0.0%)	2 (1.7%)	3 (2.5%)

1 Denominators are based on number of subjects eligible for follow-up.
2 X-rays were not required as part of the 04-04 Protocol.
Time frames for each interval are as follows: Procedure: days 0-0, Post-Procedure: days 1-14, 1 Month: days 15-60, 6 Months: days 61-244, 12 Months: days 245-548 and 24 Months: days 549-912.

Table 5. Aneurysm Diameter Change Through Two Years Post-Treatment: Imaging Core Laboratory: 04-04 Test Subjects

	6 Months		12 Months		24 Months	
	Test	Control	Test	Control	Test	Control
Number of Subjects Enrolled	139	120	139	120	139	120
Subjects with Available Data ¹	103	88	114	94	104	100
≥5 mm Increase	2 (1.9%)	2 (2.3%)	0 (0.0%)	3 (3.2%)	3 (2.9%)	19 (19.0%)
No Change	76 (73.8%)	75 (85.2%)	69 (60.5%)	76 (80.9%)	50 (48.1%)	60 (60.0%)
≥5 mm Decrease	25 (24.3%)	11 (12.5%)	45 (39.5%)	15 (16.0%)	23 (46.0%)	21 (21.0%)

1 Subjects are considered to be available for evaluation if they had both a baseline (one month) and post-baseline measurement.

Table 6. Aneurysm Volume Change Through Two Years Post-Treatment: Imaging Core Laboratory: 04-04 Test Subjects

	6 Months		12 Months		24 Months	
	Test	Control	Test	Control	Test	Control
Number of Subjects Enrolled	139	120	139	120	139	120
Overall Aneurysm Change ¹						
Subjects with Available Data ²	102	86	112	92	104	101
≥5% Increase	6 (5.9%)	20 (23.3%)	9 (8.0%)	31 (33.7%)	17 (16.3%)	47 (46.5%)
No Change	29 (28.4%)	32 (37.2%)	26 (23.2%)	26 (28.3%)	15 (14.4%)	23 (22.8%)
≥5% Decrease	67 (65.7%)	34 (39.5%)	77 (68.8%)	35 (38.0%)	72 (69.2%)	31 (30.7%)

1 Subjects are considered to be available for evaluation if they had both a baseline (one month) and post-baseline measurement.
2 Overall volume increase is computed as follows: if either distal renal to aortic bifurcation or right hypogastric to aortic bifurcation volume measurement increases ≥ 5% then the subject is considered an increaser. If either measure decreases ≥ 5% then the subject is considered a decreaser. Other non-missing combinations are considered No Change.

CLINICAL STUDY OUTCOMES (03-02)

The GORE® EXCLUDER® Device product line has been expanded to include a 31 mm Trunk-Ipsilateral Leg and 32 mm Aortic Extender to treat patients with larger proximal aortic neck diameters. The 03-02 study was initiated to evaluate the 31 mm GORE® EXCLUDER® Device compared to open surgical repair and to compare device performance to the Original GORE® EXCLUDER® Device in the treatment of AAAs with large proximal neck diameters. Enrollment of Test patients into this study was completed in January 2008 with a total of 35 patients enrolled. The larger size GORE® EXCLUDER® Devices were approved by the FDA in March 2009.

Potential Device or Procedure Related Adverse Events

Adverse events that may occur and / or require intervention or additional intraoperative procedure time include, but are not limited to:

- allergic reaction and / or anaphylactoid response to x-ray contrast dye, anti-platelet therapy, device materials
- amputation
- anesthetic complications
- aneurysm enlargement
- aneurysm rupture and death
- arterial or venous thrombosis and / or pseudoaneurysm
- arteriovenous fistula
- bleeding, hematoma, or coagulopathy
- bowel (e.g., ileus, transient ischemia, infarction, necrosis)
- cardiac (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension or hypertension)
- claudication (e.g., buttock, lower limb)
- death
- dissection, perforation, or ruptures of the aortic vessel and surrounding vasculature
- edema
- embolization (micro and macro) with transient or permanent ischemia
- endoleak
- endoprosthesis or delivery system: improper component placement; incomplete component deployment; unintentional/premature component deployment; leading end catheter component retention; component migration; separation of graft material from stent; occlusion; infection; stent fracture; graft material failure, dilatation, erosion, puncture, perigraft flow
- fever and localized inflammation
- genitourinary (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- hepatic failure
- impotence
- infection (e.g., aneurysm, device or access sites)
- lymph fistula / complications
- multi-system organ failure
- neurologic damage, local or systemic (e.g., stroke, paraplegia, paraparesis)
- occlusion of device or native vessel
- post-implant syndrome
- pulmonary complications (e.g., pneumonia, respiratory failure)
- radiation injury, late malignancy
- renal (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- surgical cut down, bypass or conversion
- tissue necrosis
- wound (e.g., infection, dehiscence)
- vascular spasm or vascular trauma (e.g., ilio-femoral vessel dissection, bleeding, rupture, death)

Device Related Adverse Event Reporting

Any adverse event involving the GORE® EXCLUDER® AAA Endoprosthesis should be reported to W. L. Gore & Associates immediately. To report an event in the US, call 800.437.8181. Outside the US, contact your local technical representative.

PATIENT SELECTION AND TREATMENT

(SEE WARNINGS AND PRECAUTIONS)

Individualization of Treatment

Gore recommends that the GORE® EXCLUDER® AAA Endoprosthesis diameter be at least 2 mm larger than the aortic inner diameter (10 – 21% oversizing) and 1 mm larger than the iliac inner diameter (7 – 25% oversizing) as described in Tables 7 and 8. The length of the GORE® EXCLUDER® AAA Endoprosthesis should be sufficient to reach from just inferior to the most distal (lowest) major renal artery to non-aneurysmal tissue in the common or external iliac arteries. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters / lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

Additional considerations for patient selection include but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- Patient's anatomical suitability for endovascular repair
- The risk of aneurysm rupture compared to the risk of treatment with the GORE® EXCLUDER® AAA Endoprosthesis
- Ability to tolerate general, regional, or local anesthesia
- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and / or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 12 Fr - 18 Fr vascular introducer sheath
- An infrarenal non-aneurysmal aortic neck length of at least 15 mm and a diameter of no greater than 32 mm
- Proximal neck angulation $\leq 60^\circ$ with minimal thrombus or calcification

- Distal segment iliac vessel lengths of at least 30 mm of which at least 10 mm must be less than or equal to 25 mm in diameter
- Freedom from significant femoral / iliac artery occlusive disease that would impede inflow or outflow of stent-grafts

The final treatment decision is at the discretion of the physician and patient.

PATIENT COUNSELING INFORMATION

The physician and patient should review the risks and benefits when discussing this endovascular device and procedure including:

- Risks and differences between endovascular repair and open surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to post-operative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term safety and effectiveness of endovascular repair has not been established.** Physicians should advise all patients that this treatment modality requires long-term, regular follow-up to assess patients' health status and stent-graft performance. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should receive enhanced follow-up. Patients should be counseled on the need for regular follow-up, even in the absence of obvious symptoms, e.g., pain, numbness, weakness (see IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP).
- Regular follow-up including imaging of the device should be performed at least every 12 months for all patients and at least every 6 to 12 months for patients with known endoleaks or aneurysm enlargement for the duration of the implant (see IMAGING GUIDELINES AND POST- OPERATIVE FOLLOW-UP).
- Physicians must advise all patients that it is important to seek prompt medical attention if he / she experiences signs of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking, or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as pain, numbness, weakness in the legs; any back, chest, abdominal, or groin pain, dizziness, fainting, rapid heartbeat, or sudden weakness.

US physicians are encouraged to refer the patient to the Patient Brochure regarding risks occurring during or after implantation of the device. Procedure related risks include cardiac, pulmonary, neurologic, bowel, and bleeding complications. Device related risks include occlusion, endoleak, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture and death (See POTENTIAL DEVICE OR PROCEDURE RELATED ADVERSE EVENTS).

US physicians are encouraged to complete the Patient Wallet Card and give it to the patient so that he / she can carry it with them at all times. The patient should refer to the wallet card any time they visit additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

HOW SUPPLIED

The GORE® EXCLUDER® AAA Endoprosthesis is preloaded on a delivery catheter and supplied sterile and non-pyrogenic.

Storage and Handling

- Do not resterilize; for single use only.
- Do not use if damaged or if sterile barrier has been compromised.
- Do not use after the "use by" (expiration) date printed on the label.
- Store in a cool, dry place.

CLINICAL USE INFORMATION

Physician Training Program

CAUTION: Always have a vascular surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The GORE® EXCLUDER® AAA Endoprosthesis should only be used by physicians experienced in vascular interventional techniques, and who have successfully completed the appropriate physician training program.

The recommended skill / knowledge requirements for physicians using the GORE® EXCLUDER® AAA Endoprosthesis are outlined below:

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair
- Knowledge of radiographic image interpretation, device selection and sizing

A multi-disciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy, and repair
- Percutaneous access and closure techniques
- Non-selective and selective guidewire and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty

- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

Recommended Materials

- 0.035" (0.89 mm) 'super stiff' guidewire, 145 cm or longer
- Angiographic radiopaque marker catheter
- Contrast media
- Snare Catheter
- Syringe
- Heparin and heparinized saline
- **Trunk-Ipsilateral Leg Endoprosthesis and Contralateral Leg Endoprosthesis:**
 - 18 Fr and 12 Fr, 14 Fr or 15 Fr introducer sheaths (Tables 7 and 8)
 - Large diameter, low pressure aortic balloon (monitor balloon volumes and pressures as recommended in balloon catheter Instructions for Use)
 - Percutaneous transluminal angioplasty (PTA) balloons (Table 8)
- **Aortic Extender Endoprosthesis:**
 - 16 Fr, 17 Fr or 18 Fr introducer sheath
 - Large diameter, low pressure aortic balloon (monitor balloon volumes and pressures as recommended in balloon catheter Instructions for Use)
- **Iliac Extender Endoprosthesis:**
 - 12 Fr, 14 Fr and 15 Fr introducer sheaths
 - PTA balloon catheters, 10 mm x 40 mm, 12 mm x 40 mm, 14 mm x 40 mm, 16 mm x 40 mm, 18 mm x 40 mm, 20 mm x 40 mm, 24 mm x 40 mm and 28 mm x 40 mm (Table 10)

Table 7. Trunk-Ipsilateral Leg Endoprosthesis Sizing Guide*

Intended Aortic Vessel Diameter (mm)	Aortic Endoprosthesis Diameter ¹ (mm)	Intended Iliac Vessel Diameter (mm)	Iliac Endoprosthesis Diameter ² (mm)	Overall Device Lengths (cm)	Recommended Introducer Sheath ³ (Fr)
19 - 21	23	10 - 11	12	12, 14, 16, 18	18
		12 - 13.5	14.5		
22 - 23	26	10 - 11	12	12, 14, 16, 18	18
		12 - 13.5	14.5		
24 - 26	28.5	10 - 11	12	12, 14, 16, 18	18
		12 - 13.5	14.5		
27 - 29	31	12 - 13.5	14.5	13, 15, 17	18
30 - 32	35	12 - 13.5	14.5	14, 16, 18	18

1 Recommended endoprosthesis oversizing relative to the aortic vessel is approximately 10-21%, and for the iliac vessel approximately 7-25%.

2 Recommended angioplasty balloon size is 12 mm and 14 mm respectively.

3 GORE® Introducer Sheaths are recommended. The GORE® EXCLUDER® AAA Endoprosthesis is not compatible with introducer sheaths longer than 40 cm (total length including the hemostatic valve).

* Note: All dimensions are nominal.

Table 8. Contralateral Leg Endoprosthesis Sizing Guide*

Intended Iliac Vessel Diameter (mm)	Iliac Endoprosthesis Diameter ¹ (mm)	Overall Device Lengths ^{2,3} (cm)	Recommended Contralateral Introducer Sheath ⁴ (Fr)	Recommended Angioplasty Balloon Size (mm x mm)
10 - 11	12	10, 12, 14	12	12 x 40
12 - 13.5	14.5	10, 12, 14	12	14 x 40
13.5 - 14.5	16	9.5, 11.5, 13.5	12	16 x 40
14.5 - 16.5	18	9.5, 11.5, 13.5	12	18 x 40
16.5 - 18.5	20	9.5, 11.5, 13.5	12	20 x 40
18.5 - 21.5	23	10, 12, 14	14	24 x 40
21.5 - 25	27	10, 12, 14	15	28 x 40

1 Recommended endoprosthesis oversizing relative to the vessel is approximately 7-25%.

2 Total treatable lengths include 4 cm of the trunk region of the Trunk-Ipsilateral Leg Endoprosthesis.

3 Labeled Contralateral Leg length includes 3 cm overlap.

4 GORE® Introducer Sheaths are recommended. The GORE® EXCLUDER® AAA Endoprosthesis is not compatible with introducer sheaths longer than 40 cm (total length including the hemostatic valve).

* Note: All dimensions are nominal.

Table 9. Aortic Extender Endoprosthesis Sizing Guide*

Intended Aortic Vessel Diameter (mm)	Aortic Extender Diameter ¹ (mm)	Endoprosthesis Length (cm)	Recommended Introducer Sheath ² (Fr)
19 - 21	23	3.3	16
22 - 23	26	3.3	16
24 - 26	28.5	3.3	16
27 - 29	32	4.5	17
30 - 32	36	4.5	18

1 Recommended endoprosthesis oversizing relative to the vessel diameter is approximately 10-21%.

2 GORE® Introducer Sheaths are recommended. The GORE® EXCLUDER® AAA Endoprosthesis is not compatible with introducer sheaths longer than 40 cm (total length including the hemostatic valve).

* Note: All dimensions are nominal.

Table 10. Iliac Extender Endoprosthesis Sizing Guide*

Intended Iliac Vessel Diameter (mm)	Distal Iliac Extender Diameter ¹ (mm)	Endoprosthesis Length ² (cm)	Recommended Introducer Sheath ⁴ (Fr)	Recommended Balloon Size (Proximal) (mm)	Recommended Angioplasty Balloon Size (Distal) (mm x mm)
8 - 9	10	7	12	14	10 x 40
10 - 11	12	7	12	14	12 x 40
12 - 13.5	14.5	7	12	14	14 x 40
13.5 - 14.5	16 ^{3,5}	9.5, 11.5, 13.5	12	14	16 x 40
14.5 - 16.5	18 ^{3,5}	9.5, 11.5, 13.5	12	14	18 x 40
16.5 - 18.5	20 ^{3,5}	9.5, 11.5, 13.5	12	14	20 x 40
18.5 - 21.5	23 ^{3,5}	10, 12, 14	14	14	24 x 40
21.5 - 25	27 ^{3,5}	10, 12, 14	15	14	28 x 40

1 Recommended endoprosthesis oversizing relative to the vessel diameter is approximately 7-25%.

2 7 cm long Iliac Extender Endoprosthesis provides a maximum extension of 4 cm when placed in the Trunk-Ipsilateral or Contralateral Leg Endoprosthesis; labeled length includes 3 cm overlap.

3 16, 18, 20, 23, and 27 mm Contralateral Legs may be used as Iliac Extenders.

4 GORE® Introducer Sheaths are recommended. The GORE® EXCLUDER® AAA Endoprosthesis is not compatible with introducer sheaths longer than 40 cm (total length including the hemostatic valve).

5 If extension of a 16, 18, 20, 23, or 27 mm Contralateral Leg is necessary, a 16 mm angioplasty balloon size is recommended for the proximal end of the extension.

* Note: All dimensions are nominal.

DIRECTIONS FOR USE

Pre-Treatment Planning

- Determine accurate size of anatomy and proper size of Trunk-Ipsilateral and Contralateral Endoprostheses (Tables 7 and 8) and Aortic and Iliac Extender Endoprostheses (Tables 9 and 10).
- Use high resolution, non-contrast and contrast enhanced computerized tomography (CT / CTA) at ≤ 3 mm acquisition and reconstruction collimation.
- Use multiple view, digital subtraction angiography with a radiopaque marker catheter or spiral CT multi-planar reconstruction.
- For angiography, use correct imaging angulation (cranial-caudal, lateral-oblique) to accurately identify origin of branch vessel anatomy.
- Consider breath-hold technique to optimize digital subtraction angiography image quality.

Anatomical Requirements

- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and / or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 12 Fr - 18 Fr vascular introducer sheath
- An infrarenal, non-aneurysmal aortic neck length of at least 15 mm and an infrarenal aortic neck treatment diameter range of 19-32 mm.
- For Trunk-Ipsilateral Leg Endoprosthesis and Aortic Extender Endoprosthesis: Proximal aortic neck angulation ≤ 60° with minimal thrombus and / or calcification.
- Key anatomic elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation, short proximal aortic neck and significant thrombus and / or calcium at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. The US clinical studies quantify significant thrombus as thrombus ≥ 2 mm in thickness and / or ≥ 25% of the vessel circumference in the intended seal zone of the aortic neck. Irregular calcium and / or plaque may compromise the fixation and sealing of the implantation sites.
- Distal segment iliac vessel lengths of at least 30 mm of which at least 10 mm must be less than or equal to 25 mm in diameter for Iliac Extender Endoprosthesis: Non-aneurysmal iliac artery length ≥ 10 mm of appropriate diameter.
- Freedom from significant femoral / iliac artery occlusive disease that would impede inflow or outflow of stent-grafts.
- Ability to tolerate general, regional, or local anesthesia.
- Patient's anatomical suitability for endovascular repair.

Arterial Access and Angiography

1. Following standard practices, access the intended contralateral side via a percutaneous diagnostic sheath, and perform marker catheter digital subtraction angiography (AP, oblique and lateral views as necessary) to confirm the correct device component sizing, and deployment locations. Consider breath-hold technique to optimize image quality. Leave marker catheter in place at the level of the renal arteries.
2. Following standard practices, perform percutaneous access and / or surgical exposure of the vessels selected to receive the Trunk-Ipsilateral and Contralateral side introducer sheaths.
3. Following the manufacturer's instructions for use, advance a 0.035" (0.89 mm) 'super stiff' guidewire, or acceptable equivalent to the level of the renal arteries.
4. Following the manufacturer's instructions for use, prepare and advance the recommended introducer sheath (Tables 7 – 10) over the guidewire, through the ilio-femoral anatomy, aortic aneurysm and up to the level of the proximal aortic neck according to standard practice.
5. **CAUTION:** Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
6. Use standard heparinized saline, pressure flush system technique to prevent thrombus formation in the introducer sheaths.
7. Use an accurate radiopaque patient marking method to assure accurate device positioning and deployment locations.

Catheter Preparation

1. Use new, sterile gloves when preparing device.

CAUTION: Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.

2. Remove the appropriately sized device delivery catheters from their packaging and examine for possible damage.
3. Remove protective packaging mandrel and packaging sheath(s) from the leading end of the delivery catheters (Figures 3A – 3F).
4. Flush with heparinized saline through the flushing port on the trailing end of the delivery catheter (Figures 3A – 3F).
5. Follow the manufacturer's recommended method for size selection, preparation and use of aortic and iliac dilation balloons. Carefully inflate the balloon to avoid complications.

Trunk-Ipsilateral Leg Endoprosthesis Positioning

1. Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
2. Advance the Trunk delivery catheter over a 0.035" (0.89 mm) 'super stiff' guidewire, through the 18 Fr introducer sheath into the aorta to the approximate approximate level of intended positioning.

WARNING: Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

WARNING: Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

3. While maintaining the delivery catheter in position, withdraw the introducer sheath (Table 7) and visually verify that the leading end of the introducer sheath is below the distal Trunk-Ipsilateral Leg radiopaque marker.
4. Magnify and center the fluoroscopic image on the proximal trunk. Reposition and rotate the Trunk-Ipsilateral delivery catheter as necessary to properly position the proximal device marker as well as orient the long contralateral, and short ipsilateral radiopaque markers and device position on the appropriate side of the anatomy. Maximize the separation between these two markers to achieve maximum lateral positioning of the iliac legs of the device. The long marker should be oriented toward the contralateral side (Figures 1A and 1B).
WARNING: Do not rotate the Trunk delivery catheter beyond 360°. Delivery system damage and / or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
5. It is recommended to view and confirm the distal position of the iliac end of the device relative to the internal iliac artery to ensure accurate and desired deployment position of the distal aspect of the device.
6. If clinically acceptable, lower the patient's mean arterial pressure to 60 – 70 mm Hg during Trunk deployment and aortic balloon inflation to decrease blood flow and reduce the risk of endoprosthesis movement.
7. Maintain a contralateral access side sheath, catheter or guidewire in position across the distal, native bifurcation to protect and ensure that contralateral access is maintained into the aneurysm sac and contralateral leg hole of the device during Trunk-Ipsilateral component deployment.
8. Re-center and magnify the image on the proximal Trunk of the device to assure final desired position of proximal device relative to anatomy. Stabilize the Trunk delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.
WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
9. When using the Original Delivery System (Figure 3A), please refer to the deployment instructions entitled, ORIGINAL DELIVERY SYSTEM - TRUNK IPSILATERAL ENDOPROSTHESIS DEPLOYMENT.
10. When using the GORE® C3® Delivery System (Figure 3B), please refer to the deployment instructions entitled, GORE® C3® DELIVERY SYSTEM - TRUNK-IPSILATERAL ENDOPROSTHESIS DEPLOYMENT.

ORIGINAL DELIVERY SYSTEM - TRUNK IPSILATERAL ENDOPROSTHESIS DEPLOYMENT

1. Loosen the deployment knob. Confirm final device position and orientation and deploy the Trunk using a steady and continuous pull of the deployment knob to release the endoprosthesis. Pull the deployment knob straight out from the catheter handle. Deployment initiates from the leading end toward the trailing end.

WARNING: Do not attempt to reposition the endoprosthesis after deployment has been initiated. Vessel damage or device misplacement may result.

WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.

CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.

2. Use fluoroscopic guidance during the withdrawal of the delivery catheter to assure safe removal from, and to avoid catching on, the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop and withdraw delivery catheter and introducer sheath together.

WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.

3. When withdrawing the device delivery system through the introducer sheath, verify all catheter components are intact.


WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred, see ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).

4. Position the aortic balloon inside the proximal region of the trunk. Avoid balloon contact with the flow splitter which is aligned with the long and short radiopaque markers. Inflate and deflate the balloon quickly with dilute contrast solution to seat the aortic end of the endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of aortic and iliac dilation balloons, carefully monitoring both volume and pressure to avoid complications.

5. Use fluoroscopic guidance to ensure the balloon is completely deflated and is safely removed from the endoprosthesis.

6. Advance and inflate the appropriate size PTA balloon catheter to seat the iliac end of the endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of PTA balloons. Carefully inflate the balloon to avoid complication.

GORE® C3® DELIVERY SYSTEM - TRUNK IPSILATERAL ENDOPROSTHESIS DEPLOYMENT

1. Loosen the White Outer Deployment Knob by turning it 90° counter-clockwise (Figure 8A). Confirm final device position and orientation and deploy the Trunk using a steady and continuous pull of the deployment knob to release the proximal endoprosthesis past the contralateral leg hole of the Trunk. Pull the deployment knob straight out from the catheter handle. Deployment initiates from the leading end toward the trailing end. The ipsilateral leg will remain constrained on the delivery catheter. If the Trunk is not properly positioned following removal of the white deployment knob, continue with Optional Repositioning of Trunk-Ipsilateral Leg Endoprosthesis, described below.

WARNING: Do not attempt to reposition the endoprosthesis after complete deployment of the device. Vessel damage or device misplacement may result.

WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.

WARNING: Do not attempt to remove a partially deployed Trunk component.

CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.

BACKUP DEPLOYMENT MECHANISM: In the event that the removal of the White Outer Deployment Knob does not initiate the deployment of the proximal end of the Trunk device, and the problem is suspected to have occurred in the deployment handle, the deployment handle has a design feature that is intended to protect the patient against this problem becoming a hazard. Remove the Deployment Line Access Hatch Lid (Figure 9) and then pull the first deployment line (labeled "1") with an appropriate dull or blunt-edged tool. Fully remove the deployment line in a steady, continuous motion. This action will cause the proximal Trunk to deploy past the contralateral leg hole.

WARNING: Do not cut the first deployment line (labeled "1"). Patient injury, catheter damage, and / or device damage may occur.

Optional Repositioning of Trunk-Ipsilateral Leg Endoprosthesis

- A. To constrain the proximal end of the Trunk, turn the Gray Constraining Dial clockwise until it stops (Figure 8B). The Black Nut within the Transparent Knob housing will move toward the trailing end of the catheter when the Trunk end of the device is constrained. Verify the distance between the three proximal radiopaque markers is reduced (Figures 1A and 1B).
- B. Magnify and center the fluoroscopic image on the proximal Trunk. Manually reposition and / or rotate the Trunk-Ipsilateral component delivery catheter as necessary to properly position the Trunk device. Maximize the separation between the long and short radiopaque markers to achieve maximum lateral iliac leg positioning of the device. The long marker should be oriented toward the contralateral side (Figures 1A and 1B).
- C. Repeat steps 4 through 8 listed above to properly position the endoprosthesis while the Trunk end of the endoprosthesis is constrained.
- D. When the Trunk is properly positioned, turn the Gray Constraining Dial counter-clockwise until it stops to reopen the proximal end of the Trunk (Figure 8C). The Black Nut within the Transparent Knob housing will move toward the leading end of the catheter when the proximal end of the Trunk is reopened. Verify that the distance between the three proximal radiopaque markers is maximized (Figures 1A and 1B).

WARNING: The constraining mechanism cannot be removed unless the proximal end of the Trunk device is fully opened.

WARNING: Do not constrain / reopen the Trunk device more than two times during a procedure. Device and / or catheter damage may occur.

WARNING: In the event that the complete counter-clockwise rotation of the Gray Constraining Dial does not initiate the reopening of the proximal end of the Trunk device, and the problem is suspected to have occurred in the deployment handle, the deployment handle has a design feature that is intended to protect the patient against this problem becoming a hazard. Verify the final desired Trunk device position relative to anatomy. Stabilize the Trunk delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site. Disengage the constraining mechanism by sliding the Red Safety Lock back while rotating the entire Transparent Knob 90° counter-clockwise (Figure 8D). Remove the constraining mechanism by pulling the Transparent Knob straight out from the catheter handle (Figure 8E). These actions will cause the constraining mechanism to be removed and the proximal end of the Trunk device to reopen.

Contralateral Leg Hole Guidewire Cannulation

1. Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
2. Following manufacturer's instructions for use, advance a 0.035" (0.89 mm) 'super stiff' guidewire into the contralateral leg hole of the Trunk according to standard practice.
3. Verify that the guidewire is within the contralateral leg hole of the Trunk by standard practice used to verify guidewire location. If multiple cannulation attempts fail, the Trunk-Ipsilateral Leg Endoprosthesis may be repositioned for better contralateral gate access by following steps A through D above.
WARNING: Do not constrain / reopen the Trunk device more than two times during a procedure. Device and / or catheter damage may occur.
4. Following manufacturer's instructions for use, introduce the recommended introducer sheath (Table 8). Advance the sheath over the guidewire and through the contralateral leg hole of the Trunk.

Ipsilateral Leg Deployment

1. Verify the proximal end of the Trunk is fully open by turning the Gray Constraining Dial counter-clockwise until it stops (Figure 8C). The Black Nut within the Transparent Knob housing will move toward the leading end of the catheter when the proximal end of the Trunk is opened. Verify that the distance between the three proximal radiopaque markers is maximized (Figures 1A and 1B).
2. Disengage the constraining mechanism by sliding the Red Safety Lock back while rotating the entire Transparent Knob 90° counter-clockwise (Figure 8D). Remove the constraining mechanism by pulling the Transparent Knob straight out from the catheter handle (Figure 8E).

WARNING: The constraining mechanism cannot be removed unless the proximal end of the Trunk device is fully opened. Verify the Gray Constraining Dial is rotated counter-clockwise until it stops, prior to removing the constraining mechanism. Attempting to remove the constraining mechanism while the Trunk is constrained may result in misplacement of the device, patient injury, and / or device damage.

BACKUP DEPLOYMENT MECHANISM: In the event, that the removal of the Transparent Knob does not initiate the removal of the constraining mechanism and the problem is suspected to have occurred in the deployment handle, the deployment handle has a design feature that is intended to protect the patient against this problem becoming a hazard. Remove the Deployment Line Access Hatch Lid (Figure 9) and then cut the silver Lock Pin line (labeled "2") with an appropriate tool. Grasp, pull, and fully remove the silver Lock Pin line (labeled "2") in a steady, continuous motion. Cut, grasp, pull, and fully remove the gray Constraining Loop line (labeled "3") in a steady, continuous motion. These actions will cause the constraining mechanism to be removed.

3. Loosen the gray deployment knob (Figure 8F). Deploy the Ipsilateral Leg of the Trunk device using a steady and continuous pull of the gray deployment knob to release the endoprosthesis. Pull the deployment line straight out from the catheter handle.

BACKUP DEPLOYMENT MECHANISM: In the event that the removal of the gray deployment knob does not initiate the deployment of the ipsilateral leg and the problem is suspected to have occurred in the deployment handle, the deployment handle has a design feature that is intended to protect the patient against this problem becoming a hazard. Remove the Deployment Line Access Hatch Lid (Figure 9) and then pull the second deployment line (labeled "4") with an appropriate dull or blunt-edged tool. Fully remove the deployment line in a steady, continuous motion. This action will cause the ipsilateral leg to deploy.

WARNING: Do not cut the second deployment line (labeled "4"). Patient injury, catheter damage, and / or device damage may occur.

4. Use fluoroscopic guidance during the withdrawal of the delivery catheter to assure safe removal from, and to avoid catching on, the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop, visually assess the cause of the resistance, and withdraw delivery catheter and introducer sheath together if appropriate.

WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.

5. When withdrawing the device delivery system through the introduce sheath, verify all catheter components are intact.

WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred, see ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).

6. Position the aortic balloon inside the proximal region of the trunk. Avoid balloon contact with the flow splitter which is aligned with the long and short radiopaque markers. Inflate and deflate the balloon quickly with dilute contrast solution to seat the aortic end of the endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of aortic and iliac dilation balloons, carefully monitoring both volume and pressure to avoid complications.
7. Use fluoroscopic guidance to ensure the balloon is completely deflated and is safely removed from the endoprosthesis.
8. Advance and inflate the appropriate size PTA balloon catheter to seat the iliac end of the endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of PTA balloons. Carefully inflate the balloon to avoid complication.

Contralateral Leg Endoprosthesis Positioning and Deployment

1. Advance the prepped Contralateral Leg Endoprosthesis delivery catheter to the level of the long radiopaque marker (Figure 1).

WARNING: Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

WARNING: Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment may occur have occurred and may result in potential patient harms, see ADVERSE EVENTS.

2. Align the radiopaque marker at the proximal end of the Contralateral Leg Endoprosthesis with the long contralateral radiopaque marker on the Trunk-Ipsilateral Leg Endoprosthesis. With the alignment of these markers, an approximate 3 cm overlap will be achieved.
3. While maintaining the delivery catheter in position, withdraw the introducer sheath (Table 8) and visually verify that the leading end of the introducer sheath is below the distal Contralateral Leg radiopaque marker.

WARNING: Do not rotate the Contralateral Leg delivery catheter during delivery, positioning or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

4. Stabilize the Contralateral Leg Endoprosthesis delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.

WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.

5. Loosen the deployment knob. Confirm final device position. Deploy the Contralateral Leg Endoprosthesis by using a steady, continuous pull of the deployment knob to release the endoprosthesis. Pull the deployment knob straight out from the catheter side-arm. Deployment initiates from the leading (aortic) end toward the trailing (iliac) end.

WARNING: Do not attempt to reposition the endoprosthesis after deployment has been initiated. Vessel damage or device misplacement may result.

WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.

CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.

6. Use fluoroscopic guidance during withdrawal of the delivery catheter to assure safe removal from the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop and withdraw delivery catheter and introducer sheath together.

WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.

7. When withdrawing the device delivery system through the introduce sheath, verify all catheter components are intact.

WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred, see ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).

8. Following manufacturer's instructions for use, advance and inflate a 14 mm PTA balloon catheter to seat the proximal end of the Contralateral Leg Endoprosthesis within the contralateral leg hole overlap region. Follow the manufacturer's recommended method for size selection, preparation and use of aortic and iliac dilation balloons, carefully monitoring both volume and pressure to avoid complications.
9. Following manufacturer's instructions for use, advance and inflate the appropriate size PTA balloon to seat the iliac end of the Contralateral Leg Endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of PTA balloons. Carefully inflate the balloon to avoid complications.

Aortic Extender Endoprosthesis Positioning and Deployment

1. Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
2. Advance the Aortic Extender Endoprosthesis delivery catheter over a 0.035" (0.89 mm) 'super stiff' guidewire, through the 16 Fr, 17 Fr or 18 Fr introducer sheath into the aorta, just proximal to the level of intended device positioning.
WARNING: Do not attempt to advance the Aortic Extender through the 12 Fr introducer sheath. The Aortic Extender is designed for a 16 Fr, 17 Fr, or 18 Fr sheath.
WARNING: Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
3. While maintaining the delivery catheter in position, withdraw the introducer sheath (Table 9) and visually verify that the leading end of the introducer sheath is below the distal Aortic Extender radiopaque marker.
WARNING: Do not rotate the Aortic Extender delivery catheter during delivery, positioning or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
4. Magnify and center the fluoroscopic image on the proximal Aortic Extender Endoprosthesis. Reposition the Aortic Endoprosthesis delivery catheter as necessary to position the proximal and distal radiopaque markers in appropriate position. The maximum recommended extension with each Aortic Extender component is approximately one-half of the Extender length inside (1.6 cm or 2.2 cm) and one-half, outside (1.6 cm or 2.2 cm), or proximal to the Trunk or Aortic Extender host component. The proximal three (3) and distal one (1) markers are visible relative to host device and anatomy pre and post-deployment (Figures 4 and 5B).
5. If clinically acceptable, lower the patient's mean arterial pressure to 60 – 70 mm Hg during Aortic Extender deployment and aortic balloon inflation to decrease blood flow and reduce the risk of endoprosthesis movement.
6. Stabilize the Extender delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.
WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
7. Loosen the deployment knob. Using fluoroscopy, confirm final device position and deploy the Aortic Extender using a steady and continuous pull of the deployment knob to release the endoprosthesis. Pull the deployment knob straight out from the catheter side-arm. Deployment initiates from the trailing end of the device toward the leading end of the device.
WARNING: Do not attempt to reposition the endoprosthesis after deployment has been initiated. Vessel damage or device misplacement may result.
WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.
8. Use fluoroscopic guidance during the withdrawal of the delivery catheter to assure safe removal from, and to avoid catching on, the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop and withdraw delivery catheter and introducer sheath together.
WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.
9. When withdrawing the device delivery system through the introduce sheath, verify all catheter components are intact.
WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred, see ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through

both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).

10. Advance the aortic dilation balloon until it is centered relative to the endoprosthesis. Inflate and deflate the balloon quickly with dilute contrast solution to seat the Aortic Extender Endoprosthesis. Follow the manufacturer's recommended method for size selection, preparation and use of aortic dilation balloons. Carefully inflate the balloon to avoid complications.
11. Use fluoroscopic guidance to ensure the balloon is completely deflated and to assure safe removal from, and to avoid catching on, the endoprosthesis.

Iliac Extender Endoprosthesis Positioning and Deployment

1. Use fluoroscopic visualization for all guidewire, sheath and device catheter manipulations.
2. Advance the Iliac Extender Endoprosthesis delivery catheter into the distal end of the host device, via the recommended introducer sheath (Table 10).
WARNING: Do not advance the device outside of the sheath while tracking it into position. Catheter breakage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not rotate the device delivery catheter while the endoprosthesis is inside the introducer sheath. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
WARNING: Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the guidewire, sheath, or catheter. Stop and assess the cause of resistance. Vessel or catheter damage or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
3. For maximum extension, align the radiopaque marker at the iliac (distal) end of the host device with the marker located 3 cm below the proximal end of the Extender component (Figures 2B, 3D, 6A and 6B).
WARNING: While using 16, 18, 20, 23, or 27 mm Contralateral Legs as an Iliac Extender, ensure that the distal end including the taper zone will not be deployed inside the previously deployed Ipsilateral Leg or Contralateral Leg of the GORE® EXCLUDER® AAA Endoprosthesis. However, when the Contralateral Leg and Iliac Extender diameters are identical, the taper zone can be deployed inside the previously deployed Contralateral Leg (Figure 7B).
WARNING: While using 16, 18, 20, 23, or 27 mm Contralateral Legs as an Iliac Extender, the 3 cm mandatory overlap must be achieved prior to the beginning of the distal taper zone of the 18, 20, 23, and 27 mm Contralateral Leg. Inadequate sealing may lead to endoleak.
4. While maintaining the delivery catheter in position, withdraw the introducer sheath (Table 10) and visually verify that the leading end of the introducer sheath is below the distal Iliac Extender radiopaque marker.
WARNING: Do not rotate the Iliac Extender delivery catheter during delivery, positioning or deployment. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
5. Stabilize the Iliac Extender delivery catheter at the level of entry into the introducer sheath and stabilize the sheath relative to the patient's access site.
WARNING: Do not attempt to withdraw any undeployed endoprosthesis through the introducer sheath. The sheath and catheter must be removed together. Catheter breakage or separation or premature deployment have occurred and may result in potential patient harms, see ADVERSE EVENTS.
6. Loosen the deployment knob. Confirm final device position. Using fluoroscopy, deploy the Iliac Extender Endoprosthesis by using a steady, continuous pull of the deployment knob to release the endoprosthesis. Pull the deployment knob straight out from the catheter side-arm. The device deploys from the leading (proximal) end toward the trailing (distal) end.
WARNING: Do not attempt to reposition the endoprosthesis after deployment has been initiated. Vessel damage or device misplacement may result.
WARNING: Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
CAUTION: Do not cover significant renal or mesenteric arteries with the endoprosthesis. Vessel occlusion may occur. During the US clinical studies, this device was not studied in patients with two occluded internal iliac arteries.
7. Use fluoroscopic guidance during the withdrawal of the delivery catheter to assure safe removal from, and avoid catching on, the endoprosthesis. If resistance is felt during removal of delivery catheter through the introducer sheath, stop and withdraw delivery catheter and introducer sheath together.
WARNING: Do not continue to withdraw the delivery catheter if resistance is felt during removal through the introducer sheath. Forcibly withdrawing the delivery catheter through the introducer sheath when resistance is encountered has resulted in adverse events including catheter breakage or separation resulting in potential patient harms, see ADVERSE EVENTS.
8. When withdrawing the device delivery system through the introduce sheath, verify all catheter components are intact.
WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal).
9. Advance and inflate an appropriate size PTA balloon catheter to seat the proximal overlap end and the distal end of the Iliac Extender Endoprosthesis. Follow the

manufacturer's recommended method for size selection, preparation and use of PTA balloons. Carefully inflate the balloon to avoid complications.

10. Use fluoroscopic guidance to ensure the balloon is completely deflated and to assure safe removal from, and to avoid catching on, the endoprosthesis.

Completion of the Procedure

1. Perform extended imaging angiography to confirm exclusion of the aneurysm. Consider breath-hold technique to optimize digital subtraction angiography image quality. Consider use of GORE® EXCLUDER® AAA Endoprosthesis Extender components as necessary. For Aortic Extenders, a minimum overlap of 1.6 cm or 2.2 cm is required, offering a maximum of 1.6 cm or 2.2 cm of extension; for Iliac Extenders, a minimum overlap of 3 cm is required.
2. Prior to removing wires and sheaths verify all catheter components have been withdrawn from the patient.
WARNING: Catheter leading end separation or breakage and related potential patient harms have occurred. See ADVERSE EVENTS. If catheter separation occurs, use best medical judgment to determine the appropriate course of action for the patient. Effective removal of the catheter component has been reported through both surgical (e.g. cut down) and endovascular techniques (e.g. snaring, sheath removal)
3. Close arterial access according to standard practice.
4. Follow-up patients as necessary to provide proper surveillance of the long-term performance of the endoprosthesis, procedure and status of the aneurysm. Minimally, annual CTs, multiple view X-rays, and ultrasound may be used for such surveillance.

IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP

General

The long-term safety and effectiveness of endovascular repair has not been established. All patients should be advised this treatment modality requires long-term, regular follow-up to assess patients' health status and stent-graft performance. Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms) should receive enhanced follow-up. Patients should be counseled on the need for regular follow-up, even in the absence of obvious symptoms (e.g., pain, numbness, weakness).

Regular and consistent follow-up is a critical part of ensuring continuing safety and efficacy of aortic endovascular repair. Physicians should tailor patient follow-up to the needs and circumstances of each individual patient. In the US clinical studies, at least one annual physician visit and the imaging schedule (Table 11) were employed.

Follow-up modalities include CT / CTA, multi-view abdominal X-ray, MRI / MRA, and ultrasound. Data from these modalities is acquired and used to compare baseline and subsequent exams to review devices and morphological changes over time and their effects on exclusion of the aneurysm.

- CT / CTA imaging provides information on aneurysm size, vascular morphological changes, proximal device-trunk fixation and migration, endoleak and patency / limb occlusion.
- Multi-view device X-ray imaging provides information on device wireform integrity (e.g., fracture, kinking) and relative component migration.
- MRI / MRA imaging provides information similar to CT / CTA and is often used as a surrogate for CT / CTA if patients are CT contrast intolerant.
- Ultrasound may be used to assess for endoleak and aneurysm size status but not for device integrity, specifically the wire form. Ultrasound is a less reliable and sensitive diagnostic method compared to CT.

Alternative imaging recommendations for patients with CT or angiography contrast intolerance issues include CO₂ angiography, MRI-MRA with or without contrast, and ultrasound. These imaging and surveillance modalities may be less sensitive and difficult to compare with diagnostic findings from previous or subsequent follow-up exams.

Table 11. Recommended Schedule for Patient Imaging Follow-Up

Visit	Schedule for Patient Imaging Follow-Up		
	Angiogram	Abdominal X-ray	CT Pre-Contrast and Contrast
Pre-Treatment	X ¹		X ¹
Treatment (Pre and Post Deployment)	X		
Discharge		X	
1 Month			X
3 Month			X ²
6 Month		X	X
12 Month (Annually Thereafter)		X	X

¹ Imaging should be performed ≤ 6 months prior to the procedure

² Recommended if endoleak reported at one month

Angiographic Imaging

Angiographic images are recommended pre-treatment to evaluate the length and tortuosity of abdominal aorta, iliac and common femoral arteries.

- Images should include an angiographic marker catheter with incremental one centimeter markers over a 10 – 20 cm length.
- The following views are recommended for optimal evaluation and case planning:
 - Abdominal aorta; Supine-AP, Lateral
 - Pelvis (to include bilateral common femorals); AP, both Obliques

Angiographic images are recommended during the treatment procedure both pre and post-deployment to evaluate device placement and orientation. Selective angiography during subsequent follow-up exams may provide useful device position and device integrity information.

CT / CTA Images

- Film sets should include all sequential images at lowest possible slice thickness (≤ 3 mm). Do NOT perform large slice thickness (> 3 mm) and / or omission of CT images / film sets (non-consecutive) as it prevents precise anatomical and device comparisons over time.
- All images should include a scale for each film / image. Images should be arranged no smaller than 20:1 images on 14" x 17" sheets if film is used.
- **If an endoleak is suspected or there is aneurysm enlargement, it is recommended that pre-contrast and contrast runs be performed.**
- Pre-contrast and contrast run slice thickness and interval must match.
- DO NOT change patient orientation or re-landmark patient between non-contrast and contrast runs.

Non-contrast and contrast enhanced baseline and follow-up exams are important for optimal patient surveillance. For the best results, use the following CT / CTA imaging guidelines listed in Table 12.

Table 12. CT / CTA Imaging Guidelines

	Pre-Contrast	CT / CTA
IV Contrast	No	Yes
Injection Volume (ml)	NA	150
Injection Rate (cc / sec)	NA	≥ 2.5
Delay	NA	Smart-Prep [®] , CARE or equivalent
Start Position	Diaphragm	1 cm above Celiac Axis
End Position	Proximal Femur	Femoral Bifurcation
Scan FOV	Large	Large
DFOV	32 cm	32 cm
Scan Type	Helical	Helical
Rotation Speed	0.8	0.8
Slice Thickness (mm)	≤ 3.0 mm	≤ 3.0 mm
Scan Mode	HS	HS
Table Speed (mm / rot)	15	15
Interval (mm)	2.0	2.0
KV / mA	120 / 300	120 / 300
Reconstruction / Algorithm	≤ 3.0 mm Soft	≤ 3.0 mm Soft
* Smart Prep ROI Loc: 1 cm Sup. to Celiac Axis Scan Phase: 3 Sec MA: 40 Monitor Delay: 6 Sec Monitor ISO: 3 Sec Enhance Thres: 100 HU		

Abdominal X-ray Film Series (plain film)

The following abdominal X-ray views are recommended for optimal visualization of the engraft:

- Supine – Frontal (AP)
- Lateral
- 30 degree LPO
- 30 degree RPO

Ensure entire device is captured on each single image format lengthwise.

If there is any concern about the device integrity (e.g., kinking, stent-wire breaks, relative component migration), it is recommended to use magnified views. The attending physician should evaluate films for device integrity (entire device length including components) using 2 – 4x magnification.

MRI Information MR Conditional

Non-clinical testing has demonstrated that the GORE® EXCLUDER® AAA Endoprosthesis is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Highest spatial gradient magnetic field of 3000-Gauss/cm or less
- Whole body-averaged specific absorption rate (SAR) of 2.0-W/kg in the Normal Operating Mode (i.e., the mode of operation of the MR equipment in which none of the outputs have a value that cause physiological stress to patients) for 15 minutes of scanning (i.e., per pulse sequence).

MR Related Heating:

In non-clinical testing, the GORE® EXCLUDER® AAA Endoprosthesis produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numeris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and during MRI performed for 15-min of scanning 3-Tesla/128-MHz (Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

Please see TABLE 1. MRI Information

Image Artifact:

For each vascular device and assembly, the artifacts that appeared on the MR images were shown as localized signal voids (i.e., signal loss) that were minor in size relative to the size and shape of these implants. The gradient echo pulse sequence produced larger artifacts than the T1 – weighted, spin echo pulse sequence for the GORE® EXCLUDER® AAA Endoprosthesis. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 10 mm relative to the size and shape of this implant. MR image quality may be

compromised if the area of interest is in the exact same area or relatively close to the position of the GORE® EXCLUDER® AAA Endoprosthesis. Therefore, it may be necessary to optimize the MR imaging parameters to compensate for the presence of this implant.

Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement ≥ 5 mm of maximum diameter (regardless of endoleak status)

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy, and the patient's personal choices. Patients should be counseled as to the possibility of subsequent reinterventions including catheter-based and open surgical conversion.

DEVICE-RELATED ADVERSE EVENT REPORTING

Any adverse event involving the GORE® EXCLUDER® AAA Endoprosthesis should be reported to W. L. Gore & Associates immediately. To report an event in the US, call 800.437.8181. Outside the US, contact your local technical representative.

PATIENT TRACKING INFORMATION

In addition to these Instructions for Use, the GORE® EXCLUDER® AAA Endoprosthesis is packaged with a Device Tracking Form which US hospital staff are required to complete and forward to Gore for the purposes of tracking all patients who receive a GORE® EXCLUDER® AAA Endoprosthesis product (as required by US Federal Regulation).

DEFINITIONS

Use By

Caution

Consult Instructions for Use

Do Not Resterilize

Do Not Reuse

Catalogue Number

Batch Code

Serial Number

Authorised Representative in the European Community

MR Conditional

CAUTION: USA Federal Law restricts the sale, distribution, or use of this device to, by, or on the order of a physician.

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Manufacturer

W. L. GORE & ASSOCIATES, INC.1505 North Fourth Street
Flagstaff, Arizona 86004**UNITED STATES**

Order Information: Tel.: +1.928.526.3030 • Tel.: +1.800.528.8763

Technical Information: Tel.: +1.928.779.2771 • Tel.: +1.800.437.8181

For international contact and additional product information,
visit **www.goremedical.com****EC REP** Authorised Representative in
the European Community**W. L. GORE & ASSOCIATES B.V.**

Ringbaan Oost 152A

5013CE Tilburg

NETHERLANDS

MADE IN USA

GORE®, EXCLUDER®, C3®, and designs are trademarks of W. L. Gore & Associates.

All other trademarks are the property of their respective owners.

© 2002-2004, 2007-2012, 2014, 2015, 2017, 2019 W. L. Gore & Associates, Inc.

Printed on recyclable paper

DECEMBER 2019

20062333